

# Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS	C1
EQUIPMENT LAYOUT (Equipment locations, heat loads, component weights, environmental specs)	A1
STRUCTURAL LAYOUT (Structural support/mounting locations for floor/wall/ceiling, wall support elevations)	S1
STRUCTURAL DETAILS (Floor and Ceiling loading information)	S2
ELECTRICAL LAYOUT (Contractor supplied wiring, interconnect methods, junction point locations and descriptions)	E1
ELECTRICAL SPECIFICATIONS (Maximum wiring run lengths, interconnect diagram, system power specifications)	E2
ELECTRICAL DETAILS	E3 THRU E4
EQUIPMENT DETAILS	D1 THRU D2

These drawings indicate the planned location and orientation of the listed equipment and the listed equipment components. These drawings are prepared on a non-warranty basis. GE Healthcare is not responsible for ultimately preparing the site for the operation of this equipment in operation. GE Healthcare complies with applicable federal, state, and/or local requirements.

## \* REQUIRED REFERENCE \*

INNOVA PLUS  
Pre Installation Manual  
5314528-4-1EN

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the pre installation manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

[www.gehealthcare.com/siteplanning](http://www.gehealthcare.com/siteplanning)

# GE Healthcare



## Cardio-Vascular Site Planning

CUSTOMER ACCEPTANCE



imagination at work

# Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

# GE Equipment Delivery Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist Rev 19				
Before using this document ensure you have the latest Rev from MyWorkshop on D00422752				
GEHC Global Order #:		Customer:		
GEHC PMI:		FE / Installer:		
The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.				
Inspection Date:	Storage is item ready?	PM is item ready?	FE is item ready?	Comments If 'N', enter comments or action plan
1				<b>MR Magnet Delivery Requirements:</b> Ensure oxygen venting system is available for magnet connection as defined by GEHC Pre-Installation Manual (PIM) requirements, exhaust fan system is installed and operational, 480V power, and chilled water supply is available 24x7 that meets system cooling requirements. External connectivity is available for magnet monitoring and phone service is available during delivery. Surface mount vitromat installed where required. Magnet room final flooring is in place.
2				<b>MR RF Screen Room Requirements:</b> RF Screen Room is tested with copy of Test Report, emailed to ISAdminCOEM@ge.com, that it is compliant with GEHC specifications. Dock Bolt and magnet anchors (if applicable) installed using 2 part anchor. For HDx systems, blower box mount bolts installed by RF vendor using 2 part anchors.
3				<b>State Regulatory Requirements:</b> Facility registration number provided for states of IL, KY, HI, RI, SC, TX. X-ray shielding plan and state acknowledgment letter provided to installer for AR, DC, NC, SC, CO & WA. <b>Site Drawing Requirements:</b> Final version of equipment network and antenna, installation drawings (including red lined versions) verified to match actual room and has been provided to installer.
4				<b>Surface Penetration Requirements:</b> Customer/Contractor scheduled to provide required drilling or cutting into floors, ceilings, and walls, OR surface penetration permit available and posted in the room when GEHC will perform the work.
5				<b>Pre-Delivery Route Requirements:</b> The equipment delivery route from the truck to the final destination within the facility has been reviewed with all key stakeholders to safely meet the minimum requirements for equipment access, and all communications/notifications have occurred. Arrangements have been made for special handling (elevator, rigging, floor protection, fork lift, rollback truck, etc).
6				<b>Finished Room Requirements:</b> Rooms that will contain equipment, including storage areas not in scan suite, are dust free. Provisions taken to maintain a dust free room. Precautions must be taken to prevent dust from entering rooms containing equipment when construction is incomplete in adjacent areas. All walls primed (final coat not needed on Day 1). Shielding, doors, and windows are to be installed. No contractor work being done during or after the installation that will cause dust in the installation areas or potential equipment damage. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of these security issues, implications and responsibility. For Storage: Room must meet PIM requirements for storage.
7				<b>Electrical Requirements:</b> Lockable (LOTO) Main Disconnect Panel (MDP) is installed per GE guidelines and system power is available. Conduits, electrical cable ducting/dividers/cable trays, and access flooring is installed in proper location and height. Surface floor duct and load-side wires can be installed at time of system installation. Validate outlet location and requirements meet specifications for device/equipment.
8				<b>HVAC Requirements:</b> The HVAC/Chilled Water systems designed to maintain the environment per spec/PM is at running state and appears to provide the desired environmental conditions including location of vents, temperature and humidity for system operation.
9				<b>Flooring Requirements:</b> Floor is clean and prepared for final floor covering. Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications. Confirm customer anchoring plan aligns with designed floor thickness. Final flooring installed where required for network racks.
10				<b>Ceiling Requirements:</b> Unistrut (or equivalent) location, levelness and spacing is measured (or vendor confirmed) and consistent with the requirement of the installation drawings. Ensure unistrut and rails are not used as mounting surfaces. Ceiling grid is installed. Remount lighting is installed and operational. HVAC diffusers are installed and connected to ductwork. Ceiling tiles installed per PMI discretion.

GE Healthcare  
[DesignCenter@gehealthcare.com](mailto:DesignCenter@gehealthcare.com)  
 Design Center  
 Milwaukee, WI  
 Copyright 2009 General Electric Company - Proprietary to GE

SHEET TITLE: SITE READINESS  
 MODALITY TYPE: INNOVA 2100/3100/4100  
 THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO ALL APPLICABLE REGULATORY REQUIREMENTS AND TO THE COMPANY'S POLICY OF ACTUAL CONSTRUCTION PURPOSES AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: HYBRID OR  
 TYPICAL FINAL LAYOUT

PROJECT	REVISION
5-94F	00
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST

REVISION HISTORY:


SHEET  
C1

PIM 901  
RQ - 140195


GE EQUIPMENT LISTING								
EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER				EQUIPMENT CROSS REFERENCE CHART				
ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	STRC PLAN	ELEC PLAN
1	1		DETECTOR CHILLER	33 lbs	706 btu	B5049F	---	DC S
2	1		WORKSTATION CART				---	---
3	1		MAC LAB CONSOLE, INCLUDES MONITORS AND KEYBOARD	181 lbs	2935 btu		---	PC S
4	1		TRAM NET RACK	8 lbs		B5047	---	TRAM S
5	1		REMOTE CONTROL FOR INJECTOR	4 lbs		B5028	---	IEC S
6	1		INJECTOR HEAD ON TABLE RAIL	15 lbs		B5030A	---	IH S
7	1		INJECTOR ELECTRONICS	37 lbs	320 btu	B5028	---	IE S
8	1		OMEGA IV/V TABLE WITH ROTATING TOP	1750 lbs	614 btu	B5049N	---	LUS C
9	1		XR BUZZER (LOCATED ABOVE CEILING)	2 lbs		B5150H	---	XR8 -
10	1		INNOVA POSITIONER (REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION)	1653 lbs	2416 btu	B5050A B5050B B5050	---	LC1 C
11	1		UPS INTERFACE BOX			E45021B	---	UTB -
12	1		ATLAS CABINET (C2)	659 lbs	1825 btu	B0558C	S100	C2 C
13	1		ATLAS CABINET (C1)	1115 lbs	3389 btu	B0558C	S100	C1 C
14	1		UPS CABINET	1170 lbs	4061 btu	E45025C	---	UPS -
15	1		WATER CHILLER	449 lbs	18716 btu	M0917B	---	CHLR C
16	1		IVUS VOLCANO S51 CONSOLE, INCLUDES FLAT PANEL MONITOR AND KEYBOARD (DESK MOUNTED)	68 lbs	1631 btu	B551	---	IVUS -
17	1		IVUS VOLCANO COLOR PRINTER					---
18	1		18 IN. MONITOR ON WALL SUPPORT	26 lbs	204 btu	C7617B	---	WBM4 C
19	1		CONTROL ROOM MONITOR WITH DL KEYPAD	22 lbs	204 btu	C7412H C7617	---	WBM2 S
20	1		OPERATORS CONSOLE	22 lbs	546 btu	C7617 C7558 B5050C	---	WBC1 C
21	1		BOLUS CHASE HANDSWITCH	2 lbs			---	WBBC -
22	1		AW WORKSTATION	81 lbs	1201 btu	M1013AM C7617	---	WBM3 C
23	1		FOUR MONITORS & SURGICAL LAMP SUSPENDED FROM ARM/ CEILING SUPPORT				---	WBM1 -
24	1		TWO MONITORS & TWO SURGICAL LAMPS SUSPENDED FROM ARM/ CEILING SUPPORT				---	WBM2 -
25	1		TWO MONITORS, RADIATION SHIELD & SURGICAL LAMP, SUSPENDED FROM ARM/ CEILING SUPPORT				---	WBM3 -

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

50	1		INNOVA MAIN DISCONNECT, REFERENCE JUNCTION POINT "PDB" ON SHEET E1 FOR DETAILED DESCRIPTION.	326 lbs	1532 btu	E4502M	---	PDB -
----	---	--	--	---------	----------	--------	-----	-------

SCALE: 1/4" = 1'-0" EQUIPMENT LAYOUT RECOMMENDED CEILING HEIGHT = 9'-6"

This equipment layout indicates the placement and interconnection of the indicated equipment components. There may be federal, state, and/or local requirements that could impact the placement of these components. It remains the Customer's responsibility for ensuring the site and final equipment placement complies with all applicable federal, state, and/or local requirements.



ANCILLARY ITEMS	
CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS	
ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
60	COUNTER TOP WITH SINK, BASE AND WALL CABINETS
61	MED GASES IN CEILING
62	LEAD APRON RACK
63	COUNTER TOP WITH BASE AND WALL CABINETS
64	COUNTER TOP FOR EQUIPMENT - MINIMUM DEPTH 30 IN. OR ADDITIONAL SHELVING MAY BE REQUIRED. PROVIDE GRIMMETTED OPENINGS AS REQUIRED TO ROUTE INTERCONNECT CABLES TO RACEWAY BELOW COUNTER TOP.
65	CONTROL WALL TO CEILING WITH LEAD GLASS WINDOW
66	CATHETER CABINETS
67	SHELF - CUSTOMER TO PROVIDE ADEQUATE WALL SUPPORT
68	X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL # 800-800-9760
69	GE CAT. NO. WX1ABW-DF-XIU MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W X 83 IN. H (118mm X 2108mm), CONTINGENT ON A 96 IN. (2438mm) CORRIDOR WIDTH
70	150-AMP DISCONNECT BREAKER (FOR LOCK-OUT/TAG-OUT CAPABILITY)
71	CUSTOMER SUPPLIED STORAGE CABINET
72	RACK FOR HYBRID NETWORK ITEMS
73	PROCEDURAL CAMERA
74	PROCEDURAL VIDEO MONITOR
THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.	
90	X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL REFERENCE JUNCTION POINT "XRLC" ON SHEET E1 FOR DETAILED DESCRIPTION - CAT. NO. E4502SS FOR WARNING LIGHT & ROOM LIGHT CONTROL.

- ### GENERAL SPECIFICATIONS
- THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.
  - CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.
  - RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.
  - THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PREDICATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC., MAKE THIS ESSENTIAL FOR A PROPER IS. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC..
  - ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.
  - DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.

- ### SITE ENVIRONMENT SPECIFICATIONS
- AMBIENT OPERATING TEMPERATURE: EQUIPMENT ROOM WITH FLUORO UPS OPTION 68° TO 77° F. (20° TO 25° C)
  - AMBIENT OPERATING TEMPERATURE: CONTROL ROOM 68° TO 77° F. (20° TO 25° C)
  - AMBIENT OPERATING TEMPERATURE: EXAM ROOM-DESIGN FOR PATIENT/OPERATOR COMFORT TARGET TEMPERATURE 64° F (18° C)
  - HUMIDITY: 30° TO 75° FOR EQUIPMENT AND CONTROL ROOMS AND 30° TO 70° FOR EXAM ROOM
  - ALTITUDE: NOT TO EXCEED 9,842 FT. (3000M) ABOVE SEA LEVEL.
  - THE ENVIRONMENT FOR THE ELECTRONICS CABINET MUST BE CONTROLLED SO THE ABOVE RESTRICTIONS ARE NOT EXCEEDED.
  - DO NOT RESTRICT THE AIR INTAKE OR AIR EXHAUST OF THE SYSTEM COMPONENTS.
  - ENVIRONMENTAL CONDITIONS LISTED ABOVE MUST BE MAINTAINED AT ALL TIMES INCLUDING FOR EXAMPLE OVERNIGHT, WEEKENDS, AND HOLIDAYS.

- ### MAGNETIC INTERFERENCE SPECIFICATIONS
- IMAGE INTENSIFIERS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 1 GAUSS TO GUARANTEE SPECIFIED IMAGING PERFORMANCE.
  - X-RAY TUBES MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE SPECIFIED PERFORMANCE.
  - SYSTEM ELECTRONICS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE DATA INTEGRITY.
  - OPERATORS CONSOLE EQUIPMENT MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO OBTAIN SPECIFIED GEOMETRIC LINEARITY.

GE Healthcare  
IS Services Design Center  
Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT LAYOUT  
MODALITY TYPE: INNOVA 2100/3100/4100  
THIS PLAN IS SUBMITTED TO ASSIST IN THE LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE LOCAL, STATE AND FEDERAL BUILDING CODES AND REGULATIONS. GE HEALTHCARE ACCEPTS NO LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
HYBRID OR  
TYPICAL FINAL LAYOUT

PROJECT	REVISION
5-94f	00
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST


REVISION HISTORY:


SHEET  
A1

TYPICAL WALL SUPPORT ELEVATIONS


S107

S115



SUPPORT FOR MAIN DISCONNECT CONTROL


(NOT TO SCALE)



SUPPORT FOR UPS INTERFACE BOX

(NOT TO SCALE)

S100




SUPPORT FOR ATLAS/SYSTEMS CABINET

(NOT TO SCALE)

SCALE: 1/4" = 1'-0"

STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-6"



STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S107, FOR MAIN DISCONNECT CONTROL.
2	AREA OCCUPIED BY GE SUPPLIED TABLE BASEPLATE
3	MOUNT XR BUZZER BRACKET ON WALL ABOVE CEILING
4	AREA OCCUPIED BY GE SUPPLIED POSITIONER BASEPLATE
5	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S115, FOR UPS INTERFACE BOX.
6	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S100, FOR ATLAS CABINET.
7	CEILING SUPPORT FOR SURGICAL LAMP/ MONITOR POST
8	AREA OF MED GASES IN CEILING

STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED TUBE HANGER OR OTHER EQUIPMENT ARE TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS. THE UNISTRUT OR EQUIVALENT STRUCTURE SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL, SQUARE AND IN THE SAME HORIZONTAL PLANE FLUSH WITH FINISHED CEILING. THE SYSTEM IS TO BE CROSS BRACED VERTICALLY, HORIZONTALLY AND DIAGONALLY TO ALLOW NO MOVEMENT AND A MAXIMUM OF 1.58mm (1/16") DEFLECTION. (10) 12.7mm (1/2") DIA. x 38.1mm (1 1/2") LONG BOLTS WITH UNISTRUT 12.7mm (1/2") NUTS WITH SPRINGS ARE TO BE PROVIDED BY CUSTOMER OR HIS CONTRACTORS FOR EACH STATIONARY AND AUXILIARY SUPPORT RAIL. CLOSURE STRIPS SHALL BE PROVIDED FOR AREAS OF UNISTRUT EXPOSED AND WITHOUT MOUNTING UNITS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 6.35mm (1/4") BELOW THE FINISHED CEILING.
- CONTROL WALLS WITH TUBE HANGER PASSAGE ABOVE SHALL BE CONSTRUCTED TO 2130mm (7'-0") HIGH.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 3.17mm (1/8") IN 3050mm (10'-0")
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.
- IT IS THE CUSTOMER'S RESPONSIBILITY TO PERFORM ANY FLOOR OR WALL PENETRATIONS THAT MAY BE REQUIRED. THE CUSTOMER IS ALSO RESPONSIBLE FOR ENSURING THAT NO SUBSURFACE UTILITIES (E.G., ELECTRICAL OR ANY OTHER FORM OF WIRING, CONDUITS, PIPING, DUCT WORK OR STRUCTURAL SUPPORTS (I.E. POST TENSION CABLES OR REBAR)) WILL INTERFERE OR COME IN CONTACT WITH SUBSURFACE PENETRATION OPERATIONS (E.G. DRILLING AND INSTALLATION OF ANCHORS/SCREWS) PERFORMED DURING THE INSTALLATION PROCESS. TO ENSURE WORKER SAFETY, GE INSTALLERS WILL PERFORM SURFACE PENETRATION OPERATIONS ONLY AFTER THE CUSTOMER'S VALIDATION AND COMPLETION OF THE "GE SURFACE PENETRATION PERMIT"

**GE Healthcare**  
 IS Services Design Center  
 Milwaukee, Wisconsin

SHEET TITLE: **STRUCTURAL LAYOUT**  
 MODALITY TYPE: **INNOVA 2100/3100/4100**  
 THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO ALL APPLICABLE REGULATORY REQUIREMENTS AND TO THE COMPANY'S POLICY OF USING ONLY ACTUAL CONSTRUCTION METHODS AND MATERIALS. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: **HYBRID OR**  
 TYPICAL FINAL LAYOUT

PROJECT	REVISION
5-94f	00

DATE: 18.Dec.13  
 DRAWN BY: JPH  
 CHECKED BY: TST


REVISION HISTORY:


SHEET  
**S1**

FLOOR MOUNTING : INNOVA 2100-3100-4100 (UNITY)/OMEGA V LONG TABLE (WITH IQ TILT TABLE BASEPLATE) INSTALLATION (TEMPLATE NO. 2360133)

B5049N


REV. DATE: 06/04/09




**WARNING!! THE RELATIONSHIP BETWEEN THE TABLE BASE AND THE POSITIONER BASEPLATE IS CRITICAL.**

PRIOR TO DRILLING MOUNTING HOLES CONTACT LOCAL GE HEALTHCARE INSTALLATION PROJECT MANAGER OR LEAD FIELD ENGINEER TO VERIFY THAT THE PROPER FULL SIZE FLOOR MOUNTING TEMPLATE IS USED.


MEDICAL GAS FLOOR EXIT LOCATIONS



THROUGH-BOLT MOUNTING OPTIONS



ANCHOR BOLT MOUNTING OPTIONS



**Customer/Contractor Alert:** It is the responsibility of the Customer or their Contractor to drill all anchor/thru-bolting holes for anchoring the positioner and table to the floor. Refer to GEHC document no. \*2290880-2-100 for installation preparation and procedures.

**NOTE: THRU BOLTING IS HIGHLY PREFERRED FOR THE INSTALLATION OF THE POSITIONER BASEPLATE AND OMEGA TABLE.** HARDENED BOLTS AND 4" x 4" [102mm x 102mm] STEEL PLATES TO BE USED ARE SUPPLIED BY GE HEALTHCARE AS INDICATED ON THE ACTUAL DETAIL DRAWING. BE ADVISED, HOWEVER, THAT ADDITIONAL SUPPORT STRUCTURES: STEEL BEAMS, PLATES, CORE BORING OF MOUNTING HOLES, ETC., ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

**NOTE: IF THRU BOLTING IS NOT POSSIBLE, FLOOR ANCHORS CAN BE USED IF APPROVED BY CUSTOMERS STRUCTURAL ENGINEER.** FOR ON GRADE INSTALLATIONS, MOUNTING KIT CAT. NO. **2286398** SHOULD BE ORDERED. ANCHORS INCLUDED IN KIT SHOULD BE APPROVED BY CUSTOMERS STRUCTURAL ENGINEER.

**NOTE: BASEPLATES MUST BE LEVEL WITHIN 1/32" [0.79mm]**

**NOTE: JOISTS MUST BE SPANNED WITH STEEL REINFORCING.** SIZE AND THICKNESS OF STEEL REINFORCING ARE DETERMINED BY THE ACTUAL PAN CONSTRUCTION ON SITE. STEEL PLATES, CHANNELS OR BEAMS MAY BE USED.

**NOTE: DETERMINE THE POSITION OF THE "REBARS" IN THE CONCRETE FLOOR SO ANCHOR HOLES WILL NOT RUN INTO THEM.**

**\* DOCUMENT FURNISHED BY GEHC AS PART OF INSTALLATION KIT**

POSITIONER BOLT FORCES FOR WORST CASE CONDITIONS

LOADS		BOLT TENSION (AT BOLT "A")	
HORIZONTAL ACCELERATION = 625 lbs. [284 Kg]	VERTICAL ACCELERATION = 209 lbs. [95 Kg]	MAXIMUM TENSION = 881 lbs. [400 Kg]	BOLT SHEAR (U-ARM LOCKED) MAXIMUM SHEAR = 120 lbs. [54 Kg]/BOLT

OMEGA TABLE BOLT FORCES FOR WORST CASE CONDITIONS

LOADS		BOLT TENSION		BOLT SHEAR	
		MAXIMUM TENSION = 1938 lbs. [880 Kg]/BOLT		MAXIMUM TENSION = 407 lbs. [185 Kg]/BOLT	MAXIMUM SHEAR = 407 lbs. [185 Kg]/BOLT

SHEET TITLE: STRUCTURAL DETAILS  
MODALITY TYPE: INNOVA 2100/3100/4100

THIS PLAN IS SUBMITTED TO SUPPORT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE MANUFACTURER'S DRAWINGS AND TO THE CODES AND STANDARDS. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
HYBRID OR  
TYPICAL FINAL LAYOUT

PROJECT	REVISION
5-94f	00

DATE: 18.Dec.13  
DRAWN BY: JPH  
CHECKED BY: TST



REVISION HISTORY:

SHEET  
S2



INTERCONNECT DIAGRAM

DATED: 10/27/09



POWER SPECIFICATIONS

INNOVA SYSTEMS  
REV. DATE: 01/04/07

VOLTAGE: PRIMARY SOURCE IS REQUIRED FOR ALL INSTALLATIONS.  
RANGE OF LINE VOLTAGES: NOMINAL LINE VOLTAGE OF 360 TO 480, 3 PHASE, 50 OR 60 Hz  
REQUIRED POWER SUPPLY: WYE DISTRIBUTION

MAXIMUM DAILY VOLTAGE VARIATION MUST FALL WITHIN ONE OF THE RANGES IN TABLE A.

**TABLE A ALLOWABLE INPUT VOLTAGES/CURRENT DEMAND**

NOMINAL VOLTAGE	NORMAL RANGE ±10 PERCENT	CURRENT (AMPS)	
		MAX. MOMENTARY	CONTINUOUS
360	324-396	304	32
380	342-418	289	31
400	360-440	274	29
420	378-462	264	28
440	396-484	249	26
460	414-506	238	25
480	432-528	228	24

ALL CALCULATIONS BASED UPON NOMINAL VOLTAGE

NOTE: LOW LINE CONDITIONS MAY INHIBIT SOME HIGH kVp TECHNIQUES. THE GENERATOR AUTOMATICALLY ESTABLISHES THESE INHIBITS BASED ON ACTUAL LINE CONDITIONS AND SYSTEM REGULATION.

PHASE-BALANCE: PHASE-TO-PHASE VOLTAGES MUST BE WITHIN +2 PERCENT OF THE LOWEST PHASE-TO-PHASE VOLTAGE. MAXIMUM ALLOWABLE TRANSIENT VOLTAGE EXCURSIONS ARE 2.5 PERCENT OF RATED LINE VOLTAGE AT A MAXIMUM DURATION OF 5 CYCLES AND FREQUENCY OF 10 TIMES PER HOUR.

POWER DEMAND: CONTINUOUS POWER DEMAND = 20KVA. (MAX DEMAND = 171 KVA)

**TABLE B MAXIMUM MOMENTARY POWER DEMAND.**

DEMAND	ADVANTX 100
kVa * POWER FACTOR AT	171 0.9
mA	1250
kVp	80

\* DEMAND INCLUDES POWER FOR ENTIRE ADVANTX SYSTEM. LINE VOLTAGE REGULATION AT MAXIMUM POWER DEMAND MUST BE LESS THAN OR EQUAL TO 6 PERCENT.

DISTRIBUTION TRANSFORMER: FOR A SINGLE UNIT INSTALLATION, THE MINIMUM TRANSFORMER SIZE IS 225 KVA.

ELECTRICAL NOTES

- NOTE 1: ALL WIRES SPECIFIED SHALL BE COPPER STRANDED, FLEXIBLE, THERMO-PLASTIC, COLOR CODED, CUT 10 FOOT LONG AT OUTLET BOXES, DUCT TERMINATION POINTS OR STUBBED CONDUIT ENDS. ALL CONDUCTORS, POWER, SIGNAL AND GROUND, MUST BE RUN IN A CONDUIT OR DUCT SYSTEM. ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS. WIRE RUNS MUST BE CONTINUOUS COPPER STRANDED AND FREE FROM SPLICES. ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.
- NOTE 2: WIRE SIZES GIVEN ARE FOR USE OF EQUIPMENT. LARGER SIZES MAY BE REQUIRED BY LOCAL CODES.
- NOTE 3: IT IS RECOMMENDED THAT ALL WIRES BE COLOR CODED, AS REQUIRED IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 4: CONDUIT SIZES SHALL BE VERIFIED BY THE ARCHITECT, ELECTRICAL ENGINEER OR CONTRACTOR, IN ACCORDANCE WITH LOCAL OR NATIONAL CODES.
- NOTE 5: CONVENIENCE OUTLETS ARE NOT ILLUSTRATED. THEIR NUMBER AND LOCATION ARE TO BE SPECIFIED BY OTHERS. LOCATE AT LEAST ONE CONVENIENCE OUTLET CLOSE TO THE SYSTEM CONTROL, THE POWER DISTRIBUTION UNIT AND ONE ON EACH WALL OF THE PROCEDURE ROOM. USE HOSPITAL APPROVED OUTLET OR EQUIVALENT.
- NOTE 6: GENERAL ROOM ILLUMINATION IS NOT ILLUSTRATED. CAUTION SHOULD BE TAKEN TO AVOID EXCESSIVE HEAT FROM OVERHEAD SPOTLIGHTS. DAMAGE CAN OCCUR TO CEILING MOUNTING COMPONENTS AND WIRING IF HIGH WATTAGE BULBS ARE USED. RECOMMEND LOW WATTAGE BULBS NO HIGHER THAN 75 WATTS AND USE DIMMER CONTROLS (EXCEPT MR). DO NOT MOUNT LIGHTS DIRECTLY ABOVE AREAS WHERE CEILING MOUNTED ACCESSORIES WILL BE PARKED.
- NOTE 7: ROUTING OF CABLE, DUCTWORK, CONDUITS, ETC., MUST RUN DIRECT AS POSSIBLE OTHERWISE MAY RESULT IN THE NEED FOR GREATER THAN STANDARD CABLE LENGTHS (REFER TO THE INTERCONNECTION DIAGRAM FOR MAXIMUM USABLE LENGTHS POINT TO POINT).
- NOTE 8: CONDUIT TURNS TO HAVE LARGE, SWEEPING BENDS WITH MINIMUM RADIUS IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 9: A SPECIAL GROUNDING SYSTEM IS REQUIRED IN ALL PROCEDURE ROOMS BY SOME NATIONAL AND LOCAL CODES. IT IS RECOMMENDED IN AREAS WHERE PATIENTS MIGHT BE EXAMINED OR TREATED UNDER PRESENT, FUTURE, OR EMERGENCY CONDITIONS. CONSULT THE GOVERNING ELECTRICAL CODE AND CONFER WITH APPROPRIATE CUSTOMER ADMINISTRATIVE PERSONNEL TO DETERMINE THE AREAS REQUIRING THIS TYPE OF GROUNDING SYSTEM.
- NOTE 10: THE MAXIMUM POINT TO POINT DISTANCES ILLUSTRATED ON THIS DRAWING MUST NOT BE EXCEEDED.
- NOTE 11: PHYSICAL CONNECTION OF PRIMARY POWER TO GE EQUIPMENT IS TO BE MADE BY CUSTOMERS ELECTRICAL CONTRACTOR WITH THE SUPERVISION OF A GE REPRESENTATIVE. THE GE REPRESENTATIVE WOULD BE REQUIRED TO IDENTIFY THE PHYSICAL CONNECTION LOCATION, AND INSURE PROPER HANDLING OF GE EQUIPMENT.

**GE Healthcare**  
IS Services Design Center  
Manufacturer  
Wisconsin













SHEET TITLE: ELECTRICAL SPECIFICATIONS  
MODALITY TYPE: INNOVA 2100/3100/4100  
THIS PLAN IS SUBMITTED TO ASSIST IN THE LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE CODES AND STANDARDS. THE USER OF THIS PLAN SHALL BE RESPONSIBLE FOR VERIFYING THE ACCURACY OF ALL INFORMATION AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: HYBRID OR TYPICAL FINAL LAYOUT

PROJECT: 5-94f  
REVISION: 00  
DATE: 18.Dec.13  
DRAWN BY: JPH  
CHECKED BY: TST


REVISION HISTORY:

SHEET E2



ELECTRICAL DETAIL  
BOX WITH COVERPLATE (TYPICAL)


ELEC-8  
REV. DATE: 09/30/94



DETAIL NOT TO SCALE


ELECTRICAL DETAIL  
POSITIONER INTERCONNECT DETAIL, UNDER FLOOR

ELEC-100  
REV. DATE: 03/30/04




ELECTRICAL DETAIL  
J.B. / WALL DUCT DETAIL (TYPICAL)

ELEC-2  
REV. DATE: 09/30/94



ELECTRICAL DETAIL  
CONDUITS THRU-FLOOR (TYPICAL)

ELEC-9  
REV. DATE: 08/08/94



GE Healthcare  
IS Services Design Center  
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL DETAILS  
MODALITY TYPE: INNOVA 2100/3100/4100

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST NATIONAL ELECTRICAL CODES AND THE NATIONAL ELECTRICAL CODE FOR CONSTRUCTION PURPOSES. HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
HYBRID OR  
TYPICAL FINAL LAYOUT

PROJECT	REVISION
5-94f	00

DATE: 18.Dec.13  
DRAWN BY: JPH  
CHECKED BY: TST













REVISION HISTORY:

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

SHEET  
E4

PIM R01  
RQ - 140195





PIM ROT

RQ - 140195


THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

EQUIPMENT DETAIL  
COOLIX 4000 RECIRCULATING CHILLER

M0917B

REV. DATE: 05/17/05


SHIPPING DIMENSIONS:  
41"(1040mm)D x 34.8"(870mm)W x 53"(1350mm)H



EQUIPMENT DETAIL  
IVUS VOLCANO S5i WORKSTATION

BS5I


REV. DATE: 04/04/08



EQUIPMENT DETAIL  
18" FLAT PANEL MONITOR ON WALL SUPPORT

C76-17B


REV. DATE: 04/29/04



EQUIPMENT DETAIL  
DLX or DL KEYPAD

C7412H


REV. DATE: 09/03/03



EQUIPMENT DETAIL  
18" FLAT PANEL MONITOR

C76-17


REV. DATE: 08/28/09



EQUIPMENT DETAIL  
RCIM WITH DL KEYBOARD CONSOLE

C75-02


REV. DATE: 08/28/09



TYPICAL CONTROL ROOM  
INNOVA SINGLE PLANE

B5050C


REV. DATE: 08/26/08



EQUIPMENT DETAIL  
WORKSTATION

M1013AW

REV. DATE: 04/25/01



SHEET TITLE: EQUIPMENT DETAILS  
MODALITY TYPE: INNOVA 2100/3100/4100

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL LOCAL, STATE AND FEDERAL REGULATIONS AND STANDARDS. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:  
HYBRID OR  
TYPICAL FINAL LAYOUT

PROJECT	REVISION
5-94f	00
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST

REVISION HISTORY:

SHEET  
D2