

Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS	C1
EQUIPMENT LAYOUT (Equipment locations, heat loads, component weights, environmental specs)	A1
STRUCTURAL LAYOUT (Structural support/mounting locations for floor/wall/ceiling, wall support elevations)	S1
STRUCTURAL DETAILS (Floor and Ceiling loading information)	S2
ELECTRICAL LAYOUT (Contractor supplied wiring, interconnect methods, junction point locations and descriptions)	E1
ELECTRICAL SPECIFICATIONS (Maximum wiring run lengths, interconnect diagram, system power specifications)	E2
ELECTRICAL DETAILS	E3 THRU E4
EQUIPMENT DETAILS	D1 THRU D4

These drawings indicate the placement and interconnection of the listed equipment components. These drawings are not construction or site preparation drawings. Customer remains ultimately responsible for preparing the site to accommodate the operation of such equipment in compliance with GE Healthcare's written specifications and all applicable federal, state, and/or local requirements.

* REQUIRED REFERENCE *

Innova IGS
Pre Installation Manual
5421046-1-1EN

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the preIS manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

www.gehealthcare.com/siteplanning

GE Healthcare



Interventional Site Planning

CUSTOMER ACCEPTANCE



imagination at work

Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

GE Equipment Delivery Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist Rev 19					
Before using this document ensure you have the latest Rev from MyWorkshop on D0C0422752					
GEHC Global Order #:		Customer:			
GEHC PMI:		FE / Installer:			
The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.					
Inspection Date:	Storage is ready?	PHI is ready?	FE is ready?	Comments	if "N", enter comments or action plan
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					

GE Healthcare
Healthcare Project Implementation - Design Center
Milwaukee, Wisconsin
Copyright 2009 General Electric Company - Proprietary to GE

SHEET TITLE: SITE READINESS
MODALITY TYPE: INNOVA IGS 520, 530, 540
THIS PLAN IS SUBMITTED TO ASSIST IN THE LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN. EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL LOCAL, STATE, AND FEDERAL REGULATIONS AND STANDARDS. THE COMPANY CANNOT ACCEPT LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY (EP)
FINAL INSTALLATION DRAWINGS

PROJECT: REVISION
5-110F 00
DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST

REVISION HISTORY:

SHEET
C1

PIM R2
RQ - 140199

GE EQUIPMENT LISTING


ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	STRC PLAN	ELEC PLAN
1	1		UPS INTERFACE BOX			E45021B	---	UIB
2	1		ATLAS CABINET (C2)	659 lbs	1825 btu	B0558C	S100	C2
3	1		ATLAS CABINET (C1)	1115 lbs	3389 btu	B0558C	S100	C1
4	1		UPS CABINET	1170 lbs	4061 btu	E45025C	---	UPS
5	1		DETECTOR CHILLER	33 lbs	706 btu	B5049F	---	DC
6	1		COOLIX 4100 AUTOTRANSFORMER	66 lbs	238 btu	B-1G505	---	AT
7	1		COOLIX 4100 WATER CHILLER	264 lbs	18730 btu	B-1G503	---	CHLR
8	2		19 IN. MONITOR ON WALL SUPPORT	26 lbs	204 btu	C7619V	---	WBM2
9	1		LARGE DISPLAY MONITOR ON SINGLE MONITOR SUSPENSION	784 lbs	1706 btu	B2004 B2015	---	LDM WBM1
10	1		LARGE DISPLAY MONITOR CABINET	253 lbs	3412 btu	B2014	---	LDC
11	1		CONTROL ROOM MONITOR WITH DL KEYPAD	22 lbs	204 btu	C7412H C7617	---	S
12	1		OPERATORS CONSOLE	22 lbs	546 btu	C7617 C7502 B5050C	---	WBC1
13	1		AW WORKSTATION	81 lbs	1201 btu	C7617 M1013AN	---	C
14	1		INNOVA IQ TABLE	1750 lbs	614 btu		B5049N	LU5
15	1		INNOVA POSITIONER (REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION)	1653 lbs	2416 btu	B5050A B5050B B5050C B5050D B5050E B5050F B5050G B5050H B5050J B5050K B5050R	---	LC1
16	1		WORKSTATION CART				---	---
17	1		COMBOLAB CONSOLE INCLUDES MONITORS AND KEYBOARD	566 lbs	2935 btu		---	PC
18	1		TRAM NET RACK	8 lbs		B5047	---	TRAM
19	1		CLAB 2 PLUS AMPLIFIER	24 lbs	204 btu	B5051	---	AMP
20	1		REMOTE MONITORING WORKSTATION WITH TWO LCD MONITORS	81 lbs	1109 btu		---	RMOT
21	1		NURSING NOTES WORKSTATION	46 lbs	682 btu		---	S
22	1		MAPPING WORKSTATION	46 lbs	682 btu		---	S
23	1		CLINICAL WORKSTATION	46 lbs	682 btu		---	S
24	1		COLOR PRINTER		1054 btu		---	S
25	1		BOLUS CHASE HANDSWITCH	2 lbs			---	WBC
26	1		COUNTERBALANCED EYE AND THYROID SHIELD WITH R96 LAMP	143 lbs		B5031E	B5031F	LMP
27	2		LONGITUDINAL STATIONARY RAIL FOR XT SUSPENSION	68 lbs			B2007B	C
28	1		REMOTE CONTROL FOR INJECTOR	4 lbs		B5028		IEC
29	1		INJECTOR HEAD ON TABLE RAIL	15 lbs		B5030A		IH
30	1		INJECTOR ELECTRONICS	37 lbs	320 btu	B5028		IE
31	1		XR BUZZER (LOCATED ABOVE CEILING)	2 lbs		B5150H		XRB
32	1		VIVID J ULTRASOUND UNIT (ON SAFEDOCK CART)			B0572		---
33	1		3 KVA UPS CABINET (LARGE DISPLAY SUBSYSTEM OPTION)	66 lbs	546 btu	B2016	---	UPS3

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

50	1		INNOVA MAIN DISCONNECT, REFERENCE JUNCTION POINT "PDB" ON SHEET E1 FOR DETAILED DESCRIPTION.	326 lbs	1532 btu	E4502M	---	PDB
51	6		VITALING SPEAKER					---
52	1		VITALING CONSOLE			B0566		---
53	2		SKYTRON LIGHTING UNIT	50 lbs	341 btu	B2063	B2063A	SL
54	1		VITALING MICROPHONE					---
55	1		VITALING MICROPHONE (ONE ON MONITOR BRIDGE IN EXAM ROOM)					---

EQUIPMENT LAYOUT RECOMMENDED CEILING HEIGHT = 9'-6"

This equipment layout indicates the placement and interconnection of the indicated equipment components. There may be federal, state, and/or local requirements that could impact the placement of these components. It remains the Customer's responsibility for ensuring the site and final equipment placement complies with all applicable federal, state, and/or local requirements.



ANCILLARY ITEMS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
60	X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL 800-800-9760 GE CAT. NO. WX1ABW-DF-XIU
61	150-AMP LOCAL SERVICE DISCONNECT FOR LOCK-OUT/TAG-OUT CAPABILITY. (MAY BE A FUSED DISCONNECT, CIRCUIT BREAKER OR SAFETY SWITCH.)
62	MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W X 82 IN. H (118mm X 2108mm), CONTINGENT ON A 96 IN. (2438mm) CORRIDOR WIDTH
63	CATHETER CABINETS
64	MED GASES IN CEILING
65	BEARING BLOCK OUTLINE, SEE S1 FOR MORE INFORMATION.
66	SHELF - CUSTOMER TO PROVIDE ADEQUATE WALL SUPPORT
67	CABLE DRAPE RAIL
68	CUSTOMER SUPPLIED STORAGE CABINET
69	LEAD APRON RACK
70	COUNTER TOP FOR EQUIPMENT - MINIMUM DEPTH 30 IN. OR ADDITIONAL SHELVING MAY BE REQUIRED. PROVIDE GROMMETTED OPENINGS AS REQUIRED TO ROUTE INTERCONNECT CABLES TO RACEWAY BELOW COUNTERTOP.
71	MOBILE RADIATION SHIELD

THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.

90	X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL REFERENCE JUNCTION POINT "XRLC" ON SHEET E1 FOR DETAILED DESCRIPTION -CAT. NO. E4502SS FOR WARNING LIGHT & ROOM LIGHT CONTROL.
----	---

GENERAL SPECIFICATIONS

- THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.
- CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.
- RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.
- THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PREDICATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC., MAKE THIS ESSENTIAL FOR A PROPER IS. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC..
- ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM

SITE ENVIRONMENT SPECIFICATIONS

- AMBIENT OPERATING TEMPERATURE: EQUIPMENT ROOM WITH FLUORO UPS OPTION 68° TO 77° F. (20° TO 25° C)
- AMBIENT OPERATING TEMPERATURE: CONTROL ROOM 68° TO 77° F. (20° TO 25° C)
- AMBIENT OPERATING TEMPERATURE: EXAM ROOM-DESIGN FOR PATIENT/OPERATOR COMFORT TARGET TEMPERATURE 64° F (18° C)
- HUMIDITY: 30° TO 75° FOR EQUIPMENT AND CONTROL ROOMS AND 30° TO 70° FOR EXAM ROOM
- ALTITUDE: NOT TO EXCEED 9,842 FT. (3000M) ABOVE SEA LEVEL.
- THE ENVIRONMENT FOR THE ELECTRONICS CABINET MUST BE CONTROLLED SO THE ABOVE RESTRICTIONS ARE NOT EXCEEDED.
- DO NOT RESTRICT THE AIR INTAKE OR AIR EXHAUST OF THE SYSTEM COMPONENTS.
- ENVIRONMENTAL CONDITIONS LISTED ABOVE MUST BE MAINTAINED AT ALL TIMES INCLUDING FOR EXAMPLE OVERNIGHT, WEEKENDS, AND HOLIDAYS.

MAGNETIC INTERFERENCE SPECIFICATIONS

- IMAGE INTENSIFIERS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 1 GAUSS TO GUARANTEE SPECIFIED IMAGING PERFORMANCE.
- X-RAY TUBES MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE SPECIFIED PERFORMANCE.
- SYSTEM ELECTRONICS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE DATA INTEGRITY.
- OPERATORS CONSOLE EQUIPMENT MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO OBTAIN SPECIFIED GEOMETRIC LINEARITY.

GE Healthcare
Healthcare Project Implementation - Design Center
Milwaukee, WI

SHEET TITLE: **EQUIPMENT LAYOUT**
MODALITY TYPE: **INNOVA IGS 520, 530, 540**
THIS PLAN IS SUBMITTED TO VERIFY LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS AND DIMENSIONS SHOWN ON DRAWINGS AND TO THE COMPANY'S STANDARD PRACTICES AND PROCEDURES. GEHC ACCEPTS NO LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: **ROOM NO. EP ELECTROPHYSIOLOGY (EP) FINAL INSTALLATION DRAWINGS**




PROJECT	REVISION
5-110F	00

DATE: **18.Dec.13**
DRAWN BY: **JPH**
CHECKED BY: **TST**

REVISION HISTORY:

SHEET **A1**


TYPICAL WALL SUPPORT ELEVATIONS



SCALE: 1/4" = 1'-0"

STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-6"



STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S115, FOR UPS INTERFACE BOX.
2	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S100, FOR ATLAS CABINET.
3	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S107, FOR MAIN DISCONNECT CONTROL.
4	AREA OCCUPIED BY GE SUPPLIED TABLE BASEPLATE
5	AREA OCCUPIED BY GE SUPPLIED POSITIONER BASEPLATE
6	HATCHED AREA INDICATES MONITOR BRIDGE BEARING BLOCK PATH.
7	UNISTRUT OR EQUIVALENT SUPPORTS FOR FASTENING THE OVERHEAD COUNTERPOISED SUSPENSION. SUPPORTS TO BE LOCATED AS SHOWN. SUPPORT SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE, FLUSH WITH FINISHED CEILING. SUSPENSION REQUIRES 100 LBS/BOLT SUPPORT. METHODS OF SUPPORT THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION.
8	UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING CEILING SUPPORTED EQUIPMENT. SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE, FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 2'-2" AND REQUIRE 350 LBS. (897 LBS. IN SEISMIC REGIONS) PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION.
9	>>COMPONENTS FLUSH WITH CEILING<< UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING CABLE DRAPE RAIL SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE, FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 2'-2" AND REQUIRE 350 LBS. PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION. TO ORDER, CALL UNISTRUT WISCONSIN AT 888-796-8710.
10	MOUNT XR BUZZER BRACKET ON WALL.
11	SKYTRON LAMP LOCATION. LIGHTS MAY NEED TO BE RECESSED TO AVOID INTERFERENCE WITH OTHER CEILING MOUNTED ITEMS. CONSULT MANUFACTURER TO VERIFY LOCATION AND MOUNTING.

STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED TUBE HANGER OR OTHER EQUIPMENT ARE TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS. THE UNISTRUT OR EQUIVALENT STRUCTURE SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL, SQUARE AND IN THE SAME HORIZONTAL PLANE FLUSH WITH FINISHED CEILING. THE SYSTEM IS TO BE CROSS BRACED VERTICALLY, HORIZONTALLY AND DIAGONALLY TO ALLOW NO MOVEMENT AND A MAXIMUM OF 1.58mm (1/16") DEFLECTION.
(10) 12.7mm (1/2") DIA. x 38.1mm (1 1/2") LONG BOLTS WITH UNISTRUT 12.7mm (1/2") NUTS WITH SPRINGS ARE TO BE PROVIDED BY CUSTOMER OR HIS CONTRACTORS FOR EACH STATIONARY AND AUXILIARY SUPPORT RAIL. CLOSURE STRIPS SHALL BE PROVIDED FOR AREAS OF UNISTRUT EXPOSED AND WITHOUT MOUNTING UNITS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 6.35mm (1/4") BELOW THE FINISHED CEILING.
- CONTROL WALLS WITH TUBE HANGER PASSAGE ABOVE SHALL BE CONSTRUCTED TO 2130mm (7'-0") HIGH.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 3.17mm (1/8") IN 3050mm (10'-0")
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.
- IT IS THE CUSTOMER'S RESPONSIBILITY TO PERFORM ANY FLOOR OR WALL PENETRATIONS THAT MAY BE REQUIRED. THE CUSTOMER IS ALSO RESPONSIBLE FOR ENSURING THAT NO SUBSURFACE UTILITIES (E.G., ELECTRICAL OR ANY OTHER FORM OF WIRING, CONDUITS, PIPING, DUCT WORK OR STRUCTURAL SUPPORTS (I.E. POST TENSION CABLES OR REBAR)) WILL INTERFERE OR COME IN CONTACT WITH SUBSURFACE PENETRATION OPERATIONS (E.G. DRILLING AND INSTALLATION OF ANCHORS/SCREWS) PERFORMED DURING THE INSTALLATION PROCESS. TO ENSURE WORKER SAFETY, GE INSTALLERS WILL PERFORM SURFACE PENETRATION OPERATIONS ONLY AFTER THE CUSTOMER'S VALIDATION AND COMPLETION OF THE "GE SURFACE PENETRATION PERMIT"

SHEET TITLE: STRUCTURAL LAYOUT

MODALITY TYPE: INNOVA IGS 520, 530, 540

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY (EP)
FINAL INSTALLATION DRAWINGS





PROJECT	REVISION
5-110F	00

DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST

REVISION HISTORY:




SHEET
S1

GE Healthcare
Healthcare Project Implementation - Design Center
Manufacturer



FLOOR MOUNTING : INNOVA 2100-3100-4100 (UNITY)/OMEGA V LONG TABLE (WITH IQ TILT TABLE BASEPLATE) INSTALLATION (TEMPLATE NO. 2360133)


B5049N
REV. DATE: 06/04/09



WARNING!! THE RELATIONSHIP BETWEEN THE TABLE BASE AND THE POSITIONER BASEPLATE IS CRITICAL.

PRIOR TO DRILLING MOUNTING HOLES CONTACT LOCAL GE HEALTHCARE INSTALLATION PROJECT MANAGER OR LEAD FIELD ENGINEER TO VERIFY THAT THE PROPER FULL SIZE FLOOR MOUNTING TEMPLATE IS USED.

MEDICAL GAS FLOOR EXIT LOCATIONS



Customer/Contractor Alert: It is the responsibility of the Customer or their Contractor to drill all anchor/thru-bolting holes for anchoring the positioner and table to the floor. Refer to GEHC document no. *2290880-2-100 for installation preparation and procedures.

NOTE: THRU BOLTING IS HIGHLY PREFERRED FOR THE INSTALLATION OF THE POSITIONER BASEPLATE AND OMEGA TABLE. HARDENED BOLTS AND 4" x 4" [102mm x 102mm] STEEL PLATES TO BE USED ARE SUPPLIED BY GE HEALTHCARE AS INDICATED ON THE ACTUAL DETAIL DRAWING. BE ADVISED, HOWEVER, THAT ADDITIONAL SUPPORT STRUCTURES: STEEL BEAMS, PLATES, CORE BORING OF MOUNTING HOLES, ETC., ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

NOTE: IF THRU BOLTING IS NOT POSSIBLE, FLOOR ANCHORS CAN BE USED IF APPROVED BY CUSTOMERS STRUCTURAL ENGINEER. FOR ON GRADE INSTALLATIONS, MOUNTING KIT CAT. NO. 2286398 SHOULD BE ORDERED. ANCHORS INCLUDED IN KIT SHOULD BE APPROVED BY CUSTOMERS STRUCTURAL ENGINEER.

NOTE: BASEPLATES MUST BE LEVEL WITHIN 1/32" [0.79mm]

NOTE: JOISTS MUST BE SPANNED WITH STEEL REINFORCING. SIZE AND THICKNESS OF STEEL REINFORCING ARE DETERMINED BY THE ACTUAL PAN CONSTRUCTION ON SITE. STEEL PLATES, CHANNELS OR BEAMS MAY BE USED.

NOTE: DETERMINE THE POSITION OF THE "REBARS" IN THE CONCRETE FLOOR SO ANCHOR HOLES WILL NOT RUN INTO THEM.

* DOCUMENT FURNISHED BY GEHC AS PART OF INSTALLATION KIT

POSITIONER BOLT FORCES FOR WORST CASE CONDITIONS		OMEGA TABLE BOLT FORCES FOR WORST CASE CONDITIONS	
LOADS		LOADS	
HORIZONTAL ACCELERATION = 625 lbs. [284 Kg]	BOLT TENSION (AT BOLT "A") MAXIMUM TENSION = 881 lbs. [400 Kg]	BOLT TENSION MAXIMUM TENSION = 1938 lbs. [880 Kg]/BOLT	BOLT SHEAR MAXIMUM SHEAR = 407 lbs. [185 Kg]/BOLT
VERTICAL ACCELERATION = 209 lbs. [95 Kg]	BOLT SHEAR (U-ARM LOCKED) MAXIMUM SHEAR = 129 lbs. [54 Kg]/BOLT		

GE Healthcare
Healthcare Project Implementation - Design Center
Milwaukee, Wisconsin

SHEET TITLE: STRUCTURAL DETAILS
MODALITY TYPE: INNOVA IGS 520, 530, 540
THIS PLAN IS SUBMITTED TO SUPPORT LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPLIANCE, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE CODES AND REGULATIONS. GE HEALTHCARE ACCEPTS NO LIABILITY FOR ANY DAMAGES RESULTING FROM THE USE OF THIS PLAN OR FOR ANY DAMAGES RESULTING FROM THE USE OF THIS PLAN OR FOR ANY DAMAGES RESULTING FROM THE USE OF THIS PLAN.

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY (EP)
FINAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-110F	00
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST

REVISION HISTORY:






SHEET
S2

SCALE: 1/4" = 1'-0"

ELECTRICAL PLAN


RECOMMENDED CEILING HEIGHT = 9'-6"

JUNCTION POINT DESCRIPTIONS




A COMPLETE REVIEW OF ELECTRICAL OPTIONS MUST BE DISCUSSED WITH YOUR GE PROJECT MANAGER OF INSTALLATIONS, BEFORE BIDDING BEGINS.

CONDUIT RUNS: INNOVA IGS 520, 530, 540			
CONDUITS REQUIRED FROM POINT "XRLC" (CONDUITS ARE LOCATED ABOVE CEILING)			
(4)	XRLC	TO XRL1	ONE 1/2" CND.
(6)	XRLC	TO C2	ONE 1/2" CND.
(7)	XRLC	TO 120-V 1 ϕ POWER	CND. AS REQ'D
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUITS REQUIRED FROM POINT "WBBC" (CONDUITS ARE LOCATED BELOW FLOOR)			
(8)	WBBC	TO LUS	ONE 2 1/2" CND.
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUITS REQUIRED FROM POINT "XRB" (CONDUITS ARE LOCATED ABOVE CEILING)			
(9)	XRB	TO POWER STRIP IN CONTROL AREA	ONE 3/4" CND.
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUITS REQUIRED FROM POINT "LMP" (CONDUITS ARE LOCATED ABOVE CEILING)			
(11)	LMP	TO 120-V 1 ϕ POWER	CND. AS REQ'D
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUITS REQUIRED FROM POINT "WBM1" (CONDUITS ARE LOCATED ABOVE CEILING)			
(12)	WBM1	TO C1	TWO 2 1/2" CNDs. (UP TO FOUR FOR SUSPENSION) USABLE CABLE LENGTH UP TO 40 FT.
(13)	WBM1	TO WBC1	ONE 2 1/2" CND. USABLE CABLE LENGTH 40 FT.
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUITS REQUIRED FROM POINT "WC" (CONDUIT IS LOCATED IN OR BELOW FLOOR)			
(14)	WC	TO LC1	ONE EMPTY 3" CND. (FOR WATER LINES) (SEE WATER LINES) (SEE WATER MAIN TO REMOVE FITTING AT LC1) (SEE WATER MAIN TO REMOVE FITTING AT LC1) USABLE CABLE LENGTH UP TO 50 FT.
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUITS REQUIRED FOR "SKYTRON LAMPS" (CONDUITS ARE LOCATED ABOVE CEILING)			
(24)	SL	TO SL	ONE 3/4" CND.
(25)	SL	TO SWC	ONE 3/4" CND.
(26)	SWC	TO 120-V 1 ϕ POWER	CND. AS REQ'D
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUITS REQ'D. FOR "LARGE DISPLAY MONITOR" (CONDUITS ARE LOCATED ABOVE CEILING)			
(34)	LDM	TO LDC	ONE 3" & ONE 3/4" CND. CABLE LENGTH 100 FT.
(36)	LDC	TO WBC1	ONE 3" CND.
(37)	LDC	TO TRAM	ONE 3" CND. (RUN TO FLOOR BOX FOR TRAM)
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUIT RUNS: PHYSIO MONITORING/ IVUS			
CONDUITS REQUIRED FOR COMBO LAB			
(52)	PC	TO WBM1	ONE 3" CND. (LOCATED ABOVE CEILING)
(53)	TRAM/AMP	TO PC	ONE 3" CND. (LOCATED IN/BELOW FLOOR)
(55)	PC	TO RMO1	ONE 3" CND. (LOCATED IN/BELOW FLOOR)
CONDUITS REQUIRED FROM POINT "PDB" (CONDUITS ARE LOCATED ABOVE CEILING)			
*CABLES THAT COULD RUN IN DUCT SYSTEM INSTEAD			
(15)	PDB	TO UPS	*TWO CNDs. AS REQ'D.
(16)	PDB	TO UIB	*ONE 1 1/2" CND. USABLE CABLE LENGTH 70 FT.
(17)	PDB	TO RDS1	ONE 1/2" CND.
(18)	PDB	TO RDS2	ONE 1/2" CND.
(19)	PDB	TO C1	*TWO 2 1/2" CNDs. FOR FOUR CUSTOMER SUPPLIED POWER/ GROUND RUNS (AND GE SUPPLIED WIRES) CABLE LENGTH 19 FT.
(20)	PDB	TO C1	*ONE 1" CND. FOR TWO GE SUPPLIED SIGNAL CABLES CABLE LENGTH 19 FT.
(21)	PDB	TO PDB1	ONE CND. AS REQ'D.
(22)	PDB1	TO 480-V 3 ϕ POWER	CND. AS REQ'D.
(23)	PDB	TO IE	(INJECTOR POWER) CONSULT MFG. (RUN IN DUCT/ CONDUIT SYSTEM)
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			
CONDUITS REQUIRED FOR AN "INJECTOR" (CONDUITS ABOVE CEILING OR BELOW FLOOR)			
(27)	IE	TO IH	ONE 2 1/2" CND.
(28)	IE	TO IEC	ONE 2 1/2" CND.
NOTE: SEE E2 PAGE FOR MAXIMUM RUN LENGTHS			



JUNCTION POINT NOTES


- ALL JUNCTION BOXES, CONDUIT, DUCT, DUCT DIVIDERS, SWITCHES, CIRCUIT BREAKERS, ETC., ARE TO BE SUPPLIED AND INSTALLED BY CUSTOMER'S ELECTRICAL CONTRACTOR.
- CONDUIT AND DUCT RUNS SHALL HAVE SWEEP RADIUS BENDS
- CONDUITS AND DUCT ABOVE CEILING OR BELOW FINISHED FLOOR MUST BE INSTALLED AS NEAR TO CEILING OR FLOOR AS POSSIBLE TO REDUCE RUN LENGTH.
- CEILING MOUNTED JUNCTION BOXES ILLUSTRATED ON THIS PLAN MUST BE INSTALLED FLUSH WITH FINISHED CEILING.
- ALL DUCTWORK MUST MEET THE FOLLOWING REQUIREMENTS:
 - DUCTWORK SHALL BE METAL WITH DIVIDERS AND HAVE REMOVABLE, ACCESSIBLE COVERS.
 - DUCTWORK SHALL BE CERTIFIED/RATED FOR ELECTRICAL PURPOSES.
 - DUCTWORK SHALL BE ELECTRICALLY AND MECHANICALLY BONDED TOGETHER IN AN APPROVED MANNER.
 - PVC AS A SUBSTITUTE MUST BE USED IN ACCORDANCE WITH ALL LOCAL AND NATIONAL CODES.
- ALL OPENINGS IN ACCESS FLOORING ARE TO BE CUT OUT AND FINISHED OFF WITH GROMMET MATERIAL BY THE CUSTOMER'S CONTRACTOR.
- GENERAL CONTRACTOR TO INSERT PULL CORDS FOR ALL CABLE RUN CONDUITS BETWEEN THE EQUIPMENT ROOM AND THE OPERATORS CONTROL ROOM.
- 10 FOOT PIGTAILS AT ALL JUNCTION POINTS.
- ALL WIRING MUST BE THIN OR TFFN STRANDED COPPER THERMOPLASTIC 600 VOLT OR EQUIVALENT INSULATION. **ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.**
- GROUNDING IS CRITICAL TO EQUIPMENT FUNCTION AND PATIENT SAFETY. SITE MUST CONFORM TO WIRING SPECIFICATIONS SHOWN ON THIS PLAN.



FEEDER TABLE REV. DATE: 12/22/10

* CALCULATIONS BASED UPON NOMINAL VOLTAGE, WIRE SIZE IN AWG.
 * RECOMMENDED FEEDER SIZES FROM DIST. TRANS. TO ROOM DISCONNECT. CALCULATIONS ARE AT NOMINAL VOLTAGE BASED UPON 1/0 WIRE SIZE FROM ROOM DISCONNECT TO POWER CABINET WITH A MAXIMUM RUN OF 25 FT.
 * NEUTRAL MUST BE TERMINATED INSIDE THE MAIN DISCONNECT PANEL, AND NOT AT ANY GE CABINET.
 * THE GROUNDING CONDUCTOR () WILL BE A 2 AWG MINIMUM OR MEET LOCAL CODE REQUIREMENTS, WHICHEVER IS LARGER. THIS GROUND WILL RUN FROM THE EQUIPMENT BACK TO THE POWER SOURCE/MAIN GROUNDING POINT AND ALWAYS TRAVEL IN THE SAME CONDUIT WITH THE FEEDERS AND NEUTRAL.
 * MINIMUM WIRE SIZE FOR CIRCUIT BREAKER, BASED ON RECOMMENDED OVERCURRENT PROTECTION.
 * FOR A FULL SYSTEM UPS, REFER TO ELECTRICAL DETAILS FOR UPS FEEDER WIRES.
 * IF THE FEEDER IS BIGGER THAN 3/0, THE HOSPITAL MUST PROVIDE AND INSTALL A REDUCTION BOX

RUN LENGTH IN FEET	POWER SUPPLY VOLTAGE							
	324-396 360	342-418 380	360-440 400	378-462 420	396-488 440	414-506 460	432-528 480	
50	FEDER/GRND	FEDER/GRND	FEDER/GRND	FEDER/GRND	FEDER/GRND	FEDER/GRND	FEDER/GRND	
100	1/0 (2)	1/0 (2)	1/0 (2)	1/0 (2)	1/0 (2)	1/0 (2)	1/0 (2)	
150	1/0 (2)	1/0 (2)	1/0 (2)	1/0 (2)	1/0 (2)	1/0 (2)	1/0 (2)	
200	3/0 (2)	2/0 (2)	2/0 (2)	2/0 (2)	2/0 (2)	2/0 (2)	2/0 (2)	
250	4/0 (2)	4/0 (2)	3/0 (2)	3/0 (2)	3/0 (2)	3/0 (2)	3/0 (2)	
300	300M (2)	300M (2)	250M (2)	4/0 (2)	3/0 (2)	3/0 (2)	3/0 (2)	
350	400M (2)	350M (2)	300M (2)	250M (2)	4/0 (2)	4/0 (2)	4/0 (2)	
400	600M (2)	500M (2)	400M (2)	300M (2)	300M (2)	250M (2)	4/0 (2)	
450	700M (2)	600M (2)	500M (2)	400M (2)	350M (2)	300M (2)	300M (2)	



POINT	DESCRIPTION	QTY.	HARDWARE	DETAIL NO., SHT. E3
AMP	AMPLIFIER	1	AMPLIFIER CABLES EXIT AT 'TRAM'	ELEC-13
AT	COOLIX 4100 AUTO TRANSFORMER	1	EXTERNALLY CONNECTED TO TRAM CHILLER	ELEC-3
C1	ATLAS CABINET	2	3 1/2 IN. DIA. CHASE NIPPLE BUSHING	ELEC-3
C2	ATLAS CABINET	2	3 1/2 IN. DIA. CHASE NIPPLE BUSHING	ELEC-3
CHLR	COOLIX 4100 WATER CHILLER	1	12 IN. OF GROMMET MATERIAL FOR 1 1/2 X 3 IN. OPENING IN DUCT COVER	ELEC-5
DC	DETECTOR CHILLER	1	2 1/2 IN. DIA. CHASE NIPPLE	ELEC-25
IE	INJECTOR ELECTRONICS	1	2 1/2 IN. DIA. BUSHING & LOCKNUT	ELEC-25
IEC	INJECTOR CONTROL	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
IH	INJECTOR HEAD	1	EXTERNALLY CONNECTED AT TABLE BASE	ELEC-6
LC1	INNOVA LC	1	FLUSH FLOOR FITTING	ELEC-101
LDC	LARGE DISPLAY CABINET	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-6
LDM	LARGE DISPLAY MONITOR	1	COVERPLATE 1 1/2 X 3 IN. CHASE NIPPLE	ELEC-8
LMP	SURGICAL LAMP	1	COVERPLATE 1 1/2 X 3 IN. CHASE NIPPLE	ELEC-8
LUS	OMEGA TABLE	1	3 1/2 IN. DIA. CHASE NIPPLE	ELEC-25
PC	WORKSTATION	1	COVERPLATE 1 1/2 X 3 IN. CHASE NIPPLE	ELEC-13
PDB	MAIN DISCONNECT	1	150-AMP PANEL INCLUDED IN ORDER	ELEC-161
PDB1	LOCAL SERVICE DISCONNECT	1	150-AMP LOCAL SERVICE DISCONNECT (CUSTOMER SUPPLIED)	ELEC-16
RDS1	EMERGENCY OFF	1	PROVIDE A SINGLE GANG, 2 1/8 IN. DEEP, FLUSH MTD. WALL BOX.	ELEC-16
RDS2	EMERGENCY OFF	1	PROVIDE A SINGLE GANG, 2 1/8 IN. DEEP, FLUSH MTD. WALL BOX.	ELEC-16
RMOT	REMOTE WORKSTATION	1	COVERPLATE 1 1/2 X 3 IN. CHASE NIPPLE	ELEC-13
SL	SKYTRON LAMP	2	3/4 IN. DIA. BUSHING & LOCKNUT	ELEC-124
SWC	SKYTRON WALL CONTROL	1	3 GANG BOX, 2 3/4 IN. DEEP	ELEC-124
TRAM	TRAM NET	1	COVERPLATE 1 1/2 X 3 IN. CHASE NIPPLE	ELEC-13
UIB	UPS INTERFACE BOX	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
UPS	UPS CABINET	1	32 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
UPS3	3 KVA UPS (LD SUBSYSTEM)	1	EXTERNALLY CONNECTED TO LARGE DISPLAY CABINET - 'LDC'	ELEC-5
WBBC	BOLUS WALLBOX	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
WBC1	OPERATORS CONSOLE	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
WBM1	TV MONITOR	2	SHARED CEILING BOX WITH 'LDM' 1 1/2 IN. DIA. CHASE NIPPLE	ELEC-8
WBM2	TV MONITOR	1	12 IN. OF GROMMET MATERIAL FOR A 3 X 3 IN. OPENING IN DUCT COVER	ELEC-5
WC	WATER CHILLER HOSE OUTLET	1	1 IN. IN. CONDUIT STUBBED 2 IN. ABOVE FLOOR	ELEC-9
XRB	XR BUZZER (LOCATED ABOVE CEILING)	1	COVERPLATE 1 1/2 X 3 IN. CHASE NIPPLE	ELEC-8
XRL1	WARNING LIGHT	1	COVERPLATE SINGLE GANG BOX 1 1/2 X 3 IN. CHASE NIPPLE	ELEC-157
XRLC	WARNING LIGHT CONTROL	1	E450SS WARNING LIGHT & DIMMER LIGHT CONTROL OR EQUIVALENT MAX 24V CONTROLLER	ELEC-157

CONTRACTOR SUPPLIED AND INSTALLED WIRING
 ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS.

WIRE RUN, FROM - TO	QUANTITY, WIRE SIZE/COLOR
<22> 3 PHASE > PDB1	3-BLACK, 1-WHITE, 1-GREEN (REFER TO FEEDER TABLE)
<21> PDB1 > PDB	3-BLACK, 1-WHITE, 1-GREEN (REFER TO FEEDER TABLE)
<19> PDB > C1 <JEDI>	3-1/0 BLACK, 1-1/0 GREEN
<19> PDB > C1 <PDU>	2-ND. 10 BLACK, 1-ND. 10 GREEN
<19> PDB > C2	3-ND. 8 BLACK, 1-ND. 8 GREEN
<19> PDB > AT	3-ND. 10 BLACK, 1-ND. 10 GREEN
<15> PDB > UPS	6-ND. 6 BLACK, 1-ND. 6 WHITE, 2-ND. 6 GREEN
<17> PDB > RDS1	2-ND. 14 BLACK, 2-ND. 14 WHITE, 1-ND. 14 GREEN
<18> PDB > RDS2	2-ND. 14 BLACK, 2-ND. 14 WHITE, 1-ND. 14 GREEN
<7> XRL > 1 PHASE	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
<6> XRLC > C2	2-ND. 14 BLACK, 2-ND. 14 WHITE, 1-ND. 14 GREEN
<4> XRL1 > XRLC	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
<11> 120-V > LMP	1-BLACK, 1-WHITE, 1-GREEN - (SIZE AS REQUIRED)
<26> 120V > SWC	1-ND. 14 BLACK, 1-ND. 14 WHITE, 1-ND. 14 GREEN
<25> SWC > SL	1-ND. 16 GREEN, 1-ND. 16 BLACK, 1-ND. 16 WHITE
<24> SWC > SL	1-ND. 16 GREEN, 1-ND. 16 BLACK, 1-ND. 16 WHITE

GE Healthcare
 Healthcare Project Implementation - Design Center
 Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL LAYOUT
 MODALITY TYPE: INNOVA IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO ASSIST IN THE DETERMINATION OF THE HEALTHCARE EQUIPMENT AND ASSOCIATED APPLIANCE, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE LOCAL, STATE AND FEDERAL REGULATIONS AND CODES. GE HEALTHCARE ACCEPTS NO LIABILITY FOR ANY DAMAGES RESULTING FROM THE ACTUAL CONSTRUCTION OF THE PROJECTS SHOWN AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY (EP)
 FINAL INSTALLATION DRAWINGS

PROJECT REVISION
 5-110F 00

DATE: 18.Dec.13
 DRAWN BY: JPH
 CHECKED BY: TST

REVISION HISTORY:


THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

PIM R2
 RQ - 140199

SHEET
E1


INTERCONNECT DIAGRAM

TYPICAL VIEWS



EQUIPMENT DESCRIPTIONS				OPTIONS					
ITEM	DESCRIPTION	WEIGHT (lb)	HEAT DISSIPATION (btu)	DRAWING DESIGNATOR	ITEM	DESCRIPTION	WEIGHT (lb)	HEAT DISSIPATION (btu)	DRAWING DESIGNATOR
1	- XR BUZZER	2		XRB	22	- ADVANTAGE WINDOWS WORKSTATION	81	1201	AW
2	- ATLAS CABINET C2	659	1825	C2	23	- IVUS VOLCANO CONSOLE	68	1631	IVUS
3	- ATLAS CABINET C1	1115	3389	C1	24	- IVUS VOLCANO COLOR PRINTER	X	X	
4	- DETECTOR CONDITIONER	33	706	DC	25	- INJECTOR HEAD	15		IH
5	- COOLIX 4100 WATER CHILLER	265	18725	CHLR	26	- INJECTOR ELECTRONICS	37	320	IE
6	- 20kva UPS CABINET	1170	4061	UPS	27	- REMOTE CONTROL FOR INJECTOR	4		IEC
7	- UPS INTERFACE BOX			UIB	28	- LAMP (RADIATION SHIELD TRACK)	143		LMP
8	- TV CEILING SUSPENSION (8 MONITOR)	557	1228	WBM1	29	- LARGE DISPLAY MONITOR CABINET	254	3412	LDC
9	- INNOVA LC POSITIONER	1653	2416	LC1	30	- LARGE DISPLAY MONITOR	784	1706	LDM
10	- OMEGA I/O TABLE	1750	614	IUS	31	- MACH 3 TRANSFORMER	70	X	M3T
11	- INNOVA VOIM FIFTH DL KEYBOARD CONSOLE	22	204		32	- MACLAB PHYSIO MONITORING	566	2935	PC
12	- VOIM OPERATOR CONSOLE	22	546		33	- PRINTER (PHYSIO)	X		
13	- ROOM LIGHTS			WBC1	34	- TRAM (PHYSIO)	8	X	TRAM
14	- XRAY WARNING LAMP			RML1	35	- REMOTE OPERATING TERMINAL (PHYSIO)	46	682	RMOT
15	- XRAY WARNING LAMP CONTROLLER			XRL1	36	- MICRO PACE STIMULATOR (PHYSIO)	X	X	MP
16	- RDS1 PUSHBUTTON			RDS1	37	- SKYTRON LIGHTING UNIT	50	341	SL
17	- RDS2 PUSHBUTTON			RDS2	38	- 150 kva UPS	2160	31802	UPS
18	- PDB MAIN DISCONNECT	326	1532	PDB	39	- UPS BATTERY CABINET	3529	X	
19	- LOTO DISCONNECT BREAKER			PDB1	40	- MAINTENANCE BYPASS PANEL	350	X	MBP
20	- 3kva UPS CABINET	81	546	UPS1	41	- COOLIX 4100 AUTOTRANSFORMER	99	239	AT
21	- BOLUS CHASE HANDSWITCH	2		WBBC					

FULL SIZE UPS (OPTIONAL)



REV DATE: 12.Mar.12

POWER SPECIFICATIONS

INNOVA SYSTEMS

REV. DATE: 01/04/07

VOLTAGE: PRIMARY SOURCE IS REQUIRED FOR ALL INSTALLATIONS. RANGE OF LINE VOLTAGES: NOMINAL LINE VOLTAGE OF 360 TO 480, 3 PHASE, 50 OR 60 HZ. REQUIRED POWER SUPPLY: WYE DISTRIBUTION. MAXIMUM DAILY VOLTAGE VARIATION MUST FALL WITHIN ONE OF THE RANGES IN TABLE A.

NOMINAL VOLTAGE	NORMAL RANGE ±10 PERCENT	CURRENT (AMPS)	
		MAX. MOMENTARY	CONTINUOUS
360	324-396	304	32
380	342-418	289	31
400	360-440	274	29
420	378-462	264	28
440	396-484	249	26
460	414-506	238	25
480	432-528	228	24

NOTE: LOW LINE CONDITIONS MAY INHIBIT SOME HIGH KVP TECHNIQUES. THE GENERATOR AUTOMATICALLY ESTABLISHES THESE INHIBITS BASED ON ACTUAL LINE CONDITIONS AND SYSTEM REGULATION. PHASE-BALANCE: PHASE-TO-PHASE VOLTAGES MUST BE WITHIN +2 PERCENT OF THE LOWEST PHASE-TO-PHASE VOLTAGE. MAXIMUM ALLOWABLE TRANSIENT VOLTAGE EXCURSIONS ARE 2.5 PERCENT OF RATED LINE VOLTAGE AT A MAXIMUM DURATION OF 5 CYCLES AND FREQUENCY OF 10 TIMES PER HOUR. POWER DEMAND: CONTINUOUS POWER DEMAND = 20KVA. (MAX DEMAND = 171 KVA)

DEMAND	ADVANTX 100
kva * POWER FACTOR AT	171 0.9
mA	1250
kvp	80

* DEMAND INCLUDES POWER FOR ENTIRE ADVANTX SYSTEM. LINE VOLTAGE REGULATION AT MAXIMUM POWER DEMAND MUST BE LESS THAN OR EQUAL TO 6 PERCENT. FOR A SINGLE UNIT INSTALLATION, THE MINIMUM TRANSFORMER SIZE IS 225 KVA.

ELECTRICAL NOTES

- NOTE 1: ALL WIRES SPECIFIED SHALL BE COPPER STRANDED, FLEXIBLE, THERMO-PLASTIC, COLOR CODED, CUT 10 FOOT LONG AT OUTLET BOXES, DUCT TERMINATION POINTS OR STUBBED CONDUIT ENDS. ALL CONDUCTORS, POWER, SIGNAL AND GROUND, MUST BE RUN IN A CONDUIT OR DUCT SYSTEM. ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS. WIRE RUNS MUST BE CONTINUOUS COPPER STRANDED AND FREE FROM SPLICES. ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.
- NOTE 2: WIRE SIZES GIVEN ARE FOR USE OF EQUIPMENT. LARGER SIZES MAY BE REQUIRED BY LOCAL CODES.
- NOTE 3: IT IS RECOMMENDED THAT ALL WIRES BE COLOR CODED, AS REQUIRED IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 4: CONDUIT SIZES SHALL BE VERIFIED BY THE ARCHITECT, ELECTRICAL ENGINEER OR CONTRACTOR, IN ACCORDANCE WITH LOCAL OR NATIONAL CODES.
- NOTE 5: CONVENIENCE OUTLETS ARE NOT ILLUSTRATED. THEIR NUMBER AND LOCATION ARE TO BE SPECIFIED BY OTHERS. LOCATE AT LEAST ONE CONVENIENCE OUTLET CLOSE TO THE SYSTEM CONTROL, THE POWER DISTRIBUTION UNIT AND ONE ON EACH WALL OF THE PROCEDURE ROOM. USE HOSPITAL APPROVED OUTLET OR EQUIVALENT.
- NOTE 6: GENERAL ROOM ILLUMINATION IS NOT ILLUSTRATED. CAUTION SHOULD BE TAKEN TO AVOID EXCESSIVE HEAT FROM OVERHEAD SPOTLIGHTS. DAMAGE CAN OCCUR TO CEILING MOUNTING COMPONENTS AND WIRING IF HIGH WATTAGE BULBS ARE USED. RECOMMEND LOW WATTAGE BULBS NO HIGHER THAN 75 WATTS AND USE DIMMER CONTROLS (EXCEPT MR). DO NOT MOUNT LIGHTS DIRECTLY ABOVE AREAS WHERE CEILING MOUNTED ACCESSORIES WILL BE PARKED.
- NOTE 7: ROUTING OF CABLE DUCTWORK, CONDUITS, ETC., MUST RUN DIRECT AS POSSIBLE OTHERWISE MAY RESULT IN THE NEED FOR GREATER THAN STANDARD CABLE LENGTHS (REFER TO THE INTERCONNECTION DIAGRAM FOR MAXIMUM USABLE LENGTHS POINT TO POINT).
- NOTE 8: CONDUIT TURNS TO HAVE LARGE, SWEEPING BENDS WITH MINIMUM RADIUS IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 9: A SPECIAL GROUNDING SYSTEM IS REQUIRED IN ALL PROCEDURE ROOMS BY SOME NATIONAL AND LOCAL CODES. IT IS RECOMMENDED IN AREAS WHERE PATIENTS MIGHT BE EXAMINED OR TREATED UNDER PRESENT, FUTURE, OR EMERGENCY CONDITIONS. CONSULT THE GOVERNING ELECTRICAL CODE AND CONFER WITH APPROPRIATE CUSTOMER ADMINISTRATIVE PERSONNEL TO DETERMINE THE AREAS REQUIRING THIS TYPE OF GROUNDING SYSTEM.
- NOTE 10: THE MAXIMUM POINT TO POINT DISTANCES ILLUSTRATED ON THIS DRAWING MUST NOT BE EXCEEDED.
- NOTE 11: PHYSICAL CONNECTION OF PRIMARY POWER TO GE EQUIPMENT IS TO BE MADE BY CUSTOMERS ELECTRICAL CONTRACTOR WITH THE SUPERVISION OF A GE REPRESENTATIVE. THE GE REPRESENTATIVE WOULD BE REQUIRED TO IDENTIFY THE PHYSICAL CONNECTION LOCATION, AND INSURE PROPER HANDLING OF GE EQUIPMENT.

---	CUSTOMER/CONTRACTOR SUPPLIED WIRING. ROUTE IN ADEQUATE CONDUIT OR RACEWAY.
---	GE FURNISHED CABLE RUNS. ROUTE IN EMPTY CONDUIT OR RACEWAY.
59' [18M]	MAXIMUM RUN LENGTH BETWEEN JUNCTION POINTS. Feet [Meters]

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

GE Healthcare
Healthcare Project Implementation - Design Center
Manufacturing

SHEET TITLE: ELECTRICAL SPECIFICATIONS
MODALITY TYPE: INNOVA IGS 520, 530, 540
THIS PLAN IS SUBMITTED TO ASSIST IN THE LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS, IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS AND SPECIFICATIONS OF THE EQUIPMENT MANUFACTURERS AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY (EP)
FINAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-110F	00

DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST

REVISION HISTORY:

SHEET
E2

PIM R2
RQ - 140199

ELECTRICAL DETAIL X-RAY WARNING LIGHT & ROOM LIGHT CONTROL PANEL

ELEC-157
REV. DATE: 04/23/09

FROM GE IMAGING SYSTEM ON SIGNAL IN "C2" CABINET
MAXIMUM 24-VAC

E4502SS
X-RAY ROOM WARNING LIGHT / ROOM LIGHTING CONTROL PANEL
X-RAY WARNING LIGHT OR ROOM LIGHT ARE NOT PART OF THIS CAT. NO.

120-VAC 20A MAXIMUM

X-RAY WARNING LIGHT

ROOM LIGHTS

THE E4502SS IS RECOMMENDED IF "X-RAY ON" WARNING LIGHT AND ROOM LIGHT CONTROL ARE UTILIZED

THE R4502RL IS RECOMMENDED IF "X-RAY ON" WARNING LIGHT ONLY

E4502RL
X-RAY ROOM WARNING LIGHT CONTROL PANEL
X-RAY WARNING LIGHT IS NOT PART OF THIS CAT. NO.

120-VAC 20A MAXIMUM

X-RAY WARNING LIGHT

CONTROL PANEL CAN BE LOCATED ABOVE THE CEILING NEAR THE WARNING LIGHT

UNLESS SPECIFIED ON SHEET A1 AS BEING INCLUDED ON EQUIPMENT ORDER, ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER'S CONTRACTOR

ELECTRICAL DETAIL HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5
REV. DATE: 03/19/04

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" (610mm)	2
18" (457mm)	2
10" (254mm)	2
6" (152mm)	1
4" (102mm)	1

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

DETAIL NOT TO SCALE

ELECTRICAL DETAIL VERTICAL WALL DUCT (TYPICAL)

ELEC-6
REV. DATE: 03/19/04

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" (610mm)	2
18" (457mm)	2
10" (254mm)	2
6" (152mm)	1
4" (102mm)	1

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

DETAIL NOT TO SCALE

ELECTRICAL DETAIL FLUSH FLOOR DUCT (TYPICAL)

ELEC-3
REV. DATE: 4/01/04

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" (610mm)	2
18" (457mm)	2
10" (254mm)	2
6" (152mm)	1
4" (102mm)	1

12" TRENCH DUCT: A = 10", B = 12"

18" TRENCH DUCT: A = 18", B = 20"

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

DETAIL NOT TO SCALE

ELECTRICAL DETAIL INSITE CONNECTION (TYPICAL)

ELEC-1
REV. DATE: 04/24/02

ONE OF THE FOLLOWING TWO SELECTIONS MUST BE INSTALLED AT THE LOCATION SHOWN ON THE ELECTRICAL PLAN (SHEET E1) FOR GE INSITE CONNECTION BASED UPON SYSTEM CONFIGURATION.

A) ONE INTERNET ACCESSIBLE VIRTUAL PRIVATE NETWORK (VPN) CONNECTION WITH A STATIC IP ADDRESS, AND ONE TELEPHONE LINE - DEDICATED-DIRECT-DIALING, VOICE GRADE.

OR

B) TWO TELEPHONE LINES - ONE DEDICATED DIRECT-DISTANCE-DIALING, VOICE GRADE AND ONE A DEDICATED DATA LINE.

FINISHED CEILING

SINGLE GANG J.B.

1" CONDUIT FROM J.B. TO ABOVE FINISHED CEILING.

TO BE DETERMINED

FINISHED FLOOR

COVERPLATE WITH TWO TELEPHONE RECEPTACLES OR ONE TELEPHONE RECEPTACLE AND ONE NETWORK RECEPTACLE

ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER OR THEIR CONTRACTOR.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL BOX WITH COVERPLATE AND NETWORK JACK

ELEC-83
REV. DATE: 10/06/98

BOX

NETWORK JACK

COVERPLATE

DETAIL NOT TO SCALE

ELECTRICAL DETAIL NETWORK CONNECTION (TYPICAL)

ELEC-84
REV. DATE: 03/06/04

FOR NUCLEAR SYSTEMS A DIRECT NETWORK CONNECTION IS TO BE MADE BETWEEN THE SYSTEM AND THE REVIEW WORKSTATION.

LOCAL AREA NETWORK

FINISHED CEILING

1/2" CONDUIT FROM J.B. TO ABOVE FINISHED CEILING.

TO BE DETERMINED

FINISHED FLOOR

COVERPLATE WITH NETWORK RECEPTACLE

DETAIL NOT TO SCALE

ELECTRICAL DETAIL INNOVA PLUS MAIN DISCONNECT PANEL

ELEC-161
REV. DATE: 09/27/10

26.78" [680.2mm]

SPACE FOR INCOMING HOSPITAL INPUT WIRES

TOP VIEW

FRONT VIEW

9.125" [231.78mm]

SPACE FOR INCOMING CONDUIT AREA BOTH SIDES

72.68" [1846.072mm]

72.72" [1847.09mm]

SIDE VIEW

8" [203.2mm]

SPACE FOR INCOMING UPS WIRES

NOTE: LINE WIRES AND UPS WIRES MUST NOT BE MIXED.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL FLUSH FLOOR DUCT (TYPICAL)

ELEC-25
REV. DATE: 4/01/04

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" (610mm)	2
18" (457mm)	2
10" (254mm)	2
6" (152mm)	1
4" (102mm)	1

12" TRENCH DUCT: A = 10", B = 12"

18" TRENCH DUCT: A = 18", B = 20"

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

DETAIL NOT TO SCALE

ELECTRICAL DETAIL CONDUITS THRU-FLOOR (TYPICAL)

ELEC-9
REV. DATE: 08/08/94

FINISHED FLOOR

HARDWARE

1.5" (38 mm) TYP.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5A
REV. DATE: 06/16/08

DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" (610mm)	2
18" (457mm)	2
10" (254mm)	2
6" (152mm)	1
4" (102mm)	1

2" X 4" OPENING CUT INTO TOP OF DUCT FOR 12" OF GROMMETED MATERIAL

REFER TO CHART FOR MINIMUM DIVIDER REQUIREMENT LOCAL CODES MAY REQUIRE ADDITIONAL DIVIDERS

DETAIL NOT TO SCALE

ELECTRICAL DETAIL BOX WITH COVERPLATE (TYPICAL)

ELEC-8
REV. DATE: 09/30/94


OUTLET BOX

HARDWARE

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
J.B. / WALL DUCT DETAIL (TYPICAL)


ELEC-2
REV. DATE: 09/30/94



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
POSITIONER INTERCONNECT DETAIL, FLUSH IN FLOOR

ELEC-101
REV. DATE: 01/07/04




NOTE: FLUSH FLOOR INTERFACE FITTING IS PART OF GE INSTALLATION KIT CAT. NO. B5079BC AND IS TO BE INSTALLED BY CUSTOMER or CUSTOMER'S CONTRACTOR

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
FLOOR BOX WITH NIPPLES (TYPICAL)


ELEC-13
REV. DATE: 09/30/94



DETAIL NOT TO SCALE


ELECTRICAL DETAIL
SKYTRON INTENSITY AND POSITIONING WALL CONTROL

ELEC-124
REV. DATE: 08/04/08



ELECTRICAL DETAIL
EMERGENCY OFF BUTTON


ELEC-16
REV. DATE: 05/14/09



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
SKYTRON LIGHTING UNIT

ELEC-98
REV. DATE: 08/04/08



DETAIL NOT TO SCALE

SHEET TITLE: ELECTRICAL DETAILS

MODALITY TYPE: INNOVA IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO SUGGEST LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE MANUFACTURER'S INSTALLATION INSTRUCTIONS AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY
(EP)
FINAL INSTALLATION DRAWINGS













PROJECT	REVISION
5-110F	00
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST

REVISION HISTORY:

SHEET
E4

PIM R2
RQ - 140199

GE Healthcare
Healthcare Project Implementation - Design Center
Manufacture, Wisconsin



GE Healthcare
Healthcare Project Implementation - Design Center
Manufacturer

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: INNOVA IGS 520, 530, 540
THIS PLAN IS SUBMITTED TO ASSIST IN THE SELECTION OF HEALTHCARE EQUIPMENT AND ASSOCIATED SUPPLIES. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE ACTUAL CONSTRUCTION PRACTICES AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY (EP)
FINAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-110F	00

DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST

REVISION HISTORY:

SHEET
D1

PIM R2
RQ - 140199

EQUIPMENT DETAIL
DLX or DL KEYPAD

C7412H
REV. DATE: 09/03/03

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
18" FLAT PANEL MONITOR

C76-17
REV. DATE: 08/28/09

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
RCIM WITH DL KEYBOARD CONSOLE

C75-02
REV. DATE: 10/25/10

DETAIL NOT TO SCALE

TYPICAL CONTROL ROOM
INNOVA SINGLE PLANE

B5050C
REV. DATE: 08/26/08

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
ADVANTAGE WINDOWS WORKSTATION

M1013AW
REV. DATE: 20.MAR.12

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA VASCULAR SYSTEM

B5050A
REV. DATE: 06/07/05

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA VASCULAR SYSTEM

B5050B
REV. DATE: 06/07/05

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INNOVA VASCULAR SYSTEM

B5050
REV. DATE: 02/22/05

SHIPPING DIMENSIONS:
110"L x 45.5"W x 77"H
(2790mmL x 1160mmW x 1950mmH)
(ON DOLLY)
WIDTH IS REDUCED TO 34" [865mm]
BY REMOVING SIDE RAILS

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
SHIPPING CREATE FOR OMEGA & INNOVA IQ TABLES

B5050E
REV. DATE: 12/07/09

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA C2 CABINET ON SHIPPING DOLLY

B5050F
REV. DATE: 12/07/09

SHIPPING WEIGHT: 821 LBS. (372 kg)(C2-HARMONY SYSTEMS)
SHIPPING WEIGHT: 785 LBS. (356 kg)(C2-BIPLANE SYSTEMS)

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
SHIPPING DOLLY FOR INNOVA LC POSITIONER

B5050G
REV. DATE: 12/07/09

NOTE:
BOTH ENDS OF THE DOLLY CAN BE REMOVED WHICH WILL SHORTEN LC GANTRY DOLLY DONE TO 86.22" (2190mm) RECOMMEND ONLY ONE SIDE BE REMOVED WHEN DELIVERY THROUGH HOSPITAL.

SHIPPING WEIGHT: 2340 lbs. (1060 kg)

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA C1 & C3 CABINETS ON SHIPPING DOLLY

B5050H
REV. DATE: 12/07/09

SHIPPING WEIGHT: 1277 lbs (579 kg) C1-HARMONY SYSTEMS
SHIPPING WEIGHT: 1052 lbs (477 kg) C1 FRONTAL - BIPLANE SYSTEMS
SHIPPING WEIGHT: 866 lbs (393 kg) C3 LATERAL - BIPLANE SYSTEMS

DETAIL NOT TO SCALE

GE Healthcare
Healthcare Project Implementation - Design Center
Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: INNOVA IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO ASSIST IN THE LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE DETAILS AND DIMENSIONS SHOWN ON THE DRAWINGS. THE USER SHALL BE RESPONSIBLE FOR VERIFYING THE ACTUAL CONSTRUCTION PURPOSES AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY (EP)
FINAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-110F	00
DATE:	18.Dec.13
DRAWN BY:	JPH
CHECKED BY:	TST

REVISION HISTORY:

SHEET
D2

PIM R2
RQ - 140199

EQUIPMENT DETAIL LP AND LC GANTRY DELIVERY PATH

B5050J
REV. DATE: 12/07/09

CORRIDOR WIDTH (FOR GANTRY DELIVERY)	DOOR SIZE (FOR GANTRY DELIVERY)
8'-0" WIDE	4'-0" OPENING
7'-0" WIDE	5'-0" OPENING
6'-0" WIDE	6'-0" OPENING
5'-0" WIDE	7'-0" OPENING
4'-0" WIDE	8'-0" OPENING

NOTE: WHEN DELIVERING GANTRY FROM CORRIDOR TO CORRIDOR THE SAME HOLDS TRUE AS CORRIDOR TO DOOR SIZE.

DETAIL NOT TO SCALE

EQUIPMENT DETAIL INNOVA RADIATION SCATTER PLOTS

B5050P
REV. DATE: 12/07/09

UNITS: RELATIVE AIR KERMA: $\mu\text{Gy}/\mu\text{Gym}^2$
DISTANCES: RADIUS AT 1, 2 AND 3 METERS

NOTE: FOR REFERENCE ONLY. PLEASE REFER TO OPERATOR'S MANUAL FOR ADDITIONAL INFORMATION.

LATERAL 1 METER
LATERAL 1.5 METERS

4-4-4 GANTRY IN LATERAL POSITION - DOSE AT 1.5 METER FROM GROUND
4-4-3 GANTRY IN LATERAL POSITION - DOSE AT 1 METER FROM GROUND

DETAIL NOT TO SCALE

EQUIPMENT DETAIL INNOVA RADIATION SCATTER PLOTS

B5050R
REV. DATE: 12/07/09

UNITS: RELATIVE AIR KERMA: $\mu\text{Gy}/\mu\text{Gym}^2$
DISTANCES: RADIUS AT 1, 2 AND 3 METERS

NOTE: FOR REFERENCE ONLY. PLEASE REFER TO OPERATOR'S MANUAL FOR ADDITIONAL INFORMATION.

VERTICAL 1 METER
VERTICAL 1.5 METERS

4-4-2 GANTRY IN VERTICAL POSITION - DOSE AT 1.5 METER FROM GROUND
4-4-3 GANTRY IN LATERAL POSITION - DOSE AT 1 METER FROM GROUND

DETAIL NOT TO SCALE

EQUIPMENT DETAIL TRAM-RAC 4A

B5047
REV. DATE: 05/26/04

TABLE RAIL MOUNT
FLOOR MOUNT

DETAIL NOT TO SCALE

EQUIPMENT DETAIL CLAB 2 PLUS AMPLIFIER

B5051
REV. DATE: 04/09/03

DETAIL NOT TO SCALE

EQUIPMENT DETAIL VITALIQ COMMUNICATION AND MUSIC SYSTEM

B0566
REV. DATE: 06/14/05

PROCEDURE ROOM

SEE MANUFACTURER'S PREINSTALL MANUAL FOR PRODUCT PLACEMENT AND INSTALLATION.

CONTROL ROOM

VIS-A-VIS, INC.
1-800-319-6014

DETAIL NOT TO SCALE

EQUIPMENT DETAIL MAVIG EYE & THYROID SHIELD WITH LAMP

B50-31E
REV. 00: 10/03/97

CEILING
CARRIAGE TRACK
65" MAX. ARC [1651mm]
65" MAX. ARC [1651mm]
LEAD GLASS SHIELD

DETAIL NOT TO SCALE

EQUIPMENT DETAIL INJECTOR REMOTE CONTROL AND ELECTRONICS

B50-28

REMOTE CONTROL
ELECTRONICS

DRAWING NOT TO SCALE

EQUIPMENT DETAIL INJECTOR ON TABLE RAIL

B50-30A

PLAN VIEW
SIDE VIEW

DRAWING NOT TO SCALE

EQUIPMENT DETAIL XR-BUZZER BRACKET

B5150H
REV. 00: 10/30/08

NOTE: XR-BUZZER BRACKET IS MOUNTED ON WALL, ABOVE CEILING. PLACE SPEAKER ABOVE GRILLED CEILING TILE FOR SOUND PENETRATION.

DETAIL NOT TO SCALE

EQUIPMENT DETAIL VIVID i ULTRASOUND UNIT

B0572
REV. DATE: 28.APR.11

TABLE RAIL MOUNT
SAFE LOCK CART

PHYSICAL SPECIFICATIONS						
Equipment	Width (in/cm)	Depth (in/cm)	Height (in/cm)	Weight (lbs/kg)	Heat (btu/hr)	
Vivid i	14.2"/35.8	12.4"/31.3	2.3"/5.9	111lb/5.4	240	

DETAIL NOT TO SCALE

EQUIPMENT DETAIL SKYTRON LIGHTING UNIT

B2063
REV. DATE: 08/04/08

PLAN VIEW
SIDE VIEW

STRUCTURAL SUPPORT (SEE SHEET "S2") (1/2" DIA. ALL-THREAD ROD)
ELECTRICAL JUNCTION BOX (SUPPLIED WITH LAMP)
FINISHED CEILING

DETAIL NOT TO SCALE

GE Healthcare

Healthcare Project Implementation - Design Center
Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: INNOVA IGS 520, 530, 540

THIS PLAN IS SUBMITTED TO SUPPORT LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE ACTUAL CONSTRUCTION PURPOSES AND THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: ROOM NO. EP
ELECTROPHYSIOLOGY (EP)
FINAL INSTALLATION DRAWINGS

PROJECT	REVISION
5-110F	00

DATE: 18.Dec.13
DRAWN BY: JPH
CHECKED BY: TST

REVISION HISTORY:

SHEET
D3