

Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS	C1
EQUIPMENT LAYOUT	A1
(Equipment locations, heat loads, component weights, environmental specs)	
STRUCTURAL LAYOUT	S1
(Structural support/mounting locations for floor/wall/ceiling, wall support elevations)	
STRUCTURAL DETAILS	S2
(Floor and Ceiling loading information)	
ELECTRICAL LAYOUT	E1
(Contractor supplied wiring, interconnect methods, junction point locations and descriptions)	
ELECTRICAL SPECIFICATIONS	E2
(Maximum wiring run lengths, interconnect diagram, system power specifications)	
ELECTRICAL DETAILS	E3 THRU E4
MECHANICAL LAYOUT	M1
(Chiller information)	
EQUIPMENT DETAILS	D1 THRU D2

These equipment IS drawings indicate the placement and interconnection of the listed equipment components. These drawings are not construction or site preparation drawings. Customer remains ultimately responsible for preparing the site to accommodate the IS and operation of such equipment in compliance with GE Healthcare's written specifications and all applicable federal, state, and/or local requirements.

*** REQUIRED REFERENCE ***
 Optima MR360/
 Brivo MR355
 Pre Installation Manual
 5338503-1EN

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the preIS manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

www.gehealthcare.com/siteplanning

GE Healthcare



MRi Site Planning



imagination at work

Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- It is the customer's responsibility to contract a vibration consultant/engineer to implement site design modifications to meet the GE vibration specification. Refer to the system preinstallation manual for the vibration specification.

GE Equipment Delivery Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist Rev 19					
<i>Before using this document ensure you have the latest Rev from IM/Workshop on D000422762</i>					
GEHC Global Order# : _____		Customer: _____			
GEHC P/N : _____		PE / Installer: _____			
The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.					
		Inspection Date:			
GEHC Minimum Requirements		Storage is item ready?	P/M is item ready?	FE is item ready?	Comments if "N", enter comments or action plan
1	MR Magnet Delivery Requirements: Ensure cryogen venting system is available for magnet connections as defined by GEHC Pre-Installation Manual (PIM) requirements, a dust fan system is installed and operational, 480V power, and chilled water supply is available 24x7 that meets system cooling requirements. External connectivity is available for magnet monitoring and phone service is available during delivery. Surface mount vibrator installed where required. Magnet room final flooring is in place.				
2	MR RF Screen Room Requirements: RF Screen Room is tested with copy of Test Report, emailed to ISAdmin@EMR-GE.com, that it is compliant with GEHC specifications. Dismantle and magnet anchors (if applicable) installed using 2 part anchor. For HDs systems, blow box mount bolts installed by RF vendor using 2 part anchors.				
3	State Regulatory Requirements: Facility registration number provided for states of IL, KY, HI, FL, SC, TX. X-ray shielding plan and state acknowledgment letter provided to installer for AR, DC, NC, SC, CO & WA. Site Drawing Requirements: Final version of equipment network and antenna installation drawings (including red lined versions) verified to match actual room and has been provided to installer.				
4	Surface Penetration Requirements: Customer/Contractor scheduled to provide required drilling or cutting into floors, ceilings, and walls; OR surface penetration permit available and posted in the room when GEHC will perform the work.				
5	Pre-Delivery Route Requirements: The equipment delivery route from the truck to the final destination within the facility has been reviewed with all key stakeholders to satisfy meet the minimum requirements for equipment access, and all communications/notifications have occurred. Arrangements have been made for special handling/elevator, rigging, floor protection, fork lift, millbacktruck, etc.				
6	Finished Room Requirements: Rooms that will contain equipment, including storage areas not in scan suite, are dust free. Provisions taken to maintain a dust free room. Precautions must be taken to prevent dust from entering rooms containing equipment when construction is incomplete in adjacent areas. All walls, primed (final coat not needed on Day 13). Shielding, doors, and windows are to be installed. No contractor work being done during or after the installation that will cause dust in the installation areas or potential equipment damage. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of their security issues, implications and responsibility. For storage, room must meet PIM requirements for storage.				
7	Electrical Requirements: Lockable (LOTO) Main Disconnect Panel (MDP) is installed per GE guidelines and system power is available. Conduits, electrical cable ducting/sliders/cable trays, and access flooring is installed in proper location and height. Surface floor duct and load-side wires can be installed at time of system installation. Validate outlet location and requirements meet specifications for device/equipment.				
8	HVAC Requirements: The HVAC/Chilled Water systems designed to maintain the environment per spec/PIM is at running state and appears to provide the desired environmental conditions including location of vents, temperature and humidity for system operation.				
9	Flooring Requirements: Floor is clean and prepared for final floor covering. Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications. Confirm customer anchoring plan aligns with designed floor thickness. Final flooring installed where required for network racks.				
10	Ceiling Requirements: Unistrut (or equivalent) location, levelness and spacing is measured (or vendor confirmed) and consistent with the requirement of the installation drawings. Ensure Unistrut and rails are not used as mounting surfaces. Ceiling grid is installed. Permanent lighting is installed and operational. HVAC diffusers are installed and connected to ductwork. Ceiling tiles installed per PIM discretion.				

GE Healthcare

IS Services Design Center
 Milwaukee, WI
 Copyright 2009 General Electric Company - Proprietary to GE

SHEET TITLE: **SITE READINESS**
 MODALITY TYPE: **OPTIMA MR360**

THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST GE DRAWING REVISIONS. IT IS THE CUSTOMER'S RESPONSIBILITY TO VERIFY ALL CONSTRUCTION PURPOSES. DURING THE CONSTRUCTION PROCESS, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE: **TYPICAL MR 8-238-C**
 TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
8-238F	03
DATE:	06.SEP.12
DRAWN BY:	PMM
CHECKED BY:	TMS

REVISION HISTORY:

SHEET
C1

PIM R3
RQ - 124576

GE EQUIPMENT LISTING								
EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER : NEITHER A QUOTE OR CON WAS ISSUED AT THE DATE OF THESE DRAWINGS				EQUIPMENT CROSS REFERENCE CHART				
NOTE: LOCAL CONDITIONS MAY DICTATE THAT ITEMS IDENTIFIED IN THIS CATEGORY BE INSTALLED BY OTHERS.								
ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	STRUC PLAN	ELEC PLAN
1	1		PATIENT TABLE - FIXED	264 lbs		B8119	-	S
2	1		SPT PHANTOM CABINET	350 lbs		M6115	-	-
3	1		MAGNET RUNDOWN UNIT	8 lbs		M1715A	MS4	C
4	1		BLOWER BOX	46 lbs	3412 btu	M5715	MSB 15	MG6 S
5	1		1.5 TESLA LCC ACTIVE SHIELD MAGNET	13115 lbs	8191 btu	M6515 M0315E M0115K M0115J	MS1	C
6	1		SYSTEM CABINET	1962 lbs	17064 btu	M0815K		MR2 C
7	1		MESH SHIELD			M0815L		L
8	1		PENETRATION PANEL			M5015G		PP1
9	1		MAGNET MONITOR	22 lbs	204 btu	M1615C		MM C
10	1		SHIELD COOLER CABINET	264 lbs	1706 btu	M30004		MS5 C
11	1		COOLING UNIT	220 lbs	19617 btu	M6015C		WC2 C
12	1		OPERATOR WORKSPACE W/COLOR LCD MONITOR	85 lbs	4948 btu	M3015R		DW C
13	1		OPERATOR WORKSPACE CABINET	198 lbs		M0615D		C
14	1		PATIENT ALERT CONTROL BOX			M4815		PA S

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.								
90	1		11KW AIRSYS CHILLER	804 lbs	50512 btu	B8117 B8116		AWC1 -
95	1		REMOTE CONTROL PANEL FOR 11KW CHILLER			B8118	O	RCP -

SCALE: 1/4" = 1'-0" EQUIPMENT LAYOUT RECOMMENDED CEILING HEIGHT = 8'-9"

MRI SITE PLANNING REMINDERS

PLEASE REFER TO PRE-INSTALLATION CHECKLIST IN PRE-INSTALLATION MANUAL LISTED ON SHEET C1 FOR ITEMS CRITICAL TO IMAGE QUALITY.

- THE LAYOUT SHOULD BE ARRANGED SO THAT THE 5G LINE IS CONTAINED TO THE MAGNET ROOM. IF NOT POSSIBLE, A BARRIER IS RECOMMENDED TO PREVENT ENTRY TO THE 5G FIELD AREA.
- THE SPACES AROUND, ABOVE, AND BELOW THE MAGNET MUST BE REVIEWED FOR EFFECTS OF THE 5G, 3G, 1G, AND 5G FIELDS. REFER TO THE PROXIMITY LIMIT CHART IN THE MR PRE-INSTALLATION MANUAL REFERENCED ON C1.
- FOR MOVING METAL, THE RESTRICTION LINES TYPICALLY EXTEND OUTSIDE OF THE MRI SPACE. PLEASE CONFIRM THERE ARE NO MOVING METAL CONCERNS WITHIN THESE AREAS. AN EMI STUDY IS RECOMMENDED IF THE RESTRICTION LINES ARE VIOLATED.
- FOR VIBRATION ANALYSIS TO BE COMPLETED AS REQUIRED PER PRE-INSTALLATION MANUAL.
- FOR EMI, REVIEW THE SITE FOR THE LOCATION OF THE MAIN ELECTRICAL FEEDERS, AC DEVICES, OR DISTRIBUTION SYSTEMS. AN EMI STUDY IS RECOMMENDED IF LARGE AC SYSTEMS ARE NEARBY.
- DETAILS OF THE FLOOR BELOW THE MAGNET MUST BE REVIEWED. THE STRUCTURAL ENGINEER MUST VERIFY THAT THE QUANTITY OF STEEL IN THE VOLUME 10FT [3.1M] X 10FT [3.1M] X 1FT [.3M] DEEP (BELOW THE MAGNET) DOES NOT EXCEED THE ALLOWABLE STEEL CONTENT AS GIVEN IN THE MR PRE-INSTALLATION MANUAL REFERENCED ON SHEET C1.

RESPONSIBILITY FOR THE COORDINATION, DESIGN, ENGINEERING, AND SITE PREPARATION RESIDES WITH THE CUSTOMER AND THEIR PROJECT ARCHITECTS AND CONTRACTORS. GE DOES NOT, BY PROVIDING REVIEWS AND FURNISHING COMMENTS AND ASSISTANCE, ACCEPT ANY RESPONSIBILITY BEYOND ITS OBLIGATIONS AS DEFINED IN THE MR SYSTEM, SALE/PURCHASE AGREEMENT.

IMAGE QUALITY CONSIDERATIONS

BROADBAND RF NOISE IS A SINGLE TRANSIENT OR CONTINUOUS SERIES OF TRANSIENT DISTURBANCES CAUSED BY AN ELECTRICAL DISCHARGE. LOW HUMIDITY ENVIRONMENTAL CONDITIONS WILL HAVE HIGHER PROBABILITY OF ELECTRICAL DISCHARGE. THE ELECTRICAL DISCHARGE CAN OCCUR DUE TO ELECTRICAL ARCING (MICRO ARCING) OR MERELY STATIC DISCHARGE. SOME POTENTIAL SOURCES CAPABLE OF PRODUCING ELECTRICAL DISCHARGE INCLUDE:

- LOOSE HARDWARE/FASTENERS VIBRATION OR MOVEMENT (ELECTRICAL CONTINUITY MUST ALWAYS BE MAINTAINED)
- FLOORING MATERIAL INCLUDING RAISED ACCESS FLOORING (PANELS & SUPPORT HARDWARE) AND CARPETING
- ELECTRICAL FIXTURES (i.e. LIGHTING FIXTURES, TRACK LIGHTING, EMERGENCY LIGHTING, BATTERY CHARGERS, OUTLETS)
- DUCTING FOR HVAC AND CABLE ROUTING
- RF SHIELD SEALS (WALLS, DOORS, WINDOWS ETC.)

FOR ADDITIONAL INFORMATION REGARDING IMAGE QUALITY, REFER TO THE PRE-INSTALLATION MANUAL LISTED ON SHEET C1.


NOTE: VERIFY DELIVERY ROUTE FOR MAGNET, EQUIPMENT, AND SERVICE EQUIPMENT PRIOR TO DELIVERY.

CRITICAL ITEMS FOR MAGNET DELIVERY

- 24/7 CHILLED WATER AND 480V POWER FOR SHIELD/CRYO COOLER
- 24/7 120V POWER FOR THE MAGNET MONITOR
- PHONE LINES FOR MAGNET MONITORING AND EMERGENCY USE
- MAGNET ROOM EXHAUST FAN
- CRYOGEN VENTING (IF ROOF HATCH, COMPLETED WITHIN 24 HRS)
- MAGNET ANCHORS INSTALLED AND TESTED

THIS IS ONLY A PARTIAL LIST OF ITEMS REQUIRED FOR DELIVERY OF THE MAGNET. FOR A COMPLETE CHECKLIST REFER TO THE PRE-INSTALLATION MANUAL REFERENCED ON SHEET C1.

* THE ISOGAUSS CONTOUR PLOTS DEPICTED ON THIS DRAWING REPRESENT MAGNETIC FRINGE FIELDS RESULTING FROM THE NORMAL OPERATION OF THE MAGNET PROVIDED WITH THE MR SYSTEM. THE ACTUAL MAGNETIC FIELD INTENSITY AT ANY POINT IN THE VICINITY OF THE MAGNET WHEN INSTALLED MAY VARY FROM THE CONTOUR PLOTS DUE TO FACTORS SUCH AS THE CONCENTRATING EFFECTS OF NEARBY FERROUS OBJECTS AMBIENT MAGNETIC FIELDS, INCLUDING THE EARTH'S MAGNETIC FIELD. THEREFORE, THE CONTOURS SHOWN ARE ONLY APPROXIMATIONS OF ACTUAL FIELD INTENSITIES FOUND AT A CORRESPONDING DISTANCE FROM THE MAGNET'S ISOCENTER.



ANCILLARY ITEMS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
60	MINIMUM 9 FT. -0 IN. [2743 mm] X 9 FT. -0 IN. [2743 mm] REMOVABLE WALL SECTION FOR MAGNET DELIVERY/REMOVAL.
61	MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 43 IN. W X 82 IN. H [1092mm X 2083mm]. CENTERING ON A 96 IN. [2438mm] CORRIDOR WIDTH.
62	NON-METAL ACCESS FLOOR WITH 2" X 2" (610 X 610mm) REMOVABLE PANELS & SUPPORT HARDWARE REQUIRED WITHIN MAGNET ROOM.
63	RF SCREEN, INCLUSIVE OF WALLS, FLOOR, DOOR, ETC. GROUND IMPEDANCE GREATER THAN 1000 OHMS. ATTENUATION 100DB AT 150MHZ +/-10MHZ PLANEWAVE.
64	COUNTERTOP WITH DRAWERS FOR MISCELLANEOUS ITEMS.
65	BASE CABINET FOR STORAGE OF: SURFACE COILS, PATIENT POSITIONING PADS, PHANTOMS, ETC.
66	AIR CONDITIONING/COOLING ISOLATION IS RECOMMENDED AT SUPPORTS OF EACH UNIT TO BE INSTALLED.
67	RF FILTERS - LOCATE WITHIN 40 IN. [1016 mm] OF THE RF COMMON GROUND STUD.
68	MAGNET ROOM EXHAUST FAN
69	VALVES AND HOSE BARBS FOR COOLING SYSTEM

THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.

90	DC LIGHTING CONTROL PANEL 135 lbs (60.8 kg) 1024 BTU/HR (300W) (CAT. NO. E4502SC/SE = BASIC SYSTEM)
91	DC LIGHTING AUTO TRANSFORMER 60 lbs (27 kg) 71 BTU/HR (20W) (PART OF VARIABLE DIMMER SYSTEM) (CAT. NO. E4502SD/SF INCLUDES BASIC SYSTEM)
92	METAL DETECTOR (HAND HELD)
93	MAIN DISCONNECT PANEL 900 BTU (264 W) 350 lbs (158kg) CAT. NO. M3088TH
94	WORKSTATION TABLE CAT. NO. M20032FK
95	OPERATOR'S CHAIR CAT. NO. E8803BE

GENERAL SPECIFICATIONS

- THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.
- CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.
- RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.
- THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PREDICATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC., MAKE THIS ESSENTIAL FOR A PROPER SERVICE. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC..
- ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.

SITE ENVIRONMENT SPECIFICATIONS

- AMBIENT OPERATING TEMPERATURE: 59-89.6 DEG (F) [15-32 (C)] FOR THE CONTROL AND EQUIPMENT AREAS, 59-69.8 DEG (F) [15-24 (C)] FOR THE MAGNET ROOM. MAXIMUM ALLOWABLE TEMPERATURE CHANGE OF 5 DEG (F) [3 (C)]/HR.
- HUMIDITY: 30 TO 75 (30-80 FOR THE MAGNET ROOM) PERCENT NON-CONDENSING, MAXIMUM ALLOWABLE CHANGE OF 5 PERCENT/HOUR.
- ALTITUDE: 100 FT [30.5M] BELOW SEA LEVEL TO 8,000 FT. [2438M] ABOVE SEA LEVEL.
- ENVIRONMENTAL RESTRICTIONS ABOVE MUST NOT BE EXCEEDED FOR THE ELECTRONICS.
- DO NOT RESTRICT THE AIR INTAKE OR AIR EXHAUST OF THE SYSTEM COMPONENTS.
- ENVIRONMENTAL CONDITIONS LISTED ABOVE MUST BE MAINTAINED AT ALL TIMES INCLUDING FOR EXAMPLE OVERNIGHT, WEEKENDS, AND HOLIDAYS.
- THE SHIELD COOLER COMPRESSOR CABINET REQUIRES WATER COOLING TO DISSIPATE THE HEAT OUTPUT. HEAT DISSIPATION TO AIR IS NEGLIGIBLE. 24 HOUR POWER AND WATER COOLING MUST BE AVAILABLE UPON MAGNET DELIVERY.
- CRYOGEN VENTING AND MAGNET ROOM EXHAUST FAN SYSTEMS MUST BE COMPLETED IN THE MAGNET ROOM PRIOR TO DELIVERY.
- FLUORESCENT LIGHTING IS NOT ALLOWED IN THE MAGNET ROOM DUE TO RF NOISE.

MAGNETIC INTERFERENCE SPECIFICATIONS

- THE CUSTOMER MUST ESTABLISH PROTOCOLS TO PREVENT PERSONS WITH CARDIAC PACEMAKERS, NEUROSTIMULATORS, AND BIOSTIMULATION DEVICES FROM ENTERING MAGNETIC FIELDS OF GREATER THAN 5 GAUSS (EXCLUSION ZONE).
- MAIN POWER TRANSFORMERS MUST REMAIN OUTSIDE THE 3 GAUSS FIELD. EMI < 40mG AC, EMI < 4.43mG DC.
- POTENTIAL EXISTS UNDER FAULT CONDITIONS THAT THE 5 GAUSS LINE MAY EXPAND RADIALLY TO 16.4 FT. [5.0 m] AND AXIALLY TO 22.96 FT. [7.0 m] FOR 2 SECONDS OR LESS. IT SHOULD BE NOTED THAT NORMAL RAMPDOWNS OR MRI (MAGNET RUNDOWN UNIT) INITIATED QUENCHES WILL NOT CAUSE THE MAGNETIC FIELD TO EXPAND.
- IT IS RECOMMENDED EVERY SITE CONSIDER THE EVENT OF A QUENCH AND PLAN ACCORDINGLY (SUCH AS PLACING 5 GAUSS WARNING SIGNS AT EXPANDED LOCATIONS).
- THE FERROUS METAL OBJECTS LISTED BELOW MUST NOT MOVE INTO OR INSIDE OF THE MOVING METAL SENSITIVITY LINE DURING SCANS.

TYPICAL MOVING MAGNETIC MASS	DISTANCE RADIALLY	DISTANCE AXIALLY
CARTS, GURNEYS 100-400 lbs [45-182 kg]	3 GAUSS LINE	3 GAUSS LINE
FORKLIFTS, SMALL ELEVATOR, CARS, MINIVANS VANS, PICKUP TRUCKS, AMBULANCES (OBJECTS GREATER THAN 400 lbs [182 kg])	15.5 ft. [4.72 m]	21.0 ft. [6.4 m]
BUSES AND TRUCKS (DUMP, TRACTOR TRAILER, UTILITY, FIRE TRUCKS)	18.1 ft. [5.52 m]	24.5 ft. [7.47 m]

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED

GE Healthcare

Healthcare Project Implementation - Design Center

SHEET TITLE: EQUIPMENT LAYOUT MODALITY TYPE: OPTIMA MR360

PROJECT TITLE: TYPICAL MR 8-238-C TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
8-238F	03


DATE: 06.SEP.12
DRAWN BY: PMM
CHECKED BY: TMS

REVISION HISTORY:

SHEET
A1


TYPICAL WALL SUPPORT ELEVATIONS

S63




SUPPORT FOR
MAGNET RUNDOWN UNIT
(NOT TO SCALE)

S86




SUPPORT FOR
MAGNET MONITOR
(NOT TO SCALE)

S62



SUPPORT FOR
DC LIGHTING CONTROLLER
(NOT TO SCALE)

S60





SUPPORT FOR
MAIN DISCONNECT CONTROL
(NOT TO SCALE)

SCALE: 1/4" = 1'-0"

STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 8'-9"



STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S63, FOR MAGNET RUNDOWN UNIT.
2	FLOOR MOUNTING AREA, FOR BLOWER BOX. SEE DETAIL MSB-15 ON SHEET S2.
3	MAGNET FLOOR MOUNTING, SEE DETAIL M6615A2 ON SHEET S2 FOR MORE INFORMATION.
4	FLOOR LEVELNESS - FLOOR AREA MUST BE HARD. FLOOR SLOPE < 1/8" IN. (3.2mm) PER 10' (3.0m). FLOOR SURFACE FINISH SHALL BE 0.157 IN. (4mm).
5	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S86, FOR MAGNET MONITOR.
6	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S62, FOR DC LIGHTING CONTROLLER.
7	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S60, FOR MAIN DISCONNECT CONTROL.
8	SEE DETAIL B5031A ON THE S2 SHEET.
9	SEE DETAIL M05-15K ON SHEET S2 FOR FLOOR MOUNTING OF OPERATOR WORKSPACE.

STRUCTURAL NOTES

- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CERTAIN MR PROCEDURES REQUIRE AN EXTREMELY STABLE ENVIRONMENT TO ACHIEVE HIGH RESOLUTION IMAGE QUALITY. VIBRATION IS KNOWN TO INTRODUCE FIELD INSTABILITIES INTO THE IMAGING SYSTEM. THE VIBRATION EFFECTS ON IMAGE QUALITY CAN BE MINIMIZED DURING THE INITIAL SITE PLANNING OF THE MR SUITE BY MINIMIZING THE VIBRATION ENVIRONMENT. SEE MOUNTING DETAIL ON SHEET S2 FOR ADDITIONAL INFORMATION.
- STANDARD STEEL STUDS, NAILS, SCREWS, CONDUIT, PIPING, DRAINS AND OTHER HARDWARE ARE ACCEPTABLE IF PROPERLY SECURED. ANY LOOSE STEEL OBJECTS CAN BE VIOLENTLY ACCELERATED INTO THE BORE OF THE MAGNET. CAREFUL THOUGHT SHOULD BE GIVEN TO THE SELECTION OF LIGHT FIXTURES, CABINETS, WALL DECORATIONS, ETC. TO MINIMIZE THIS POTENTIAL HAZARD. FOR SAFETY, ALL REMOVABLE ITEMS WITHIN THE MAGNET ROOM SUCH AS FAUCET HANDLES, DRAIN COVERS, SWITCH BOX COVER PLATES, LIGHT FIXTURE COMPONENTS, MOUNTING SCREWS, ETC. MUST BE NON-MAGNETIC. IF YOU HAVE A SPECIFIC QUESTION ABOUT MATERIAL, BRING IT TO THE ATTENTION OF YOUR GE PROJECT MANAGER OF INSTALLATIONS.
- FLOOR LEVELNESS IN THE MAGNET ROOM SHOULD NOT EXCEED 0.3125 in. (8 mm) WHEN MEASURING BETWEEN DEPRESSIONS AND HIGH SPOTS OVER ANY 120 in. (3048 mm) DISTANCE WITHIN THE 87.5 in. (2178 mm) BY 139.3 in. (3539 mm) AREA OF THE MAGNET ENCLOSURE AND THE AREA IN FRONT OF THE ENCLOSURE. THIS FLOOR LEVELNESS REQUIREMENT IS IMPORANT FOR ACCURATE PATIENT TABLE DOCKING.
- NON-MOVABLE STEEL SUCH AS WALL STUDS OR HVAC COMPONENTS WILL PRODUCE NEGLIGIBLE EFFECT ON THE ACTIVE SHIELD MAGNET.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.
- CUSTOMERS CONTRACTOR TO PROVIDE AND INSTALL APPROPRIATE SUPPORTS FOR THE STORAGE OF EXCESS CABLES.
- IT IS THE CUSTOMER'S RESPONSIBILITY TO PERFORM ANY FLOOR OR WALL PENETRATIONS THAT MAY BE REQUIRED. THE CUSTOMER IS ALSO RESPONSIBLE FOR ENSURING THAT NO SUBSURFACE UTILITIES (E.G., ELECTRICAL OR ANY OTHER FORM OF WIRING, CONDUITS, PIPING, DUCT WORK OR STRUCTURAL SUPPORTS (I.E. POST TENSION CABLES OR REBAR)) WILL INTERFERE OR COME IN CONTACT WITH SUBSURFACE PENETRATION OPERATIONS (E.G. DRILLING AND INSTALLATION OF ANCHORS/SCREWS) PERFORMED DURING THE INSTALLATION PROCESS. TO ENSURE WORKER SAFETY, GE INSTALLERS WILL PERFORM SURFACE PENETRATION OPERATIONS ONLY AFTER THE CUSTOMER'S VALIDATION AND COMPLETION OF THE "GE SURFACE PENETRATION PERMIT"

SHEET TITLE: STRUCTURAL LAYOUT

MODALITY TYPE: OPTIMA MR360

THIS PLAN IS SUBMITTED TO SUPPORT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE DETAILS AND REQUIREMENTS OF THE GE HEALTHCARE PROJECT IMPLEMENTATION DESIGN CENTER. GE HEALTHCARE PROJECT IMPLEMENTATION DESIGN CENTER CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
TYPICAL MR
8-238-C
TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
8-238F	03





DATE:	06.SEP.12
DRAWN BY:	PMM
CHECKED BY:	TMS

REVISION HISTORY:

SHEET
S1

GE Healthcare
Healthcare Project Implementation - Design Center
Milwaukee, Wisconsin

PIM R3
RQ - 124576



ACOUSTICS AND VIBRATION GUIDELINES: SIGNA LCC MAGNET - (CXK4)

M66-15G
REV. DATE: 08/22/05

SYSTEM ACOUSTIC NOISE LEVELS

ANY GE FACTORY-INSTALLED PROTOCOL CAN BE MODIFIED BY OPERATORS, WHICH CAN INCREASE OR DECREASE ACOUSTIC SPL (SOUND PRESSURE LEVEL); OR OPERATORS MAY CREATE THEIR OWN PROTOCOL WHICH COULD PRODUCE A HIGHER OR LOWER ACOUSTIC SPL AS STATED UNDER OPERATING CONDITIONS CONDITION 1 BELOW. TYPICAL SCANS GENERATE ACOUSTIC LEVELS AS STATED UNDER OPERATING CONDITIONS CONDITION 2 BELOW. IN ADDITION, THE EXPOSURE TIMES ARE COMPLETELY UNDER OPERATOR CONTROL. CONSEQUENTLY, HEARING PROTECTION IS REQUIRED FOR ALL PEOPLE IN THE MAGNET ROOM DURING SCANS TO PREVENT HEARING IMPAIRMENT. ACOUSTIC LEVELS MAY EXCEED 99 dBA. AGAIN, FOR MORE INFORMATION ABOUT RECOMMENDED SAFETY PROCEDURES REGARDING PATIENT EXPOSURE TO MR-GENERATED ACOUSTIC NOISE, SEE THE MR SAFETY GUIDE INCLUDED IN THE USER MANUAL.

AMBIENT CONDITIONS

TO REDUCE ANY BACKGROUND NOISE DUE TO CABINET BLOWERS, ETC., ACOUSTICAL CEILINGS, WALLS, AND FLOORS ARE RECOMMENDED. THE FOLLOWING ARE TYPICAL NOISE LEVEL READINGS:

- OPERATOR AREA55 dBA
- EQUIPMENT ROOM75 dBA
- MRCC (MR COMMON CHILLERS).....69.1 dBA

OPERATING CONDITIONS

CONDITION 1

MR SCANNERS UNDER "WORST-CASE" OPERATING CONDITIONS, COULD GENERATE ACOUSTIC LEVELS (AS MEASURED AT THE MAGNET ISO-CENTER) AS FOLLOWS:

AVERAGE SPL 118 dBA SPL = SOUND PRESSURE LEVEL
PEAK 128 dB
FREQUENCY RANGE 20 TO 20k Hz

CONDITION 2

MR SCANNERS FOR MANY TYPICAL CLINICAL SCANNING SCENARIOS THOUGH, GENERATE ACOUSTIC LEVELS (AS MEASURED AT THE MAGNET ISO-CENTER) SOMEWHAT LOWER AS FOLLOWS:

AVERAGE SPL 95 TO 110 dBA
PEAK 110 TO 120 dB
FREQUENCY RANGE 20 TO 20k Hz

AS RECENT HISTORY HAS SHOWN AN EVOLUTION TOWARDS MORE POWERFUL (AND HENCE LOUDER) GRADIENT SUBSYSTEMS, ARCHITECTS SHOULD CONSIDER THE ACOUSTIC LEVELS STATED IN THE "WORST CASE" CONDITION 1, MENTIONED ABOVE. NOTE THAT HIGH-FIELD SIGNA SYSTEMS HAVE THE ABILITY TO RUN SCANNING PROTOCOLS WHICH CAN GENERATE ACOUSTIC LEVELS OVER THE ENTIRE HUMAN PERCEPTIBLE FREQUENCY RANGE (20 TO 20k Hz), THEREFORE ATTENUATION OVER THIS ENTIRE RANGE MUST BE CONSIDERED FOR SITE DESIGN.

TEST MEASUREMENTS (1.1)

VIBRATION MEASUREMENTS ARE IN THE RANGE OF 10⁻⁶g. TEST EQUIPMENT MUST HAVE THE REQUIRED SENSITIVITY TO THESE LEVELS.

INSTRUMENTATION IS RECOMMENDED TO HAVE A LOW TOLERANCE TO TEMPERATURE EFFECTS AS MANY TIMES THE LOW FREQUENCY THERMAL DRIFT MAY INFLUENCE THE MEASUREMENTS.

IT IS HIGHLY RECOMMENDED ALL MEASURED DATA IS REAL TIME DATA ACQUISITION. RECORDING THE VIBRATION DATA WILL NOT ALLOW FOR A PROPER SITE SURVEY, SPECIFICALLY WHEN STUDYING TRANSIENT VIBRATION AND WHEN SEARCHING FOR SPECIFIC VIBRATION SOURCES.

ALL ANALYSES ARE TO BE NARROWBAND FAST FOURIER TRANSFORMS (FFT'S) OVER THE FREQUENCY BANDS LISTED BELOW:

FREQUENCYBAND	FREQUENCY RESOLUTION
0.2 TO 50 HZ	Δf = 0.125 HZ

TIME HISTORIES OF THE VIBRATION MUST BE RECORDED AS ACCELERATION LEVELS VS. TIME. THE RESOLUTION OF THE TIME HISTORY MUST BE ADJUSTED TO CLEARLY CAPTURE THE TRANSIENT EVENT. THE ANALYZER SET-UP WILL BE SITE DEPENDENT AND, IN SPECIAL CASES, VIBRATION RESPONSE DEPENDENT. IT IS THE RESPONSIBILITY OF THE VIBRATION CONSULTANT TO STUDY THE TRANSIENT ENVIRONMENT, CAPTURE DATA TO CONFIRM TRANSIENT ACTIVITY EXCEEDS THE TRIGGER LEVEL, THEN EXPAND THE TIME HISTORY DATA TO EXHIBIT THE STRUCTURAL RESPONSE.

EQUIPMENT (SPECTRAL ANALYZER) SET-UP (1.2)

- FREQUENCY AVERAGE A MINIMUM OF 20 LINEAR AVERAGES. DO NOT USE PEAK HOLD OR 1/3 OCTAVE ANALYSIS.
- AVERAGE AND STORE A MINIMUM OF 10 PLOTS TO SUPPORT THE SITE VIBRATIONS CONSISTENCY.
- HANNING WINDOW MUST BE APPLIED TO THE ENTIRE SPECTRA

SPECTRUM ANALYZERS CAPABLE OF THESE MEASUREMENTS ARE READILY AVAILABLE FOR PURCHASE OR RENTAL. MODELS SUCH AS THE HP 3560A, NICOLET PHASZER, B&K PULSE, AND HP 35670 ARE ALL CAPABLE OF MAKING THE SITE VIBRATION MEASUREMENTS. ACCELEROMETERS MUST HAVE THE CAPABILITY TO MEASURE FROM 0.2 HZ BEYOND 50 HZ. TIME HISTORIES CAN BE RECORDED USING ANY OF THE ANALYZERS LISTED ABOVE. PLEASE NOTE THAT THE EQUIPMENT MENTIONED ARE FOR EXAMPLE ONLY. IT IS THE RESPONSIBILITY OF THE ENGINEERING TEST FIRM TO PROVIDE EQUIPMENT THAT WILL ALLOW MEASUREMENTS COMPLIANT WITH THIS GUIDELINE.

DATA COLLECTIONS (1.3)

AMBIENT BASELINE CONDITION:

ALL OF THE MEASUREMENTS DEFINED IN 1.1 AND 1.2 (ABOVE) MUST BE MADE IN A 'QUIET' ENVIRONMENT, THAT IS, IN AREAS WHERE EXCESSIVE TRAFFIC, SUBWAY TRAINS, ETC. EXISTS. A VIBRATION MEASUREMENT MUST ALSO BE MADE DURING PERIODS WITHOUT TRAFFIC OR DURING PERIODS OF LIGHT TRAFFIC. MEASUREMENTS MUST DEFINE THE LOWEST LEVELS OF VIBRATION POSSIBLE AT THE SITE.

THE SOURCE OF ANY STEADY STATE VIBRATION WHOSE LEVELS EXCEED THE SPECIFICATIONS MUST BE IDENTIFIED AS TO THE SOURCE OF THE VIBRATION DISTURBANCE. A SECOND MEASUREMENT SHOULD BE MADE WITH ALL OF THE IDENTIFIED CONTRIBUTORS POWERED DOWN IF POSSIBLE. IN SITUATIONS WHERE IT IS NOT POSSIBLE TO POWER DOWN EQUIPMENT, VIBRATION DATA MUST BE COLLECTED TO IDENTIFY SPECIFIC SOURCE OF THE VIBRATION CONCERN. THE MAJORITY OF STEADY STATE VIBRATION PROBLEMS CAN BE NEGATED BY ISOLATING THE VIBRATION SOURCE.

NORMAL CONDITION

ALL OF THE VIBRATION MEASUREMENTS LISTED ABOVE MUST BE REPEATED DURING PERIODS OF 'NORMAL' ENVIRONMENTAL CONDITIONS INCLUDING THE FFT'S AND TIME HISTORIES. THE TRANSIENT MEASUREMENTS MUST BE PROVIDED TO DEFINE THE DYNAMIC DISTURBANCES THE MR SYSTEM MIGHT BE EXPOSED TO. TRANSIENT ANALYSIS IS REQUIRED FOR A TRUE ASSESSMENT OF THE SITE.

SPECIAL ATTENTION MUST BE PAID TO THE SITE ASSESSMENT DURING THE ENTIRE ANALYSIS. SINCE TRANSIENT VIBRATION IS NOT EASILY ADDRESSED ONCE THE MR SUITE IS FULLY CONSTRUCTED, THE TEST CONSULTANT MUST FULLY UNDERSTAND THE NEEDS FOR THIS ANALYSIS. THE SOURCE OF ANY TRANSIENT MUST BE IDENTIFIED AND SUPPORTED WITH VIBRATION PLOTS. IF THE SOURCE OF ANY TRANSIENT IS NOT ABLE TO BE LOCATED, IT IS RECOMMENDED THAT THE CUSTOMER SHOULD HAVE AN ALTERNATE LOCATION IDENTIFIED AND VIBRATION STUDIED.

TRANSIENT VIBRATION IS DIFFICULT TO ASSESS IF THE DETAILS OF THE TRANSIENT VIBRATION IS NOT UNDERSTOOD. THE 0.0005 g, ZERO TO PEAK TRIGGER LEVEL IS A STARTING POINT TO BEGIN UNDERSTANDING THE VIBRATION STABILITY. THE TRANSIENT VIBRATION PEAK AMPLITUDE, STRUCTURAL (TIME VARIANT) RESPONSE, DECAY RATE AND AN ESTIMATE OF THE NUMBER OF EVENTS PER UNIT TIME WOULD CONSTITUTE A COMPLETE TRANSIENT ANALYSIS. ALL TRANSIENT FAILURES MUST BE SUPPORTED BY TIME HISTORY PLOTS. THE PLOTS MUST CLEARLY SHOW THE STRUCTURAL RESPONSE, THE FREQUENCY OF THE SIGNATURE AND THE DECAY RATE. FROM THIS DATA, GE CAN HELP DETERMINE COMPLIANCE TO THE VIBRATION GUIDELINES.

TEST CONSULTANT MUST PROVIDE DESIGN RECOMMENDATIONS FOR ALL SITES/BUILDING STRUCTURES WHICH ARE FOUND TO EXCEED THE SPECIFICATIONS.

PRESENTATION/INTERPRETATION OF RESULTS (1.4)

THE RECOMMENDED FORMAT FOR SITE VIBRATION DATA COLLECTION, PRESENTATION, AND ANALYSIS IS ILLUSTRATED IN THE EXAMPLES SHOWN IN ILLUSTRATIONS 1.1 THROUGH 1.4. IN THE PRE-INSTALLATION MANUAL. PRESENTATION OF THE DATA IN ANY OTHER FORMAT (LINEAR UNITS ONLY) MAY RESULT IN AN INCORRECT INTERPRETATION AND DIAGNOSIS OF THE SITE. ADDITIONAL DATA COLLECTION OR PRESENTATION METHODS IS AT THE OPTION OF THE VIBRATION TESTING SERVICE.

IT IS THE RESPONSIBILITY OF THE CUSTOMER'S VIBRATION TESTING SERVICE TO INTERPRET THE RESULTS AND DETERMINE IF THAT SITE MEETS SPECIFICATIONS. ILLUSTRATIONS A-1 AND A-2 ARE EXAMPLES PROVIDED TO ASSIST A TEST CONSULTANT IN THE USE OF GE STEADY STATE SPECIFICATIONS (VIBRATION SPECIFICATIONS ABOVE AMBIENT BASELINE). IF THE VIBRATION LEVELS ARE TOO HIGH, ADDITIONAL DATA ACQUISITION MAY BE NECESSARY TO:



- DETERMINE THE SOURCE OF THE VIBRATION
- PROPOSE A SOLUTION TO THE PROBLEM
- FIND AN ALTERNATE SITE LOCATION.

ILLUSTRATIONS A-3 AND A-4 IN THE PRE-INSTALLATION MANUAL ARE EXAMPLES PROVIDED TO ASSIST A TEST CONSULTANT IN THE USE OF GE TRANSIENT SPECIFICATIONS. THE 500 MICRO-G, ZERO TO PEAK TRIGGER LEVELS DATA COLLECTION TO BEGIN ASSESSMENT OF THE SITE VIBRATION ANALYSIS. THE RESPONSE OF THE TRANSIENT MUST BE ASSESSED RELATIVE TO THE STEADY STATE VIBRATION SPECIFICATIONS IN SECTION SPECIFICATIONS.

ANY QUESTIONS REGARDING TEST EQUIPMENT REQUIREMENTS, TEST PARAMETERS, OR GENERAL QUESTIONS SHOULD BE DISCUSSED WITH YOUR GE INSTALLATION SPECIALIST.

FLOOR MOUNTING DETAIL: SIGNA LCC MAGNET (CXK4)

REFER TO SHEET A1 FOR ACTUAL MAGNET ORIENTATION
M6615A2
REV. DATE: 10/03/08



VIBRATION

- THE MAGNET MAY BE SENSITIVE TO VIBRATIONS IN THE FREQUENCY RANGE OF 0.5 TO 45 HZ DEPENDING ON THE AMPLITUDE OF THE VIBRATION. IN THE PHYSICAL AREA WHERE THE MR SYSTEM IS TO BE LOCATED, EVERY PRECAUTION MUST BE TAKEN TO ENSURE THAT THE VIBRATION IS MINIMIZED. IN THE MAGNET SITING AREA, THE STRUCTURAL STABILITY AND BEHAVIOR CHARACTERISTICS CAN BE ASSESSED. THE VIBRATION TESTS OUTLINED CAN BE USED TO ASSESS THE VIBRATION ENVIRONMENT. SITES WHICH CURRENTLY PASS THE VIBRATION STABILITY CRITERIA MAY PROCEED WITH INSTALLATION. SITES WHICH HAVE MARGINAL VIBRATION STABILITY REQUIRE SOURCE ISOLATION OR STRUCTURAL MODIFICATIONS. THEN IT IS THE CUSTOMER'S RESPONSIBILITY TO CONTRACT A VIBRATION CONSULTANT OR QUALIFIED ENGINEER TO IMPLEMENT DESIGN MODIFICATIONS TO MEET THE SPECIFIED LIMITS. WITH THE VIBRATION CONSULTANT PRESENT, LOCAL GE FIELD SERVICE AND/OR INSTALLATION SPECIALIST MUST VERIFY THE ELIMINATION/REDUCTION OF ALL IDENTIFIED SOURCES DO IMPROVE THE VIBRATION ENVIRONMENT. GE CAN ASSIST IN INTERPRETING MARGINAL SITE TEST RESULTS AND PREDICTING THE IMPACT ON SYSTEM PERFORMANCE. HOWEVER IT IS ULTIMATELY THE CUSTOMER/ARCHITECT/ENGINEER RESPONSIBILITY TO DESIGN SITE SOLUTION.
- TO MINIMIZE THE INTERFERENCE, THE MAGNET SHOULD BE PLACED ON A SOLID FLOOR, LOCATED AS FAR AS POSSIBLE FROM THE VIBRATION SOURCES, SUCH AS PARKING LOTS, ROADWAYS, SUBWAYS, TRAINS, HALLWAYS, ELEVATORS, HELIPORTS AND HOSPITAL PHYSICAL PLANTS CONTAINING PUMPS, MOTORS, AIR HANDLING EQUIPMENT, OR AIR CONDITIONING EQUIPMENT.

PLEASE NOTE THAT OTHER ITEMS NOT LISTED COULD ALSO BE POTENTIAL SOURCES OF VIBRATION.

VIBRATION ISOLATION IS RECOMMENDED AT FLOOR CONNECTION POINTS OF THE AIR CONDITIONING UNIT(S) TO BE INSTALLED FOR THE PURPOSE OF COOLING THE MR SUITE.

ISOLATION OF THE MR MAGNET IS NOT A RECOMMENDED SOLUTION FOR REDUCING ENVIRONMENTAL VIBRATION.

- VIBRATION MEASUREMENTS SHOULD BE MADE WHEN THE PROPOSED SITE IS LOCATED NEAR ANY OF THE SOURCES LISTED HERE. MEASUREMENTS SHOULD BE MADE USING A SPECTRUM ANALYZER CAPABLE OF PERFORMING THE TEST GUIDELINES.

MAGNET SITING REQUIREMENT

- THE MAGNET MUST BE RIGIDLY BOLTED TO THE FLOOR. VIBRATION MEASUREMENTS ON THE MAGNET SUPPORT MUST MEET THE GUIDELINES BELOW. CUSTOMER/CONTRACTOR IS RESPONSIBLE FOR THE PROPER MAGNET ANCHORING.

TRANSIENT VIBRATION

- TIME HISTORY VIBRATION LEVELS (WITH ALL STEADY STATE VIBRATION SOURCES POWERED DOWN) EXCEEDING TRIGGER OF 0.0005 g, ZERO TO PEAK MUST BE FULLY ANALYZED TO ASSESS THE POTENTIAL IMPACT TO THE BUILDING STRUCTURE. THE BUILDING (SPECTRAL) RESPONSE IMMEDIATELY FOLLOWING THE 0.0005 g, ZERO TO PEAK TRIGGER LEVEL (ENDING AT THE DECAY OF THE VIBRATION SIGNAL) MUST NOT CAUSE THE SITE ENVIRONMENT TO EXCEED THE STEADY STATE VIBRATION LEVELS DEFINED BELOW.

STEADY STATE VIBRATION

- THE MAXIMUM STEADY STATE VIBRATION TRANSMITTED THROUGH THE FLOOR MUST NOT EXCEED THE FOLLOWING MAXIMUM SINGLE FREQUENCY COMPONENTS ABOVE AMBIENT BASELINE:
 - 5 x 10⁻⁵ g rms at 0 Hz ramping to 10 x 10⁻⁵ g at 20 Hz
 - 10 x 10⁻⁵ g rms 20-40 Hz
 - 45 x 10⁻⁵ g rms 40-50 Hz


IN ORDER TO ENSURE THAT ANY DISCRETE SIGNAL REPRESENTS A REAL MECHANICAL VIBRATION SOURCE, THE SIGNAL MUST HAVE A BANDWIDTH THAT TYPICALS DYNAMIC SYSTEM RESPONSE.

PROJECT	REVISION
8-238F	03

DATE: 06.SEP.12
DRAWN BY: PMM
CHECKED BY: TMS

REVISION HISTORY:

INTERCONNECT DIAGRAM



POWER SPECIFICATIONS

SIGNA OPTIMA MR360 MR355/BRIVO (TYPE C)
 (REV. DATE 29Apr.11)

VOLTAGE
 PRIMARY SOURCE IS REQUIRED FOR ALL INSTALLATIONS. RANGE OF LINE VOLTAGES: NOMINAL LINE VOLTAGE OF 380 TO 480, 3 PHASE, 50 OR 60 HZ.
 RECOMMENDED POWER SUPPLY: WYE-CONNECTED OR DELTA-CONNECTED (GROUNDING DELTA).
 MAXIMUM DAILY VOLTAGE VARIATION MUST FALL WITHIN ONE OF THE RANGES IN TABLE A.

TABLE A ALLOWABLE INPUT VOLTAGES/CURRENT DEMAND

NOMINAL VOLTAGE	ABSOLUTE RANGE	CURRENT (AMPS)		MINIMUM STANDARD OVERCURRENT PROTECTION **
		MAX MOMENTARY	CONTINUOUS	
380	342-418	83	68	150-A
400	360-440	79	64	150-A
415	374-456	76	62	150-A
480	432-528	66	53	150-A

** OVERCURRENT PROTECTION SIZED FOR 125% CONTINUOUS CURRENT. (CALCULATIONS BASED UPON NOMINAL VOLTAGE).

PHASE-BALANCE.

PHASE-TO-PHASE VOLTAGES MUST BE WITHIN 2 PERCENT OF THE LOWEST PHASE-TO-PHASE VOLTAGE. MAXIMUM ALLOWABLE TRANSIENT VOLTAGE EXCURSIONS ABOVE OR BELOW NOMINAL WAVESHAVE FORM NOT TO EXCEED 200V AT A MAXIMUM DURATION OF 1 CYCLE AND FREQUENCY OF 10 TIMES PER HOUR. VOLTAGE TRANSIENT OR IMPULSE ON THE INCOMING POWER MUST BE HELD TO A MINIMUM. TRANSIENTS CAUSED BY LIGHTNING, SURGES, LOAD SWITCHING, STATIC ELECTRICITY ETC. CAN CAUSE SCAN ABORTS OR IN EXTREME INSTANCES, COMPONENT FAILURE IN THE COMPUTER SUBSYSTEM.

POWER DEMAND

MAXIMUM POWER DEMAND AVERAGED OVER 5 SECONDS = 54.4 KVA. 64.7 KVA CONSISTING OF 30 KVA FOR PDU + 10.9 KVA (CONTINUOUS OPERATION) FOR MRCC + 9 KVA (CONTINUOUS OPERATION) FOR SHIELD/CRYO COOLER + 4.5 KVA FOR MAGNET MONITOR EQUIPMENT.

TABLE B MAXIMUM POWER DEMAND.

DEMAND	SIGNA TWINSPED
kVA*	54.4
POWER FACTOR AT	0.9

* DEMAND INCLUDES POWER FOR ENTIRE MR SYSTEM. LINE VOLTAGE REGULATION AT MAXIMUM POWER DEMAND MUST BE LESS THAN OR EQUAL TO 2 PERCENT OR 4 PERCENT FROM POWER SOURCE.

DISTRIBUTION TRANSFORMER

FOR A SINGLE UNIT INSTALLATION, THE MINIMUM TRANSFORMER SIZE IS 112.5 KVA. REGULATED TRANSFORMER IS NOT REQUIRED UNLESS VOLTAGE CHANGES EXCEED ±10% OVER A PERIOD OF 1 HOUR OR LONGER.

REFER TO DIRECTION LISTED ON C1 FOR ADDITIONAL INFORMATION.

ELECTRICAL NOTES

- NOTE 1: ALL WIRES SPECIFIED SHALL BE COPPER STRANDED, FLEXIBLE, THERMO-PLASTIC, COLOR CODED, CUT 10 FOOT LONG AT OUTLET BOXES, DUCT TERMINATION POINTS OR STUBBED CONDUIT ENDS. ALL CONDUCTORS, POWER, SIGNAL AND GROUND, MUST BE RUN IN A CONDUIT OR DUCT SYSTEM. ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS. WIRE RUNS MUST BE CONTINUOUS COPPER STRANDED AND FREE FROM SPLICES. ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.
- NOTE 2: WIRE SIZES GIVEN ARE FOR USE OF EQUIPMENT. LARGER SIZES MAY BE REQUIRED BY LOCAL CODES.
- NOTE 3: IT IS RECOMMENDED THAT ALL WIRES BE COLOR CODED, AS REQUIRED IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 4: CONDUIT SIZES SHALL BE VERIFIED BY THE ARCHITECT, ELECTRICAL ENGINEER OR CONTRACTOR, IN ACCORDANCE WITH LOCAL OR NATIONAL CODES.
- NOTE 5: CONVENIENCE OUTLETS ARE NOT ILLUSTRATED. THEIR NUMBER AND LOCATION ARE TO BE SPECIFIED BY OTHERS. LOCATE AT LEAST ONE CONVENIENCE OUTLET CLOSE TO THE SYSTEM CONTROL, THE POWER DISTRIBUTION UNIT AND ONE ON EACH WALL OF THE PROCEDURE ROOM. USE HOSPITAL APPROVED OUTLET OR EQUIVALENT.
- NOTE 6: GENERAL ROOM ILLUMINATION IS NOT ILLUSTRATED. CAUTION SHOULD BE TAKEN TO AVOID EXCESSIVE HEAT FROM OVERHEAD SPOTLIGHTS. DAMAGE CAN OCCUR TO CEILING MOUNTING COMPONENTS AND WIRING IF HIGH WATTAGE BULBS ARE USED. RECOMMEND LOW WATTAGE BULBS NO HIGHER THAN 75 WATTS AND USE DIMMER CONTROLS (EXCEPT MR). DO NOT MOUNT LIGHTS DIRECTLY ABOVE AREAS WHERE CEILING MOUNTED ACCESSORIES WILL BE PARKED.
- NOTE 7: ROUTING OF CABLE DUCTWORK, CONDUITS, ETC., MUST RUN DIRECT AS POSSIBLE OTHERWISE MAY RESULT IN THE NEED FOR GREATER THAN STANDARD CABLE LENGTHS (REFER TO THE INTERCONNECTION DIAGRAM FOR MAXIMUM USABLE LENGTHS POINT TO POINT).
- NOTE 8: CONDUIT TURNS TO HAVE LARGE, SWEEPING BENDS WITH MINIMUM RADIUS IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 9: A SPECIAL GROUNDING SYSTEM IS REQUIRED IN ALL PROCEDURE ROOMS BY SOME NATIONAL AND LOCAL CODES. IT IS RECOMMENDED IN AREAS WHERE PATIENTS MIGHT BE EXAMINED OR TREATED UNDER PRESENT, FUTURE, OR EMERGENCY CONDITIONS. CONSULT THE GOVERNING ELECTRICAL CODE AND CONFER WITH APPROPRIATE CUSTOMER ADMINISTRATIVE PERSONNEL TO DETERMINE THE AREAS REQUIRING THIS TYPE OF GROUNDING SYSTEM.
- NOTE 10: THE MAXIMUM POINT TO POINT DISTANCES ILLUSTRATED ON THIS DRAWING MUST NOT BE EXCEEDED.
- NOTE 11: PHYSICAL CONNECTION OF PRIMARY POWER TO GE EQUIPMENT IS TO BE MADE BY CUSTOMERS ELECTRICAL CONTRACTOR WITH THE SUPERVISION OF A GE REPRESENTATIVE. THE GE REPRESENTATIVE WOULD BE REQUIRED TO IDENTIFY THE PHYSICAL CONNECTION LOCATION, AND INSURE PROPER HANDLING OF GE EQUIPMENT.

DIAGRAM KEY

- CUSTOMER/CONTRACTOR SUPPLIED WIRING. ROUTE IN ADEQUATE CONDUIT OR RACEWAY.
- GE FURNISHED CABLE RUNS. ROUTE IN EMPTY CONDUIT OR RACEWAY.
- 59' [18M] MAXIMUM RUN LENGTH BETWEEN JUNCTION POINTS. Feet, [Meters]

GE Healthcare
 Healthcare Project Implementation - Design Center
 Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL SPECIFICATIONS
 MODALITY TYPE: OPTIMA MR360
 THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS AND SPECIFICATIONS. THIS PLAN IS NOT TO BE USED FOR CONSTRUCTION PURPOSES. GE HEALTHCARE SHALL NOT BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
**TYPICAL MR
 8-238-C**
 TYPICAL INSTALLATION DRAWINGS













PROJECT	REVISION
8-238F	03

DATE: 06.SEP.12
 DRAWN BY: PMM
 CHECKED BY: TMS

REVISION HISTORY:

SHEET
E2

PIM R3
 RQ - 124576




ELECTRICAL DETAIL
NETWORK CONNECTION (TYPICAL)

ELEC-84

REV. DATE: 03/06/04

• FOR NUCLEAR SYSTEMS, A DIRECT NETWORK CONNECTION IS TO BE MADE BETWEEN THE SYSTEM AND THE REVIEW WORKSTATION.



DETAIL NOT TO SCALE



GE Healthcare
Healthcare Project Implementation - Design Center
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL DETAILS
MODALITY TYPE: OPTIMA MR360

THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE MANUFACTURER'S RECOMMENDATIONS. THE COMPANY CANNOT BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
**TYPICAL MR
8-238-C**
TYPICAL INSTALLATION DRAWINGS

PROJECT	REVISION
8-238F	03
DATE:	06.SEP.12
DRAWN BY:	PMM
CHECKED BY:	TMS


REVISION HISTORY:

SHEET
E4

RQ - 124576 PIM R3

TYPICAL CRYOGEN VENT PIPE DETAIL

MECH-01



REV. DATE: 05/04/10

CUSTOMER SUPPLIED WARNING SIGN TO READ:

CAUTION
FREEZING GASES AND SMALL OBJECTS MAY BE DISCHARGED WITHOUT NOTICE. STAY AT LEAST 20 FT (6.1 m) AWAY. *35FT (10.7M) FOR 3I

THIS SIGN MUST BE PLACED AT THE EXTERIOR EXIT POINT OF THE CRYOGEN VENT FOR THIS FACILITY. SEE TYPICAL CRYOGEN VENT PIPE DETAIL FOR POSSIBLE WARNING SIGN LOCATIONS.

THE FOLLOWING ARE MATERIALS THAT MUST BE USED FOR CONSTRUCTION OF THE VENT INSIDE THE MAGNET ROOM.

SS 304
AL 6061-T6
CU DWVM OR L

NOTE:
VENTGLASS AND CLAMPS FOR 8 IN. (203mm) DIAMETER PIPE SUPPLIED BY GEHC.

NOTE:
THE VENT GLASS ISOLATION JOINT INSIDE THE MAGNET ROOM MUST BE A MAXIMUM OF 114" (2.95m) ABOVE THE FINISHED FLOOR.

NOTE:
GE SUPPLIES VENTGLASS & CLAMPS WHICH CAN BE USED FOR 8 IN. (203mm) DIAMETER PIPE ONLY. THESE MATERIALS MAY BE USED FOR ISOLATION JOINT OUTSIDE OF ROOM AT THE CONTRACTOR'S OPTION IF THE MATERIALS MEETS THE CONTRACTOR'S DESIGN REQUIREMENTS.

- THE MATING DIAMETERS MUST MATCH WITHIN ±0.125 IN (3mm)
- THE VENTGLASS MUST NOT BE USED FOR STRUCTURAL SUPPORT

CRYOGENIC VENT SYSTEM PRESSURE DROP MATRIX (A)

MECH-04

(THIS TABLE MUST BE USED FOR CRYOGENIC VENT SYSTEM DESIGN)

INSIDE DIAMETER OF VENT PIPE in.(mm)	DISTANCE OF VENT SYSTEM COMPONENT FROM MAGNET SURFACE ft.(m)	PRESSURE DROP PER STRAIGHT VENT PIPE WITH SMOOTH INSIDE SURFACE		PRESSURE DROP PER ELBOW USED ANYWHERE WITHIN 20 FT. VENT SEGMENT			
		psi (kPa/m)	psi (kPa/m)	STANDARD SWEEP ELBOW psi (KPa)	STANDARD SWEEP 90° ELBOW psi (KPa)	LONG SWEEP 90° ELBOW psi (KPa)	LONG SWEEP ELBOW psi (KPa)
8(203)	0-20 (0-6.1)	0.10 (2.26)	1.10 (7.58)	2.06 (14.20)	0.55 (3.79)	1.03 (7.10)	1.85 (12.76)
	20-40 (6.1-12.2)	0.21 (4.75)	2.10 (14.48)	3.70 (25.51)	1.03 (7.10)	1.85 (12.76)	3.36 (23.17)
	40-60 (12.2-18.3)	0.30 (6.79)	2.88 (19.80)	5.21 (35.92)	1.44 (10.33)	2.60 (17.92)	4.65 (32.17)
	60-80 (18.3-24.4)	0.38 (8.60)	3.70 (25.51)	6.71 (46.27)	1.85 (12.76)	3.36 (23.17)	5.85 (40.68)
	80-100 (24.4-30.5)	0.47 (10.63)	4.52 (31.17)	8.22 (56.66)	2.26 (15.58)	4.11 (28.34)	7.15 (49.15)
10(254)	0-20 (0-6.1)	0.03 (0.68)	0.55 (3.79)	0.82 (5.56)	0.27 (1.86)	0.41 (2.83)	0.75 (5.17)
	20-40 (6.1-12.2)	0.07 (1.98)	0.82 (5.56)	1.51 (10.41)	0.41 (2.83)	0.75 (5.17)	1.30 (9.01)
	40-60 (12.2-18.3)	0.10 (2.26)	1.23 (8.48)	2.19 (15.10)	0.62 (4.27)	1.10 (7.58)	1.95 (13.57)
	60-80 (18.3-24.4)	0.12 (2.71)	1.51 (10.41)	2.74 (18.89)	0.75 (5.17)	1.37 (9.45)	2.45 (16.98)
	80-100 (24.4-30.5)	0.16 (3.62)	1.92 (13.24)	3.43 (23.65)	0.96 (6.62)	1.71 (11.79)	3.15 (21.70)
12(305)	0-20 (0-6.1)	0.013 (0.29)	0.27 (1.86)	0.41 (2.83)	0.14 (0.97)	0.21 (1.45)	0.35 (2.41)
	20-40 (6.1-12.2)	0.027 (0.51)	0.41 (2.83)	0.82 (5.56)	0.27 (1.86)	0.41 (2.83)	0.65 (4.48)
	40-60 (12.2-18.3)	0.041 (0.93)	0.55 (3.79)	1.10 (7.58)	0.41 (2.83)	0.65 (4.48)	1.05 (7.24)
	60-80 (18.3-24.4)	0.054 (1.22)	0.69 (4.76)	1.37 (9.45)	0.54 (3.79)	0.82 (5.56)	1.30 (9.01)
	80-100 (24.4-30.5)	0.069 (1.56)	0.86 (5.95)	1.51 (10.41)	0.68 (4.76)	1.04 (7.10)	1.65 (11.37)
14(356)	0-20 (0-6.1)	0.008 (0.22)	0.18 (1.25)	0.27 (1.86)	0.10 (0.68)	0.16 (1.10)	0.25 (1.70)
	20-40 (6.1-12.2)	0.016 (0.44)	0.36 (2.51)	0.54 (3.79)	0.20 (1.37)	0.32 (2.21)	0.50 (3.41)
	40-60 (12.2-18.3)	0.022 (0.61)	0.54 (3.79)	0.82 (5.56)	0.30 (2.07)	0.46 (3.17)	0.70 (4.83)
	60-80 (18.3-24.4)	0.029 (0.81)	0.72 (5.00)	1.09 (7.58)	0.40 (2.74)	0.62 (4.27)	0.95 (6.53)
	80-100 (24.4-30.5)	0.037 (1.03)	0.90 (6.23)	1.36 (9.45)	0.50 (3.41)	0.75 (5.17)	1.15 (7.92)
16(406)	0-20 (0-6.1)	0.005 (0.14)	0.11 (0.75)	0.17 (1.17)	0.07 (0.48)	0.11 (0.75)	0.17 (1.17)
	20-40 (6.1-12.2)	0.010 (0.29)	0.22 (1.51)	0.34 (2.33)	0.14 (0.97)	0.21 (1.45)	0.32 (2.21)
	40-60 (12.2-18.3)	0.015 (0.44)	0.33 (2.33)	0.51 (3.41)	0.21 (1.45)	0.32 (2.21)	0.48 (3.28)
	60-80 (18.3-24.4)	0.020 (0.58)	0.44 (3.11)	0.68 (4.65)	0.28 (1.91)	0.41 (2.83)	0.60 (4.11)
	80-100 (24.4-30.5)	0.026 (0.72)	0.56 (3.92)	0.86 (5.95)	0.36 (2.41)	0.54 (3.79)	0.80 (5.56)

NOTE 1: ELBOWS WITH ANGLES GREATER THAN 90° MUST NOT BE USED.



NOTE 2: THE TABLE DATA IS BASED ON THE FOLLOWING:
A. INITIAL FLOW CONDITIONS AT MAGNET INTERFACE.
B. GAS TEMPERATURE STARTING AT 4.5 KELVIN (-452 F OR -268 F C).
C. HELIUM GAS FLOW RATE OF 2.737 CUBIC FEET (77.5 CUBIC METERS) PER MINUTE.
D. 45° STANDARD SWEEP ELBOW K = 15 F.
E. 90° STANDARD SWEEP ELBOW K = 30 F.
F. 45° LONG SWEEP ELBOW K = 7.5 F.
G. 90° LONG SWEEP ELBOW K = 15 F.

NOTE 3: THE TOTAL PRESSURE DROP OF THE ENTIRE CRYOGENIC VENT SYSTEM MUST BE LESS THAN 17 PSI (117.2 KPa). THE CALCULATION STARTS AT THE MAGNET VENT INTERFACE AND ENDS AT THE TERMINATION POINT OUTSIDE THE BUILDING.

NOTE 4: FOR 14 IN. (356mm) AND 16 IN. (406mm) VENT PIPE DIAMETERS REFER TO PRE-INSTALLATION MANUAL REFERENCED ON SHEET C1.

SYSTEM CHILLER PIPING (COOLING TYPE C)



MECH-05



SCALE: 1/4" = 1'-0"

MECHANICAL/PLUMBING LAYOUT

RECOMMENDED CEILING HEIGHT = 8'-9"



MECHANICAL/PLUMBING ITEMS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	REFER TO PRE-INSTALLATION MANUAL LISTED ON SHEET C1 FOR CRYOGEN VENT REQUIREMENTS. SEE SHEET S-2 FOR CRYOGEN VENT LOCATION. 8" (203 mm) CRYOGEN VENT - TOLERANCE FOR VENT LOCATION +/- 0.25" (6 mm). SEE DETAILS MECH-04 AND MECH-01. THE CUSTOMER'S DESIGNER IS RESPONSIBLE FOR SELECTING VENT MATERIALS AND HARDWARE CAPABLE OF SAFELY HANDLING THE PRESSURES AND COLD TEMPERATURE GENERATED WITHIN THE VENT AT EACH WEJ SITE. THE CUSTOMER'S CONTRACTOR IS RESPONSIBLE FOR PROVIDING AND INSTALLING THE CRYOGEN VENT FROM THE MAGNET VENT ADAPTER TO THE BUILDING'S EXTERIOR. FOR NON-STANDARD VENT CONFIGURATIONS (I.E. OFFSET CEILING EXITS, WALL EXITS, AND GEODESIC DOMES) THE CUSTOMER'S CONTRACTOR IS RESPONSIBLE FOR THE DESIGN AND INSTALLATION OF THE CRYOGENIC VENT SYSTEM AND VENT SUPPORTS WITHIN THE MAGNET ROOM.
2	MINIMUM 2 FT. X 2 FT. (0.61m X 0.61m) PRESSURE EQUALIZING WAVEGUIDE VENT IN THE MAGNET ROOM CEILING.
3	MINIMUM CEILING HEIGHT REQUIREMENT AREA REFER TO MAGNET EQUIPMENT DETAILS FOR MORE INFORMATION.
4	EXHAUST FAN AND AIR INLET MUST BE SIZED FOR A MINIMUM OF 1800 CFM (84 M3/MINUTE) AND A MINIMUM OF 18 AIR EXCHANGES PER HOUR. SEE DETAIL ELEC-55 ON THE ELECTRICAL DETAIL SHEET(S). MAGNET ROOM EXHAUST FAN INTAKE VENT MUST BE LOCATED AT THE HIGHEST CEILING PLANE NEAR THE MAGNET CRYOGEN VENT.
5	SEE PRE-INSTALLATION MANUAL FOR RECOMMENDED BACK-UP WATER SPECIFICATIONS.
6	CUSTOMER/CONTRACTOR RESPONSIBLE FOR RIGGING AND INSTALLATION OF SYSTEM COOLING CABINET. THERE IS A MAXIMUM OF 100 FEET (30.5 M) VERTICAL DIFFERENCE ABOVE OR 10 FEET (3.0 M) BELOW BETWEEN THE OUTDOOR CHILLER CABINET (MRCC) AND THE CRYO COMPRESSOR. A TOTAL MAXIMUM DISTANCE OF 800 FEET (241 M) EXISTS BETWEEN THE OUTDOOR CHILLER CABINET (MRCC) AND CRYO COMPRESSOR OR THE MAGNET. PLEASE REFER TO THE PRE-INSTALLATION MANUAL FOR COMPLETE SITE PREPARATION REQUIREMENTS.
7	TWO (2) 3/4 IN. (19MM) COPPER LINES (INSULATED). FOUR (4) 3/4 IN. (19MM) HOSE BARBS. TWO (2) 1/2 IN. (13MM) HOSE BARBS. FOUR (4) 3/4 IN. (19MM) BALL VALVES. TWO (2) 3/4 IN. (19MM) TO 1/2 IN. (13MM) REDUCERS. ONE (1) 150 MICRON FILTER. ONE (1) SHUT OFF VALVE. ONE (1) BY-PASS VALVE. REFER TO DETAIL MECH-50. PLEASE REFER TO THE PRE-INSTALLATION MANUAL FOR COMPLETE FACILITY WATER REQUIREMENTS.

MECHANICAL/PLUMBING NOTES

- ALL PIPING, FITTINGS, SUPPORTS, HOSES, CLAMPS, VENTILATION SYSTEMS, ETC. ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS.
- FOR COMPLETE DESIGN AND IS REQUIREMENTS, SPECIFICATIONS AND GUIDELINES REFER TO THE PRE-IS MANUAL REFERENCED ON SHEET C1 FOR:
MR SYSTEMS - SYSTEM COOLING, CRYOGEN VENTING, WAVEGUIDES AND EXHAUST VENTING.
CYCLOTRON SYSTEMS - CHEMISTRY LINES, GAS LINES, AND SYSTEM COOLING.

GE Healthcare
Healthcare Project Implementation - Design Center
Milwaukee, Wisconsin

SHEET TITLE: MECHANICAL LAYOUT
MODALITY TYPE: OPTIMA MR360
THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE CODES AND REGULATIONS. THE USER SHALL BE RESPONSIBLE FOR VERIFYING ALL LOCAL AND STATE REQUIREMENTS. GE HEALTHCARE SHALL NOT BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
TYPICAL MR
8-238-C
TYPICAL INSTALLATION DRAWINGS













PROJECT	REVISION
8-238F	03

DATE: 06.SEP.12
DRAWN BY: PMM
CHECKED BY: TMS

REVISION HISTORY:

SHEET
M1

PIM R3
RQ - 124576



PROJECT	REVISION
8-238F	03
DATE:	06.SEP.12
DRAWN BY:	PMM
CHECKED BY:	TMS

REVISION HISTORY:

NO.	DATE	DESCRIPTION

PIM R3
RQ - 124576

