

Drawing Index

These sheets are a document set and should not be separated. Electrical information and references are contained on all sheets.

SITE READINESS	C1
EQUIPMENT LAYOUT	A1
(Equipment locations, heat loads, component weights, environmental specs)	
STRUCTURAL LAYOUT	S1
(Structural support/mounting locations for floor/wall/ceiling, wall support elevations)	
STRUCTURAL DETAILS	S2
(Floor and Ceiling loading information)	
ELECTRICAL LAYOUT	E1
(Contractor supplied wiring, interconnect methods, junction point locations and descriptions)	
ELECTRICAL SPECIFICATIONS	E2
(Maximum wiring run lengths, interconnect diagram, system power specifications)	
ELECTRICAL DETAILS	E3 THRU E4
EQUIPMENT DETAILS	D1 THRU D3

These equipment IS drawings indicate the placement and interconnection of the listed equipment components. These drawings are not construction or site preparation drawings. Customer remains ultimately responsible for preparing the site to accommodate the IS and operation of such equipment in compliance with GE Healthcare's written specifications and all applicable federal, state, and/or local requirements.

* REQUIRED REFERENCE *


Innova Optima
Pre Installation Manual
5400243-2-1EN

A mandatory component of this drawing set is the GE Healthcare Pre Installation manual. Failure to reference the preIS manual will result in incomplete documentation required for site design and preparation.

Pre Installation documents for GE Healthcare products can be accessed on the web at:

www.gehealthcare.com/siteplanning

GE Healthcare



Cardio-Vascular Site Planning

CUSTOMER ACCEPTANCE



imagination at work

Customer Site Readiness Requirements

- Any deviation from these drawings must be communicated in writing to and reviewed by your local GE Healthcare Installation Project Manager prior to making changes.
- Make arrangements for any rigging, special handling, or facility modifications that must be made to deliver the equipment to the installation site. If desired, your local GE Healthcare Installation Project Manager can supply a reference list of rigging contractors.
- New construction requires the following; 1. Secure area for equipment, 2. Power for drills and other test equipment, 3. Capability for image analysis, 4. Restrooms.
- Provide for refuse removal and disposal (e.g. crates, cartons, packing)
- Contact a radiation physicist or consultant to specify radiation containment requirements.

GE Equipment Delivery Requirements

The items on the GE Healthcare Site Readiness Checklist are REQUIRED to facilitate equipment delivery to the IS site. Equipment will not be delivered if these requirements are not satisfied.

GE Healthcare Site Readiness Checklist Rev 18					
<i>Before using this document ensure you have the latest Rev from MyWorkshop on DCCC429762</i>					
GEHC Global Order # :		Customer:			
GEHC PMI :		FE / Installer:			
The customer is responsible for proper site preparation regardless of any GEHC measurements/inspections/assessments.					
Item #	Inspection Date:	Storage Is Item Ready?	PM Is Item Ready?	FE Is Item Ready?	Comments If "N", enter comments or action plan
1					Pre-Project Delivery Requirements: Ensure oxygen venting system is designed and installed with objective evidence that it is compliant with the GEHC Pre-Installation Manual (PIM) requirements, exhaust system is installed and operational, 480V power, and chilled water supply is available that meets system cooling requirements. Electrical connectivity is available for magnet monitoring and phone service is available during delivery.
2					Pre-Install Screen Room Requirements: If Screen Room is tested with objective evidence that it is compliant with GEHC specifications. Do not install using 2-part anchors for pipe system, blowers must be installed by SE vendor using 2-part anchors.
3					Site Regulatory Requirements: Site Drawing Requirements: final version of equipment installation drawings (including red lined version) verified to match actual room and has been provided to installer. X-ray shielding plan and state acknowledgment letter provided to installer for AR, DC, NC, SC, CO & WA.
4					Site Drawing Requirements: Final version of equipment installation drawings (including red lined version) verified to match actual room and has been provided to installer.
5					Surface Penetration Requirements: Customer/contractor scheduled to provide required drilling or cutting into floors, ceilings, and walls. DR surface penetrations permit available and posted in the room when GEHC will perform the work.
6					Delivery Route Requirements: The equipment delivery route from the truck to the final destination within the facility has been reviewed with all key stakeholders to safely meet the minimum requirements for equipment access, and all communications/notifications have occurred. Arrangements have been made for special handling (elevator, rigging, floor protectors, fork lift, rollback truck, etc).
7					Finished Room Requirements: Rooms that will contain equipment, including storage areas, and clean rooms, are dust free. Precautions taken to maintain dust free room. Precautions must be taken to prevent dust from entering rooms containing equipment when construction is incomplete in adjacent areas. All walls primed (final coat not needed on Day 1). No contractor work being done during or after the installation that will cause dust in the installation areas or potential equipment damage. Room security to prevent unauthorized access and theft has been discussed with customer. The customer is aware of these security issues, implications and responsibilities. For Storage Room must meet PIM requirements for storage.
8					Electrical Requirements: Main Disconnect Panel (MDP) is installed and system power is available. Conduits, electrical cable ducting/dividers/cable trays, and access flooring is installed in proper location and height. Surface floor duct and lead-side wires can be installed at time of system installation.
9					HVAC Requirements: The HVAC/Chilled Water systems designed to maintain the environment are running and appear to provide the desired environmental conditions (temperature and humidity) for system operation.
10					Flowing Requirements: Floor is clean and prepared for final floor covering. Floor levelness/flatness is measured and within tolerance, and there are no visible defects per GEHC specifications.
11					Ceiling Requirements: Units for equivalent location, levelness and spacing is measured per vendor confirmed and consistent with the requirements of the installation drawings. Ceiling grids installed. Permanent lighting is installed and operational. HVAC diffusers are installed and connected to ductwork. Ceiling tiles installed per PIM direction.

GE Healthcare

IS Services Design Center

Minneapolis, WI
Copyright © 2009, General Electric Company. Proprietary to GE.

SHEET TITLE: SITE READINESS

MODALITY TYPE: INNOVA 3100-IQ OPTIMA

THIS PLAN IS SUBMITTED TO CURRENT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE LATEST CONSTRUCTION PRACTICES. HOWEVER, THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

**INTERVENTIONAL
CARDIOLOGY – OPTIMA**

TYPICAL FINAL LAYOUT

PROJECT	REVISION
4-64F	00
DATE: 07 Apr 11	
DRAWN BY: LLM	
CHECKED BY: TST	

REVISION HISTORY:

SHEET

C1

PIM R1
RQ - 121212

GE EQUIPMENT LISTING

EQUIPMENT ON ORDER FROM GE HEALTHCARE, INSTALLED BY GE HEALTHCARE, PER : NEITHER A QUOTE OR GON WAS ISSUED AT THE DATE OF THESE DRAWINGS

NOTE: LOCAL CONDITIONS MAY DICTATE THAT ITEMS IDENTIFIED IN THIS CATEGORY BE INSTALLED BY OTHERS.


ITEM NO.	QUANTITY ORDERED	REFER TO SHEET "D"	ITEM DESCRIPTION (* = EXISTING/REINSTALL)	WEIGHT	HEAT OUTPUT (PER HOUR)	DETAIL NO.	STRC PLAN	ELEC PLAN
1	1		WORKSTATION CART				---	---
2	1		MAC LAB CONSOLE, INCLUDES MONITORS AND KEYBOARD	566 lbs	2935 btu		---	PC S
3	1		TRAM NET RACK	8 lbs		B5047	---	TRAM S
4	1		COLOR PRINTER		1054 btu		---	S
5	1		COUNTERBALANCED EYE AND THYROID SHIELD WITH 996 LAMP	143 lbs		B5031E	---	LMP S
6	1		REMOTE CONTROL FOR INJECTOR	4 lbs		B5028	---	IEC S
7	1		INJECTOR HEAD ON TABLE RAIL	15 lbs		B5030A	---	IH S
8	1		INJECTOR ELECTRONICS	37 lbs	320 btu	B5028	---	IE S
9	1		ERGO PRE-CABLED THREE LCD MONITOR SUSPENSION (SELF WEIGHT WITHOUT MONITOR AND ACCESSORIES. MAXIMUM MONITOR WEIGHT IS 4 LBS EACH)			B8301 B8302	---	WBM1 S
10	2		ERGO STATIONARY RAILS (10'-10" (INCLUDES CABLE TRACK))	81 lbs			CBB 300	C
11	1		INNOVA POSITIONER (REFERENCE TABLE BASE-PLATE DETAIL FOR FLOOR MOUNTING INFORMATION)	1653 lbs	2416 btu	B5050A B5050B B5050C B5050E B5050F B5050G B5050H B5050J B5050P B5050R	---	LC1 C
12	1		OMEGA IV/V TABLE WITH ROTATING TOP	1300 lbs	614 btu	B5037	B50 49M	LUS C
13	1		XR BUZZER (LOCATED ABOVE CEILING)	2 lbs		B5150H	---	XR B
14	1		UPS INTERFACE BOX			E4502IB	---	UI B
15	1		ATLAS CABINET (C2)	659 lbs	1825 btu	B0558C	\$100	C2 C
16	1		ATLAS CABINET (C1)	1115 lbs	3389 btu	B0558C	\$100	C1 C
17	1		UPS CABINET	1170 lbs	4061 btu	E4502SC	---	UPS
18	1		DETECTOR CHILLER	33 lbs	706 btu	B5049F	---	DC S
19	1		WATER CHILLER	449 lbs	18716 btu	M0917B	---	CHLR S
20	1		CONTROL ROOM MONITOR WITH DL KEYPAD	22 lbs	204 btu	C7412H C7617	---	S
21	1		OPERATORS CONSOLE	22 lbs	546 btu	C7617 C7502 B5050C	---	WBC1 C
22	1		CAJODD WORKSTATION WITH TWO LCD MONITORS	50 lbs	1040 btu	B8136	---	WS
23	1		BOLUS CHASE HANDSWITCH	2 lbs			---	WBC

THE FOLLOWING ITEMS, WHICH HAVE BEEN ORDERED FROM GE HEALTHCARE, ARE TO BE INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

6	6		VITALING SPEAKER					
1	1		INNOVA MAIN DISCONNECT, REFERENCE JUNCTION POINT "PDB" ON SHEET E1 FOR DETAILED DESCRIPTION.	386 lbs	1532 btu	E4502M	---	PDB
2	2		VITALING MICROPHONE					
1	1		VITALING CONSOLE			B0566		

EQUIPMENT LAYOUT SCALE: 1/4" = 1'-0" RECOMMENDED CEILING HEIGHT = 9'-6"

This equipment layout indicates the placement and interconnection of the indicated equipment components. There may be federal, state, and/or local requirements that could impact the placement of these components. It remains the Customer's responsibility for ensuring the site and final equipment placement complies with all applicable federal, state, and/or local requirements.



ANCILLARY ITEMS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
60	MED GASES IN CEILING
61	CUSTOMER SUPPLIED STORAGE CABINET
62	CATHETER CABINETS
63	LEAD APRON RACK
64	COUNTER TOP WITH SINK, BASE AND WALL CABINETS
65	COUNTER TOP FOR EQUIPMENT - MINIMUM DEPTH 30 IN. OR ADDITIONAL SHELVING MAY BE REQUIRED. PROVIDE GROMMETTED OPENINGS AS REQUIRED TO BE INTERCONNECT CABLES TO RACEWAY BELOW COUNTERTOP.
66	CONTROL WALL TO CEILING WITH LEAD GLASS WINDOW
67	SHELF - CUSTOMER TO PROVIDE ADEQUATE WALL SUPPORT
68	X-RAY ON WARNING LIGHT - AVAILABLE FROM GE SUPPLY CALL 800-200-9760 GE CAT. NO. WX1ABW-DF-XIU
69	MINIMUM DOOR OPENING FOR EQUIPMENT DELIVERY IS 44 IN. W x 83 IN. H. (118mm x 2109mm). CONTINGENT ON A 96 IN. (2438mm) CORRIDOR WIDTH
70	150-AMP LOCAL SERVICE DISCONNECT FOR LOCK-OUT/TAG-OUT CAPABILITY (MAY BE A FUSED DISCONNECT, CIRCUIT BREAKER OR SAFETY SWITCH.)
71	COUNTERTOP WITH BASE AND WALL CABINETS

THE FOLLOWING ITEMS ARE AVAILABLE FROM GE HEALTHCARE TECHNOLOGIES. CONTACT YOUR LOCAL GE HEALTHCARE SERVICE REPRESENTATIVE FOR PRICING AND AVAILABILITY.

90	X-RAY ROOM WARNING LIGHT/ROOM LIGHTING CONTROL PANEL. REFERENCE JUNCTION POINT "XRL" ON SHEET "E1" FOR DETAILED DESCRIPTION -CAT. NO. E4502SS FOR WARNING LIGHT & ROOM LIGHT CONTROL.
----	---

GENERAL SPECIFICATIONS

- THE REQUIRED CEILING HEIGHT INDICATED ON THESE PLANS IS TO ENSURE EQUIPMENT FUNCTION IS NOT INHIBITED. CONSULT WITH YOUR LOCAL GEHC SPECIALIST REGARDING ACCEPTABILITY OF OTHER CEILING HEIGHTS.
- CHECK ALL DOOR OPENINGS AND HALLWAYS FROM DELIVERY LOCATION TO WHERE EQUIPMENT IS TO BE INSTALLED TO ENSURE THE ROUTE PHYSICALLY AND STRUCTURALLY WILL ACCOMMODATE THE EQUIPMENT AS SHIPPED.
- RADIATION PROTECTION REQUIREMENTS ARE NOT INDICATED ON THIS PLAN. WHERE NEEDED PER NATIONAL OR LOCAL CODE THEY SHALL BE SPECIFIED BY A QUALIFIED RADIOLOGICAL PHYSICIST.
- THE DEVELOPMENT OF THE EQUIPMENT LAYOUT, ROOM DIMENSIONS, MECHANICAL AND ELECTRICAL SUGGESTIONS IS PREDICATED UPON THE BEST INFORMATION OBTAINABLE FROM THE SITE, COUPLED WITH THE CUSTOMER'S KNOWN DESIRES. ARCHITECTURAL OR ELECTRICAL CHANGES INCLUDING RELOCATION OF EQUIPMENT ILLUSTRATED ON THIS DRAWING IS ALLOWED ONLY WITH NOTIFICATION, IN WRITING, AND REVIEW BY GEHC SERVICE DEPARTMENT. EQUIPMENT OPERATION, SERVICEABILITY, AND RESTRICTING CABLE LENGTHS, ETC. MAKE THIS ESSENTIAL FOR A PROPER IS. GEHC RESERVES THE RIGHT TO MAKE ON THE JOB CHANGES BECAUSE OF CUSTOMER REQUIREMENTS AND/OR OBSTACLES IN CONSTRUCTION, ETC..
- ALL WORK TO BE IN COMPLIANCE WITH NATIONAL AND LOCAL BUILDING SAFETY CODES.
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.

SITE ENVIRONMENT SPECIFICATIONS

- EQUIPMENT ROOM AMBIENT OPERATING TEMPERATURE: 58 TO 75 DEGREES (F). MAXIMUM ALLOWABLE TEMPERATURE CHANGE OF 15 DEGREES (F)/HOUR, WITH 20% - 75% HUMIDITY.
- EXAM ROOM AMBIENT OPERATING TEMPERATURE: 55 TO 75 DEGREES (F). MAXIMUM ALLOWABLE TEMPERATURE CHANGE OF 15 DEGREES (F)/HOUR, HUMIDITY: 10% - 70%
- CONTROL ROOM AMBIENT OPERATING TEMPERATURE: 59 TO 75 DEGREES (F). MAXIMUM ALLOWABLE TEMPERATURE CHANGE OF 15 DEGREES (F)/HOUR, HUMIDITY: 30% - 80%
- ALTITUDE: NOT TO EXCEED 8,000 FT. ABOVE SEA LEVEL.
- DO NOT RESTRICT THE AIR INTAKE AT THE LOWER FRONT OR AIR EXHAUST AT THE TOP OF THE ELECTRONICS CABINETS.

MAGNETIC INTERFERENCE SPECIFICATIONS

IMAGE INTENSIFIERS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 1 GAUSS TO GUARANTEE SPECIFIED IMAGING PERFORMANCE.

X-RAY TUBES MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE SPECIFIED PERFORMANCE.

SYSTEM ELECTRONICS MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO GUARANTEE DATA INTEGRITY.

OPERATORS CONSOLE EQUIPMENT MUST BE LOCATED IN AMBIENT STATIC MAGNETIC FIELDS OF LESS THAN 10 GAUSS TO OBTAIN SPECIFIED GEOMETRIC LINEARITY.

GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT LAYOUT
MODALITY TYPE: INNOVA 3100-IQ OPTIMA
THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE REGULATORY REQUIREMENTS. THE USER SHALL BE RESPONSIBLE FOR OBTAINING ALL NECESSARY PERMITS AND APPROVALS. GE HEALTHCARE ASSUMES NO LIABILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
**INTERVENTIONAL
CARDIOLOGY - OPTIMA**
TYPICAL FINAL LAYOUT




PROJECT	REVISION
4-64F	00

DATE: 07.Apr.11
DRAWN BY: LLM
CHECKED BY: TST

REVISION HISTORY:

SHEET
A1


TYPICAL WALL SUPPORT ELEVATIONS



SCALE: 1/4" = 1'-0"

STRUCTURAL LAYOUT

RECOMMENDED CEILING HEIGHT = 9'-6"



STRUCTURAL SUPPORT METHODS

CUSTOMER/CONTRACTOR SUPPLIED AND INSTALLED ITEMS

ITEM NO.	ITEM DESCRIPTION (* INDICATES EXISTING)
1	UNISTRUT OR EQUIVALENT SUPPORTS FOR FASTENING THE OVERHEAD COUNTERSUSPENSION. SUPPORTS TO BE LOCATED AS SHOWN. SUPPORT SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE. FLUSH WITH FINISHED CEILING. SUSPENSION REQUIRES 102 LBS/BOLT SUPPORT. METHODS OF SUPPORT THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION.
2	UNISTRUT OR EQUIVALENT SUPPORT IN CEILING FOR FASTENING CEILING SUPPORTED EQUIPMENT. SUPPORTS TO RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL. RUN WALL TO WALL, BE PARALLEL, SQUARE, AND IN THE SAME HORIZONTAL PLANE. FLUSH WITH THE FINISHED CEILING. RAILS ARE MOUNTED TO THESE SUPPORTS EVERY 2'-2" AND REQUIRE 425 LBS. (997 LBS. IN SEISMIC REGIONS) PER BOLT LOAD. METHODS OF SUPPORT THAT PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE SHOULD BE FAVORED. DO NOT USE SCREW ANCHORS IN DIRECT TENSION.
3	AREA OCCUPIED BY GE SUPPLIED POSITIONER BASEPLATE
4	AREA OCCUPIED BY GE SUPPLIED OMEGA TABLE BASE
5	MOUNT XR BUZZER BRACKET ON WALL ABOVE CEILING
6	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S115, FOR UPS INTERFACE BDX.
7	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S100, FOR ATLAS CABINET.
8	SUPPORT BACKING, LOCATE AS SHOWN, REFER TO ELEVATION DETAIL S107, FOR MAIN DISCONNECT CONTROL.

STRUCTURAL NOTES

- ALL STEEL WORK AND PARTS NECESSARY TO SUPPORT CEILING MOUNTED TUBE HANGER OR OTHER EQUIPMENT ARE TO BE SUPPLIED BY THE CUSTOMER OR HIS CONTRACTORS. THE UNISTRUT OR EQUIVALENT STRUCTURE SHOULD RUN CONTINUOUS WITH NO FITTINGS EXTENDING BELOW FACE OF UNISTRUT CHANNEL, RUN WALL TO WALL, BE PARALLEL, SQUARE AND IN THE SAME HORIZONTAL PLANE FLUSH WITH FINISHED CEILING. THE SYSTEM IS TO BE CROSS BRACED VERTICALLY, HORIZONTALLY AND DIAGONALLY TO ALLOW NO MOVEMENT AND A MAXIMUM OF 1.58mm (1/16") DEFLECTION. (10) 12.7mm (1/2") DIA. x 38.1mm (1 1/2") LONG BOLTS WITH UNISTRUT 12.7mm (1/2") NUTS WITH SPRINGS ARE TO BE PROVIDED BY CUSTOMER OR HIS CONTRACTORS FOR EACH STATIONARY AND AUXILIARY SUPPORT RAIL. CLOSURE STRIPS SHALL BE PROVIDED FOR AREAS OF UNISTRUT EXPOSED AND WITHOUT MOUNTING UNITS.
- METHODS OF SUPPORT FOR THE STEELWORK THAT WILL PERMIT ATTACHMENT TO STRUCTURAL STEEL OR THROUGH BOLTS IN CONCRETE CONSTRUCTION SHOULD BE FAVORED. DO NOT USE CONCRETE OR MASONRY ANCHORS IN DIRECT TENSION.
- ALL UNITS THAT ARE WALL MOUNTED OR WALL SUPPORTED ARE TO BE PROVIDED WITH SUPPORTS WHERE NECESSARY. WALL SUPPORTS ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTORS. SEE PLAN AND DETAIL SHEETS FOR SUGGESTED LOCATIONS AND MOUNTING HOLE LOCATIONS.
- ALL CEILING MOUNTED FIXTURES, AIR VENTS, SPRINKLERS, ETC. TO BE FLUSH MOUNTED, OR SHALL NOT EXTEND MORE THAN 6.35mm (1/4") BELOW THE FINISHED CEILING.
- CONTROL WALLS WITH TUBE HANGER PASSAGE ABOVE SHALL BE CONSTRUCTED TO 2130mm (7'-0") HIGH.
- FLOOR SLABS ON WHICH EQUIPMENT IS TO BE INSTALLED MUST BE LEVEL TO 3.17mm (1/8") IN 3050mm (10'-0").
- DIMENSIONS ARE TO FINISHED SURFACES OF ROOM.
- CUSTOMERS CONTRACTOR MUST PROVIDE ALL PENETRATIONS IN POST TENSION FLOORS.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL ANY NON-STANDARD ANCHORING. DOCUMENTS FOR STANDARD ANCHORING METHODS ARE INCLUDED WITH GE EQUIPMENT DRAWINGS FOR GEOGRAPHIC AREAS THAT REQUIRE SUCH DOCUMENTATION.
- CUSTOMERS CONTRACTOR MUST PROVIDE AND INSTALL HARDWARE FOR "THROUGH THE FLOOR" ANCHORING AND/OR ANY BRACING UNDER ACCESS FLOORS. THIS CONTRACTOR MUST ALSO PROVIDE FLOOR DRILLING THAT CANNOT BE COMPLETED BECAUSE OF AN OBSTRUCTION ENCOUNTERED WHILE DRILLING BY THE GE INSTALLER SUCH AS REBAR ETC.
- IT IS THE CUSTOMER'S RESPONSIBILITY TO PERFORM ANY FLOOR OR WALL PENETRATIONS THAT MAY BE REQUIRED. THE CUSTOMER IS ALSO RESPONSIBLE FOR ENSURING THAT NO SUBSURFACE UTILITIES (E.G., ELECTRICAL OR ANY OTHER FORM OF WIRING, CONDUITS, PIPING, DUCT WORK OR STRUCTURAL SUPPORTS (I.E. POST TENSION CABLES OR REBAR)) WILL INTERFERE OR COME IN CONTACT WITH SUBSURFACE PENETRATION OPERATIONS (E.G. DRILLING AND INSTALLATION OF ANCHORS/SCREWS) PERFORMED DURING THE INSTALLATION PROCESS. TO ENSURE WORKER SAFETY, GE INSTALLERS WILL PERFORM SURFACE PENETRATION OPERATIONS ONLY AFTER THE CUSTOMER'S VALIDATION AND COMPLETION OF THE "GE SURFACE PENETRATION PERMIT"

SHEET TITLE: STRUCTURAL LAYOUT
MODALITY TYPE: INNOVA 3100-IQ OPTIMA

PROJECT TITLE:
INTERVENTIONAL
CARDIOLOGY - OPTIMA
TYPICAL FINAL LAYOUT

PROJECT	REVISION
4-64F	00

DATE: 07.Apr.11
DRAWN BY: LLM
CHECKED BY: TST

REVISION HISTORY:

SHEET
S1

PIM R1
RQ - 121212

THIS SHEET IS PART OF THE DOCUMENT SET LISTED ON SHEET C1 AND SHOULD NOT BE SEPARATED


GE Healthcare
IS Services Design Center
Minneapolis, Wisconsin

THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF GE HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE LOCAL, STATE AND FEDERAL REGULATIONS AND CODES. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

SUPPORT DETAIL
Mavig Ceiling Track Mounting

B50-31F

REV. 00: 05/09/05




SUPPORT DETAIL

XT Radiographic Suspension, Optima Mounting

CB8-300


REV. DATE: 12/29/10



FLOOR MOUNTING : ALL INNOVA (UNITY AND HARMONY) SYSTEMS/OMEGA V LONG TABLE (NO IQ TILT TABLE BASEPLATE) INSTALLATION (TEMPLATE NO. 2127792)

B5049M


REV. DATE: 06/04/09




WARNING!! THE RELATIONSHIP BETWEEN THE TABLE BASE AND THE POSITIONER BASEPLATE IS CRITICAL.

PRIOR TO DRILLING MOUNTING HOLES CONTACT LOCAL GE HEALTHCARE INSTALLATION PROJECT MANAGER OR LEAD FIELD ENGINEER TO VERIFY THAT THE PROPER FULL SIZE FLOOR MOUNTING TEMPLATE IS USED.


MEDICAL GAS FLOOR EXIT LOCATIONS



THROUGH-BOLT MOUNTING OPTIONS



ANCHOR BOLT MOUNTING OPTIONS



Customer/Contractor Alert: It is the responsibility of the Customer or their Contractor to drill all anchor/thru-bolting holes for anchoring the positioner and table to the floor. Refer to GEHC document no. *2290880-2-100 for installation preparation and procedures.

NOTE: THRU BOLTING IS HIGHLY PREFERRED FOR THE INSTALLATION OF THE POSITIONER BASEPLATE AND OMEGA TABLE. HARDENED BOLTS AND 4" x 4" [102mm x 102mm] STEEL PLATES TO BE USED ARE SUPPLIED BY GE HEALTHCARE AS INDICATED ON THE ACTUAL DETAIL DRAWING. BE ADVISED, HOWEVER, THAT ADDITIONAL SUPPORT STRUCTURES: STEEL BEAMS, PLATES, CORE BORING OF MOUNTING HOLES, ETC., ARE TO BE SUPPLIED AND INSTALLED BY THE CUSTOMER OR HIS CONTRACTOR.

NOTE: IF THRU BOLTING IS NOT POSSIBLE, FLOOR ANCHORS CAN BE USED IF APPROVED BY CUSTOMERS STRUCTURAL ENGINEER. FOR ON GRADE INSTALLATIONS, MOUNTING KIT CAT. NO. **2286398** SHOULD BE ORDERED. ANCHORS INCLUDED IN KIT SHOULD BE APPROVED BY CUSTOMERS STRUCTURAL ENGINEER.

NOTE: BASEPLATES MUST BE LEVEL WITHIN 1/32" [0.79mm]

NOTE: JOISTS MUST BE SPANNED WITH STEEL REINFORCING. SIZE AND THICKNESS OF STEEL REINFORCING ARE DETERMINED BY THE ACTUAL PAN CONSTRUCTION ON SITE. STEEL PLATES, CHANNELS OR BEAMS MAY BE USED. **NOTE:** DETERMINE THE POSITION OF THE "REBARS IN THE CONCRETE FLOOR SO ANCHOR HOLES WILL NOT RUN INTO THEM.

POSITIONER BOLT FORCES FOR WORST CASE CONDITIONS

LOADS	BOLT TENSION (AT BOLT "A")
HORIZONTAL ACCELERATION = 625 lbs. [284 Kg]	MAXIMUM TENSION = 881 lbs. [400 Kg]
VERTICAL ACCELERATION = 209 lbs. [95 Kg]	BOLT SHEAR (U-ARM LOCKED)
	MAXIMUM SHEAR = 120 lbs. [54 Kg]/BOLT

OMEGA TABLE BOLT FORCES FOR WORST CASE CONDITIONS

LOADS	BOLT TENSION	BOLT SHEAR
	MAXIMUM TENSION = 1938 lbs. [880 Kg]/BOLT	MAXIMUM SHEAR = 407 lbs. [185 Kg]/BOLT

SHEET TITLE: STRUCTURAL DETAILS
MODALITY TYPE: INNOVA 3100-IQ OPTIMA

THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN. EVERY EFFORT HAS BEEN MADE TO CONFORM TO ALL APPLICABLE CODES AND REGULATIONS. THIS PLAN IS TO BE USED FOR CONSTRUCTION PURPOSES. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
**INTERVENTIONAL
CARDIOLOGY - OPTIMA**
TYPICAL FINAL LAYOUT

PROJECT REVISION
4-64F 00

DATE: 07.Apr.11
DRAWN BY: LLM
CHECKED BY: TST

REVISION HISTORY:


SHEET

S2

PIM R1
RQ - 121212

INTERCONNECT DIAGRAM

TYPICAL VIEWS



EQUIPMENT DESCRIPTIONS

ITEM	DESCRIPTION	WEIGHT (lb)	HEAT DISSIPATION (btu)	DRAWING DESIGNATOR
1	XR BUZZER	2		XRB
2	ATLAS CABINET C2	659	1825	C2
3	ATLAS CABINET C1	1115	3389	C1
4	DETECTOR CONDITIONER	33	706	DC
5	WATER CHILLER	449	18716	CHLR
6	20kva UPS CABINET (OPTIONAL)	1170	4061	UPS
7	UPS INTERFACE BOX (OPTIONAL)			UIB
8	ERGO TV CEILING SUSPENSION (3 MONITOR)	208	614	WBM1
9	INNOVA V CEM POSITIONER	1653	2416	LC1
10	OMEGA V LONG TABLE	1750	614	LU5
11	INNOVA VCIM HITH DL KEYBOARD CONSOLE	22	204	
12	VCIM OPERATOR CONSOLE	22	546	WBC1
13	ROOM LIGHTS			RML1
14	XRAY WARNING LAMP			XRL1
15	XRAY WARNING LAMP CONTROLLER			XRLC
16	RDS1 PUSHBUTTON			RDS1
17	RDS2 PUSHBUTTON			RDS2
18	PDB MAIN DISCONNECT	326	1532	PDB
19	LOTO DISCONNECT BREAKER			PDB1

OPTIONS

ITEM	DESCRIPTION	WEIGHT (lb)	HEAT DISSIPATION (btu)	DRAWING DESIGNATOR
21	BOLUS CHASE HANDSWITCH	2		WBBC
22	ADVANTAGE WINDOWS WORKSTATION	81	1201	AW
23	IVUS VOLCANO CONSOLE	68	1631	IVUS
24	IVUS VOLCANO COLOR PRINTER		X	
25	INJECTOR HEAD	15		IH
26	INJECTOR ELECTRONICS	37	320	IE
27	REMOTE CONTROL FOR INJECTOR	4		IEC
28	LAMP (RADIATION SHIELD TRACK)	143		LMP
29	NOT USED			
30	NOT USED			
31	MACH 3 TRANSFORMER	70	X	M3T
32	MACLAB PHYSIO. MONITORING	566	2935	PC
33	PRINTER (PHYSIO.)	X	309	
34	TRAM (PHYSIO.)	8	X	TRAM
35	REMOTE OPERATING TERMINAL (PHYSIO.)	46	682	ROT
36	MICRO PACE (PHYSIO.)	X	X	MP
37	SKYTRON LIGHTING UNIT	50	341	SL
38	150 KVA UPS	2160	31802	UPS
39	UPS BATTERY CABINET	3529	X	
40	MAINTENANCE BYPASS PANEL	350	X	MBP

POWER SPECIFICATIONS

INNOVA SYSTEMS
REV. DATE: 01/04/07

VOLTAGE PRIMARY SOURCE IS REQUIRED FOR ALL INSTALLATIONS.
RANGE OF LINE VOLTAGES
NOMINAL LINE VOLTAGE OF 360 TO 480, 3 PHASE, 50 OR 60 Hz

REQUIRED POWER SUPPLY: WYE DISTRIBUTION

MAXIMUM DAILY VOLTAGE VARIATION MUST FALL WITHIN ONE OF THE RANGES IN TABLE A.

TABLE A ALLOWABLE INPUT VOLTAGES/CURRENT DEMAND

NOMINAL VOLTAGE	NORMAL RANGE ±10 PERCENT	CURRENT (AMPS)	
		MAX. MOMENTARY	CONTINUOUS
360	324-396	304	32
380	342-418	289	31
400	360-440	274	29
420	378-462	264	28
440	396-484	249	26
460	414-506	238	25
480	432-528	228	24

ALL CALCULATIONS BASED UPON NOMINAL VOLTAGE

NOTE LOW LINE CONDITIONS MAY INHIBIT SOME HIGH KVp TECHNIQUES. THE GENERATOR AUTOMATICALLY ESTABLISHES THESE INHIBITS BASED ON ACTUAL LINE CONDITIONS AND SYSTEM REGULATION.

PHASE-BALANCE. PHASE-TO-PHASE VOLTAGES MUST BE WITHIN +2 PERCENT OF THE LOWEST PHASE-TO-PHASE VOLTAGE. MAXIMUM ALLOWABLE TRANSIENT VOLTAGE EXCURSIONS ARE 2.5 PERCENT OF RATED LINE VOLTAGE AT A MAXIMUM DURATION OF 5 CYCLES AND FREQUENCY OF 0 TIMES PER HOUR.

POWER DEMAND CONTINUOUS POWER DEMAND = 20KVA. (MAX DEMAND = 171 KVA)

TABLE B MAXIMUM MOMENTARY POWER DEMAND.

DEMAND	ADVANTX 100
kva * POWER FACTOR AT	171 0.9
mA	1250
kvP	80



* DEMAND INCLUDES POWER FOR ENTIRE ADVANTX SYSTEM. LINE VOLTAGE REGULATION AT MAXIMUM POWER DEMAND MUST BE LESS THAN OR EQUAL TO 6 PERCENT.

DISTRIBUTION TRANSFORMER FOR A SINGLE UNIT INSTALLATION, THE MINIMUM TRANSFORMER SIZE IS 225 KVA.

ELECTRICAL NOTES

- NOTE 1: ALL WIRES SPECIFIED SHALL BE COPPER STRANDED, FLEXIBLE, THERMO-PLASTIC, COLOR CODED, CUT 10 FOOT LONG AT OUTLET BOXES, DUCT TERMINATION POINTS OR STUBBED CONDUIT ENDS. ALL CONDUCTORS, POWER, SIGNAL AND GROUND, MUST BE RUN IN A CONDUIT OR DUCT SYSTEM. ELECTRICAL CONTRACTOR SHALL RING OUT AND TAG ALL WIRES AT BOTH ENDS. WIRE RUNS MUST BE CONTINUOUS COPPER STRANDED AND FREE FROM SPLICES. ALUMINUM OR SOLID WIRES ARE NOT ALLOWED.
- NOTE 2: WIRE SIZES GIVEN ARE FOR USE OF EQUIPMENT. LARGER SIZES MAY BE REQUIRED BY LOCAL CODES.
- NOTE 3: IT IS RECOMMENDED THAT ALL WIRES BE COLOR CODED, AS REQUIRED IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 4: CONDUIT SIZES SHALL BE VERIFIED BY THE ARCHITECT, ELECTRICAL ENGINEER OR CONTRACTOR, IN ACCORDANCE WITH LOCAL OR NATIONAL CODES.
- NOTE 5: CONVENIENCE OUTLETS ARE NOT ILLUSTRATED. THEIR NUMBER AND LOCATION ARE TO BE SPECIFIED BY OTHERS. LOCATE AT LEAST ONE CONVENIENCE OUTLET CLOSE TO THE SYSTEM CONTROL, THE POWER DISTRIBUTION UNIT AND ONE ON EACH WALL OF THE PROCEDURE ROOM. USE HOSPITAL APPROVED OUTLET OR EQUIVALENT.
- NOTE 6: GENERAL ROOM ILLUMINATION IS NOT ILLUSTRATED. CAUTION SHOULD BE TAKEN TO AVOID EXCESSIVE HEAT FROM OVERHEAD SPOTLIGHTS. DAMAGE CAN OCCUR TO CEILING MOUNTING COMPONENTS AND WIRING IF HIGH WATTAGE BULBS ARE USED. RECOMMEND LOW WATTAGE BULBS NO HIGHER THAN 75 WATTS AND USE DIMMER CONTROLS (EXCEPT MR). DO NOT MOUNT LIGHTS DIRECTLY ABOVE AREAS WHERE CEILING MOUNTED ACCESSORIES WILL BE PARKED.
- NOTE 7: ROUTING OF CABLE DUCTWORK, CONDUITS, ETC., MUST RUN DIRECT AS POSSIBLE OTHERWISE MAY RESULT IN THE NEED FOR GREATER THAN STANDARD CABLE LENGTHS (REFER TO THE INTERCONNECTION DIAGRAM FOR MAXIMUM USABLE LENGTHS POINT TO POINT).
- NOTE 8: CONDUIT TURNS TO HAVE LARGE, SWEEPING BENDS WITH MINIMUM RADIUS IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES.
- NOTE 9: A SPECIAL GROUNDING SYSTEM IS REQUIRED IN ALL PROCEDURE ROOMS BY SOME NATIONAL AND LOCAL CODES. IT IS RECOMMENDED IN AREAS WHERE PATIENTS MIGHT BE EXAMINED OR TREATED UNDER PRESENT, FUTURE, OR EMERGENCY CONDITIONS. CONSULT THE GOVERNING ELECTRICAL CODE AND CONFER WITH APPROPRIATE CUSTOMER ADMINISTRATIVE PERSONNEL TO DETERMINE THE AREAS REQUIRING THIS TYPE OF GROUNDING SYSTEM.
- NOTE 10: THE MAXIMUM POINT TO POINT DISTANCES ILLUSTRATED ON THIS DRAWING MUST NOT BE EXCEEDED.
- NOTE 11: PHYSICAL CONNECTION OF PRIMARY POWER TO GE EQUIPMENT IS TO BE MADE BY CUSTOMERS ELECTRICAL CONTRACTOR WITH THE SUPERVISION OF A GE REPRESENTATIVE. THE GE REPRESENTATIVE WOULD BE REQUIRED TO IDENTIFY THE PHYSICAL CONNECTION LOCATION, AND INSURE PROPER HANDLING OF GE EQUIPMENT.

FULL SIZE UPS (OPTIONAL)



GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin

SHEET TITLE: ELECTRICAL SPECIFICATIONS
MODALITY TYPE: INNOVA 3100-IQ OPTIMA

THIS PLAN IS SUBMITTED TO SUBMIT TO THE ELECTRICAL WIRING OF THE HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE REQUIREMENTS OF THE NATIONAL ELECTRICAL CODE AND LOCAL CODES. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
INTERVENTIONAL
CARDIOLOGY - OPTIMA
TYPICAL FINAL LAYOUT

PROJECT	REVISION
4-64F	00

DATE: 07.Apr.11
DRAWN BY: LLM
CHECKED BY: TST

REVISION HISTORY:

SHEET
E2

ELECTRICAL DETAIL
INSITE CONNECTION (TYPICAL)


ELEC-1
REV. DATE: 04/24/02

ONE OF THE FOLLOWING TWO SELECTIONS MUST BE INSTALLED AT THE LOCATION SHOWN ON THE ELECTRICAL PLAN (SHEET E1) FOR GE INSITE CONNECTION BASED UPON SYSTEM CONFIGURATION.

A) ONE INTERNET ACCESSIBLE VIRTUAL PRIVATE NETWORK (VPN) CONNECTION WITH A STATIC IP ADDRESS, AND ONE TELEPHONE LINE - DEDICATED-DIRECT-DIALING, VOICE GRADE.

OR

B) TWO TELEPHONE LINES - ONE DEDICATED DIRECT-DISTANCE-DIALING, VOICE GRADE AND ONE A DEDICATED DATA LINE.




ALL ITEMS ILLUSTRATED ARE TO BE FURNISHED AND INSTALLED BY CUSTOMER OR THEIR CONTRACTOR.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
FLOOR BOX WITH NIPPLES (TYPICAL)


ELEC-13
REV. DATE: 09/30/94



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
BOX WITH COVERPLATE AND NETWORK JACK


ELEC-83
REV. DATE: 10/06/98



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
NETWORK CONNECTION (TYPICAL)

ELEC-84
REV. DATE: 03/06/04




DETAIL NOT TO SCALE

FOR NUCLEAR SYSTEMS A DIRECT NETWORK CONNECTION IS TO BE MADE BETWEEN THE SYSTEM AND THE REVIEW WORKSTATION.

ELECTRICAL DETAIL
HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5
REV. DATE: 03/19/04




DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" [610mm]	2
18" [457mm]	2
10" [254mm]	2
6" [152mm]	1
4" [102mm]	1

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
VERTICAL WALL DUCT (TYPICAL)

ELEC-6
REV. DATE: 03/19/04




DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" [610mm]	2
18" [457mm]	2
10" [254mm]	2
6" [152mm]	1
4" [102mm]	1

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
BOX WITH COVERPLATE (TYPICAL)


ELEC-8
REV. DATE: 09/30/94



DETAIL NOT TO SCALE


ELECTRICAL DETAIL
X-RAY WARNING LIGHT & ROOM LIGHT CONTROL PANEL

ELEC-157
REV. DATE: 04/23/09



ELECTRICAL DETAIL
EMERGENCY OFF BUTTON


ELEC-16
REV. DATE: 05/14/09



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
HORIZONTAL WALL DUCT (TYPICAL)

ELEC-5A
REV. DATE: 06/16/08




DUCT WIDTH	MINIMUM DIVIDERS REQUIRED
24" [610mm]	2
18" [457mm]	2
10" [254mm]	2
6" [152mm]	1
4" [102mm]	1

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
POSITIONER INTERCONNECT DETAIL, UNDER FLOOR

ELEC-100
REV. DATE: 03/30/04




NOTE: PIPE, JUNCTION BOX AND DUCT OR CONDUIT ARE TO BE SUPPLIED AND INSTALLED BY CUSTOMER OR CUSTOMER'S CONTRACTOR.

DETAIL NOT TO SCALE

ELECTRICAL DETAIL
CONDUITS THRU-FLOOR (TYPICAL)

ELEC-9
REV. DATE: 08/08/94



DETAIL NOT TO SCALE

SHEET TITLE: ELECTRICAL DETAILS
MODALITY TYPE: INNOVA 3100-iQ OPTIMA
THIS PLAN IS SUBMITTED TO SUBJECT LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS IN PREPARING THIS PLAN. HOWEVER, THE USER SHALL BE RESPONSIBLE FOR VERIFYING ALL CONDITIONS, DIMENSIONS, AND MATERIALS TO BE USED FOR EACH INSTALLATION. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
INTERVENTIONAL
CARDIOLOGY - OPTIMA
TYPICAL FINAL LAYOUT

PROJECT	REVISION
4-64F	00

DATE: 07.Apr.11
DRAWN BY: LLM
CHECKED BY: TST


REVISION HISTORY:

SHEET
E3

GE Healthcare
IS Services Design Center
Minneapolis, Wisconsin

ELECTRICAL DETAIL
TABLE INTERCONNECTION - BOX BELOW FLOOR


ELEC-48
REV. DATE: 01/04/98



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
J.B. / WALL DUCT DETAIL (TYPICAL)


ELEC-2
REV. DATE: 09/30/94



DETAIL NOT TO SCALE

ELECTRICAL DETAIL
INNOVA PLUS MAIN DISCONNECT PANEL

ELEC-161
REV. DATE: 09/27/10



DETAIL NOT TO SCALE

SHEET TITLE: ELECTRICAL DETAILS
MODALITY TYPE: INNOVA 3100-IQ OPTIMA

THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO DETAILS AND SPECIFICATIONS OF THE MANUFACTURER. THE COMPANY CANNOT BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
INTERVENTIONAL
CARDIOLOGY - OPTIMA
TYPICAL FINAL LAYOUT

PROJECT	REVISION
4-64F	00
DATE:	07.Apr.11
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:

SHEET
E4

RQ - 121212 PIM R1

IS Services Design Center
Minneapolis, Wisconsin

EQUIPMENT DETAIL
TRAM-RAC 4A

B5047
REV. DATE: 05/26/04

TABLE RAIL MOUNT

FLOOR MOUNT

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
MAVIG EYE & THYROID SHIELD WITH LAMP

B50-31E
REV. 00: 10/03/97

CEILING

CARRIAGE TRACK

65" MAX. ARC [1651mm]

65" MAX. ARC [1651mm]

LEAD GLASS SHIELD

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INJECTOR REMOTE CONTROL AND ELECTRONICS

B50-28

REMOTE CONTROL

7.5" [191mm]

14.5" [368mm]

3.25" [83mm]

ELECTRONICS

10.5" [267mm]

17" [432mm]

12.75" [324mm]

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INJECTOR ON TABLE RAIL

B50-30A

PLAN VIEW

SIDE VIEW

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INNOVA OPTIMA - ERGO SUSPENSION PLACEMENT

B8301
REV. DATE: 12/29/10

HEAD

LEFT

RIGHT

FOOT

ISO-CENTER

51.2" [1300mm]

10'-10" (2000mm)

43.3" [1100mm]

NOTE: CLINICAL APPLICATIONS CAN BE MOSTLY COVERED BY THE ABOVE SUGGESTED INSTALLATION POSITION. THE TRAVEL RANGE OF THE PULLEY ALONG THE RAILS IS ABOUT 2.5M (8'-2").

PLAN VIEW

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA OPTIMA - ERGO SUSPENSION DETAIL

B8302
REV. DATE: 12/29/10

FINISHED CEILING

26" (660mm)

330°

330°

Min. 5.9" (149mm)

Max. 34.2" (869mm)

Min. 50.7" (1287mm)

Max. 75" (2007mm)

60" (1512mm)

PLAN VIEW

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA 2100/3100/4100 VASCULAR SYSTEM

B5050A
REV. DATE: 06/07/05

FINISHED CEILING

50.1" [1273mm]

88" [2234mm]

9'-6" [2896mm] RECOMMENDED

46" [1168mm]

42.1" [1070mm]

27.7" [705mm]

4.72" [120mm]

10.7" [272mm]

50°

45°

ISO-CENTER

FOCAL SPOT

FINISHED FLOOR

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA 2100/3100/4100 VASCULAR SYSTEM

B5050B
REV. DATE: 06/07/05

ISO-CENTER

117°

105°

88" [2234]

42.1" [1070mm]

FRONT VIEW

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
INNOVA 3100/4100 VASCULAR SYSTEM

B5050
REV. DATE: 02/22/05

95°

95°

PLAN VIEW

SHIPPING DIMENSIONS:
110"L x 45.5"W x 77"H
[2790mmL x 1160mmW x 1950mmH]
(ON DOLLY)
WIDTH IS REDUCED TO 34" [865mm]
BY REMOVING SIDE RAILS

DRAWING NOT TO SCALE

EQUIPMENT DETAIL
SHIPPING CREATE FOR OMEGA & INNOVA IQ TABLES

B5050E
REV. DATE: 12/07/09

OMEGA IV/V TABLE

76.75" (1950mm)

OMEGA TABLE BASE ASSEMBLY ON PALLET

39" (810mm)

OMEGA 4/5 AND INNOVA IQ TILT TABLE

137" (3470mm)

OMEGA 4/5 and INNOVA IQ TABLE TOP ASSEMBLY ON PALLET

9" (220mm)

93" (2365 mm)

47.20" (1200mm)

33.50" (850mm)

INNOVA IQ TILT TABLE

INNOVA IQ TILT TABLE BASE ASSEMBLY ON PALLET

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
INNOVA C2 CABINET ON SHIPPING DOLLY

B5050F
REV. DATE: 12/07/09

INNOVA C2 CABINET (ON DOLLY)

APPROX. 77" (1950 mm)

52" (1320 mm)

29" (740 mm)

SHIPPING WEIGHT: 821 LBS. (372 kg)(C2-HARMONY SYSTEMS)
SHIPPING WEIGHT: 785 LBS. (356 kg)(C2-BIPLANE SYSTEMS)

DETAIL NOT TO SCALE

EQUIPMENT DETAIL
SHIPPING DOLLY FOR INNOVA LC POSITIONER

B5050G
REV. DATE: 12/07/09

INNOVA LC POSITIONER (ON DOLLY)

110" (2790mm) (maximum)

34" (865mm)

11.5" (290mm)

77" (1950mm)

10" (250mm)

86.22" (2190mm) (minimum)

13.78" (350mm)

NOTE:
BOTH ENDS OF THE DOLLY CAN BE REMOVED WHICH WILL SHORTEN LC GANTRY DOLLY DONE TO 86.22" (2190mm) RECOMMEND ONLY ONE SIDE BE REMOVED WHEN DELIVERY THROUGH HOSPITAL.

SHIPPING WEIGHT: 2340 lbs. (1060 kg)

DETAIL NOT TO SCALE

GE Healthcare
IS Services Design Center
Minneapolis, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: INNOVA 3100-IQ OPTIMA













THIS PLAN IS SUBMITTED TO SUPPORT LOGGING OF HEALTHCARE EQUIPMENT AND ASSOCIATED ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE INSTRUCTIONS AND DIMENSIONS PROVIDED. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:
**INTERVENTIONAL
CARDIOLOGY - OPTIMA**
TYPICAL FINAL LAYOUT

PROJECT	REVISION
4-64F	00
DATE:	07.Apr.11
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:

SHEET
D1



GE Healthcare
IS Services Design Center
Milwaukee, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS
MODALITY TYPE: INNOVA 3100-iQ OPTIMA

THIS PLAN IS SUBMITTED TO ASSIST IN THE CORRECT LOCATION OF THE HEALTH CARE EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS IN PREPARING THIS PLAN, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE DETAILS AND DIMENSIONS SHOWN. HOWEVER, THE USER SHALL BE RESPONSIBLE FOR VERIFYING THE ACTUAL CONSTRUCTION DIMENSIONS. THE COMPANY SHALL NOT BE HELD RESPONSIBLE FOR ANY DAMAGES RESULTING THEREFROM.






PROJECT TITLE:
INTERVENTIONAL
CARDIOLOGY - OPTIMA
TYPICAL FINAL LAYOUT

PROJECT	REVISION
4-64F	00

DATE: 07.Apr.11
DRAWN BY: LLM
CHECKED BY: TST

REVISION HISTORY:

SHEET
D2



GE Healthcare

IS Services Design Center

Minneapolis, Wisconsin

SHEET TITLE: EQUIPMENT DETAILS

MODALITY TYPE: INNOVA 3100-IQ OPTIMA

THIS PLAN IS SUBMITTED TO SUBMIT LOCATION OF HEALTHCARE EQUIPMENT AND ASSOCIATED APPARATUS. EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS IN PREPARING THIS PLAN. HOWEVER, IT IS NOT TO BE USED FOR CONSTRUCTION PURPOSES WITHOUT THE CONSULTATION OF THE COMPANY. THE COMPANY CANNOT ACCEPT RESPONSIBILITY FOR ANY DAMAGES RESULTING THEREFROM.

PROJECT TITLE:

**INTERVENTIONAL
CARDIOLOGY - OPTIMA**

TYPICAL FINAL LAYOUT

PROJECT	REVISION
4-64F	00
DATE:	07.Apr.11
DRAWN BY:	LLM
CHECKED BY:	TST

REVISION HISTORY:

SHEET

D3

RQ - 121212 PIM R1