



Tec 850 Vaporizer

Deliver with confidence



The Tec 850 vaporizer from GE Healthcare combines clinical performance with ergonomic design to help you efficiently deliver anesthesia for every patient. Enhanced design delivers consistent, reliable operation while reducing overall operating costs.

Efficient and effective anesthesia delivery

Designed for clinical confidence

The Tec 850 vaporizer is designed to deliver consistent performance with improved operating excellence throughout the product life so clinicians can focus on patients, not the vaporizer.

New design features include:

- An ergonomic dial release allows either left- or right-hand operation
- Wide, centered liquid level indicator provides a clear indication of the fill level of the vaporizer
- Easy-to-turn dial and fine graduations with a larger dial font for enhanced readability and easily control agent delivery

Lower costs and flexible service

- The Tec 850 vaporizer requires an inspection once every five years and no planned factory maintenance reducing your total cost of ownership
- It comes with a comprehensive five-year warranty. Plus, you get access to a wide range of flexible service offerings provided at your site or at one of our service centers*

**Contact your service representative for a comprehensive list of offerings in your region.*



Uncompromised anesthesia delivery

The only vaporizer validated by GE to integrate seamlessly for GE anesthesia machine.

They feature a non-spill system that limits the movement of liquid agent if the vaporizer is tilted or inverted, protecting internal components and helping maintain output within clinically acceptable limits. Plus, they can be mounted without tools on GE Healthcare anesthesia systems equipped with the Selectatec manifold and without taking the anesthesia system out of service.



Partner with an experienced leader

With 100 years of experience in anesthesia delivery, GE Healthcare has built a reputation for pioneering design and patient safety.

Vaporizer sold in more than

150 countries



Over 100 years of anesthesia innovation



Estimated over 250,000 GE Tec **vaporizers in use**



GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter – great people and technologies taking on tough challenges.

From medical imaging, software & IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

Imagination at work

© 2017 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information. GE, GE Monogram, Easy-Fil, Tec and Selectatec are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company. GE Medical Systems, Inc., doing business as GE Healthcare.

All other product names and logos are trademarks or registered trademarks of their respective companies.

JB51513XX